

# Siluron®

Ultrapurified silicone oils  
for intraocular use

new!



**Geuder®**  
Precision made in Germany

**FLUORON®**  
Leading in purity and variety

# The innovative silicone oils with the special molecular design.

- High resistance to emulsification
- Short injection time
- Exceptional long-term tolerance
- Excellent chemical purity

Innovative molecular design for a new generation of silicone oils

Due to their special molecular structure, the Siluron® Xtra and Siluron® 2000 have, compared to the standard silicone oils, in vitro a high emulsification resistance. Both innovative silicone oils consist of a mixture of ultra-long molecular chains having a viscosity in the range of 2,500,000 mPas and of short molecular chains with a viscosity in the range of 1,000 mPas.

Thanks to this special molecular design, Siluron® Xtra and Siluron® 2000 can modify their viscous properties depending on the permanent high shear forces, such as they appear in the eye due to its constant movement: the greater the applied shear force is, the more viscous behaves the silicone oil, i.e. the

more resistant it is to emulsification. By contrast, the viscosity of the conventional oils decreases continuously as a result of permanently acting external forces, which causes a higher tendency for emulsification.

The new generation of silicone oils – Siluron® 2000 and Siluron® Xtra – is characterized by its special property of a significantly higher emulsification resistance. This is based on an intelligent mixture of different long chains of molecules and the resulting dynamic viscosity. The good injectability in cases of small incisions is a further advantage of innovative silicone oils.

Literature: Chan YK., Ng CO., Knox PC., Garvey MJ., Williams RL., Wong D.: Emulsification of silicone oil and eye movements; Invest Ophthalmol Vis Sci. 2011; 52: 9721-9727

**Siluron® 2000**    **G-80740** Siluron® 2000 syringe 10 ml, sterile  
**Siluron® XTRA**    **G-80750** Siluron® Xtra syringe 10 ml, sterile



Siluron® 1000 | Siluron® 5000

# The proven standard silicone oil tamponades.

Exceptional long-term tolerance  
Excellent chemical purity



Standard  
silicone oils  
also available  
as vial

**Siluron® 1000** G-80710 Siluron® 1000 vial 10 ml, sterile | G-80720 Siluron® 1000 syringe 10 ml, sterile  
**Siluron® 5000** G-80810 Siluron® 5000 vial 10 ml, sterile | G-80820 Siluron® 5000 syringe 10 ml, sterile



The proven standard  
silicone oil tamponade



The proven standard silicone oil  
long-term tamponade

# Overview of Properties

## Physicochemical properties of Siluron® oils

Property	Siluron® 1000	Siluron® 5000	Siluron® 2000	Siluron® Xtra
Density [g/cm³] 25 °C	0.97	0.97	0.97	0.97
Viscosity [mPas] 25 °C	900 - 1200	4800 - 5500	2000 - 2400	4,100 - 4,800
Refractive index	1.404	1.404	1.404	1.404
Solubility in water	non miscible	non miscible	non miscible	non miscible
Composition [w %]	100 % Poly-dimethylsiloxan	100 % Poly-dimethylsiloxan	95 % Siluron 1,000 + 5 % 2.5 Mio. mPas	90 % Siluron 1,000 + 10 % 2.5 Mio. mPas
Elasticity [Je°] [Pas]	2 x 10 <sup>-5</sup>	1 x 10 <sup>-5</sup>	6.5 x 10 <sup>-4</sup>	1.4 x 10 <sup>-3</sup>
Shear viscosity (at 8,37 s <sup>-1</sup> , 37 °C) [mPas]	931	4,303	1,800	4,377
Volatile components (200 °C, 24 h) [%]	≤ 0.2 %	≤ 0.2 %	≤ 0.2 %	≤ 0.2 %


Source: Caramoy A., Hagedorn N., Fauser S., Kugler W., Gross T., Kirchhof B.: Development of emulsification-resistant silicone oils: can we go beyond 2000 mPas silicone oil? Invest Ophthalmol Vis Sci. 2011; 52: 5432-5436

## Comparison of emulsification rate



Source: Caramoy A., Hagedorn N., Fauser S., Kugler W., Gross T., Kirchhof B.: Development of emulsification-resistant silicone oils: can we go beyond 2000 mPas silicone oil? Invest Ophthalmol Vis Sci 2011; 52: 5432-5436


## Comparison of emulsification resistance



Source: Wong D., Stappler T. et al. in preparation

## Comparison of injection time

5.5 bar injection pressure, 20 gauge injection cannula



Source: Williams RL., Day MJ., Garvey MJ., Morphis G., Irigoyen C., Wong D., Stappler T.: Br J Ophthalmol. 2011; 95: 273-276



# The perfect accessories for Siluron

**G-34497**

## Single-use cannula

to inject silicone oil  
20 Gauge / 0.9 mm x 8 mm  
5 pcs. per box, sterile



**G-34498**

## Single-use cannula


to inject silicone oil  
23 Gauge / 0.6 mm x 8 mm  
5 pcs. per box, sterile



**G-28766**

## Single-use oil injection system

to inject silicone oil pneumatically,  
with protective cover for glass  
syringe, pressure tube fits  
MEGATRON S3 / S4 HPS, sterile



**G-31891**

## Single-use syringe

Luer-Lock, 10 ml, sterile



## Adapter to connect G-28766 silicone oil injection systems

**G-28791** for Geuder Megatron,  
ALCON Accurus, AMO Gemini


**G-28792** for Örtli systems

**G-28793** for Bausch & Lomb  
Millennium

**G-28794** Luer-Lock female

**G-28795** Luer-Lock male

**G-28796** for DORC Associate



**G-32696**

## Single-use pressure tube

for injection of viscous fluid,  
Luer-Lock female / male  
10 pcs. per box, sterile




## Stopper for viscous fluid aspiration

**G-33060** with tube connection  
for disposable syringe 5 ml

**G-33065** with tube connection  
for disposable syringe 10 ml

**G-33066** with tube connection  
for disposable syringe 20 ml



**G-33032**

## Hoerauf spreading forceps

for silicone oil removal



**G-26230**

## Sclera spreading forceps

for silicone oil removal



Manufacturing of silicone oils:

**FLUORON®**

Fluoron GmbH Magirus-Deutz-Strasse 10 89077 Ulm Germany  
Phone: +49 731 205 5997 0 Fax: +49 731 205 5997 28  
info@fluoron.de www.fluoron.de

Distribution of silicone oils / Manufacturing of instruments:

**Geuder®**  
Precision made in Germany

GEUDER AG Hertzstrasse 4 69126 Heidelberg Germany  
Phone: +49 6221 3066 Fax: +49 6221 303122  
info@geuder.de www.geuder.com

GEUDER AG reserves the right to make changes to technical details in response to recent developments.  
Geuder does not assume liability for the accuracy of each individual statement. Illustrations not drawn to scale.