


Geuder®

OPHTHALMIC SURGICAL PRODUCTS



INSTRUMENTS
MAIN CATALOGUE
06 / 2018

Geuder®
Precision made in Germany

FOCUSED ON PRECISION.

GEUDER - PRECISION MADE IN GERMANY.

GEUDER AG was founded by Hans Geuder Sr. in Heidelberg 1951, since then it has become a worldwide leading manufacturer of ophthalmic surgical instruments and systems as well as provider for highly purified biomaterials. The broad portfolio offered by GEUDER AG includes a wide range of surgical equipment for the anterior and posterior segments of the eye, along with highly innovative biomaterials for retinal and cataract surgery developed and produced by our affiliated company FLUORON GmbH.

As a result of intensive contact and cooperations with leading experts in the ophthalmology world, GEUDER AG is continuously developing and

manufacturing the most innovative products that fulfill the international standards.

A numerous amount of granted patents in more than 65 years demonstrates our intensive involvement in research and development, as well as our forward-looking philosophy. Always keeping our goal to provide the highest level of precision and reliability while complying with the requirements of ophthalmologists to serve the patient's expectations and well being.

With over 3.000 ophthalmic surgical products being distributed in more than 100 countries, GEUDER AG covers a large range of instruments



required for surgical procedures in the anterior and posterior segment of the eye. Thanks to innovative production technologies our instruments offer a unique level of precision and longevity. The constant extension of our product range – particularly single-use instruments – is following the current and especially prospective regulatory requirements and statutory regulations in the medical sector.

Qualified and specialized employees in the field of service and repair ensure that the product's performance and quality is maintained over years – an important feature for the ophthalmic surgeon who is aware of quality and environment.

Our catalogue, in addition to classic standard instruments also displays the most representative current instruments as well as novelties and new developments for innovative surgical methods.

Please contact us for more detailed information about our surgical instruments, systems and biomaterials.

Certified to ISO 9001/13485

EYE SPECULA	01
RETRACTORS	02
CALIPERS · MARKERS · MEASURING INSTRUMENTS · TONOMETER · OCULOPRESSOR	03
CURETTES · SPOONS	04
FOREIGN BODY INSTRUMENTS	05
KNIVES · DIAMOND KNIVES	06
TREPHINES	07
GLAUCOMA PUNCHES · RONGEURS	08
SCISSORS	09
FORCEPS	10
NEEDLE HOLDERS	11
MANIPULATORS · CHOPPERS · FIXATION RINGS · HOOKS · SPATULAS · LENS LOOPS	12
SCLERAL DEPRESSORS · LOCALIZERS	13
PROBES · DILATORS	14
CANNULAS · I/A INSTRUMENTS	15
SINGLE USE INSTRUMENTS	16
MISCELLANEOUS	17
STERILIZATION CONTAINERS AND ADAPTERS	18
SETS	19
INDEX	IN

CHAPTERS AND SUBCHAPTERS

01.	EYE SPECULA	13
Wire	15 – 17	
Wire · Adjustable	18 – 19	
Adjustable · Locking	20 – 26	
Aspiration · Adjustable · Locking	27 – 30	
Glaucoma · Adjustable	31	
Glaucoma · Adjustable · Aspiration	32	
02.	RETRACTORS	33
Lid	35 – 36	
Lacrimal Sac	37 – 38	
Adjustable	38	
03.	CALIPERS · MARKERS · MEASURING INSTRUMENTS	
TONOMETER · OCULOPRESSOR	39	
Corneal Markers	41 – 42	
Muscle Tucker	42	
Calipers	43	
Rulers	43 – 44	
Trabeculectomy Marker · Strabismus Marker	44	
LASIK Markers · Keratometer	45	
Degree Gauges · Axis Markers	46	
Toric Markers	47	
Tonometer · Oculopressor	48	
04.	CURETTES · SPOONS	49
Chalazion	51 – 52	
Enucleation	52	
05.	FOREIGN BODY INSTRUMENTS	53
Foreign Body Spuds · Needles · Gouges	55 – 56	
06.	KNIVES · DIAMOND KNIVES	57
Diamond Knives: Paracentesis · Sideport · Universal	59	
Diamond Knives: Tri-Facet · Clear Cornea · Trapezoid · Phaco	60	
Diamond Knives: Clear Cornea · Crescent · Step Knife · LRI	61	
Diamond Knives: Micrometers · Viscocanalostomy · Retinal Surgery	62	
Keratome · Chalazion · Viscocanalostomy · Scleral	63	
Corneal Scarifier · DMEK Scraper	64	
SATO	65	
CRVO	66	

07.	TREPHINES	67
	Trephines	69 – 70
	Punch	70
	Silicone Blocks	71 – 72
08.	GLAUCOMA PUNCHES · RONGEURS	73
	Rongeurs	75
	Glaucoma Punches	76 – 77
09.	SCISSORS	79
	Corneal · Corneal Section	81 – 82
	Stitch	82
	Synechia · Trabeculum	83
	Capsulotomy	83 – 85
	Capsule	84 – 86
	Capsulotomy · Rotatable · Interchangeable	85
	Capsule (Interchangeable and Rotatable)	86
	Iris	86 – 87
	IOL	87 – 88
	Scissors with Needle Holder	88
	Tenotomy · Conjunctival · Ring Handle	89
	Tenotomy · Strabismus · Ring Handle	90
	Utility · Ring Handle	91
	Utility · Stitch · Ring Handle	92
	Utility · Probe Tips · Ring Handle	93
	Utility · Enucleation · Ring Handle	94
	Vitreoretinal 20 Gauge	95
	Vitreoretinal 23 & 25 Gauge	96
	Vitreoretinal (Interchangeable and Rotatable)	97
10.	FORCEPS	99
	Cilia · Foreign Body	101
	Conjunctiva	102
	Suturing · Fixation	103
	Suturing · Colibri Style	104
	Corneal Fixation	104 – 105
	PAUFIQUE Suturing/Tying	105
	Tying	106 – 108
	Iris suture · Fragment	108
	Jewelers · Fixation	109
	Dressing · Tissue	110

CHAPTERS AND SUBCHAPTERS

Iris Fixation	110 – 112
Capsulorhexis	113 – 117
Capsulorhexis · Micro Incision	115
Capsulorhexis · Micro Incision · Interchangeable	116
Capsulorhexis · DMEK · Descemetorhexis	117
DSAEK	118
LASIK Flap	118 – 119
Lenticule Extraction	119
Nucleus Cracker	120
Nucleus Prechoppers	120 – 121
IOL Implantation	121
IOL Implantation · Silicone Coated	122
IOL Implantation · Iris Claw	123
IOL Folders	123 – 124
IOL Holders	124
IOL Explantation	125
Superior Rectus	126
Muscle	126 – 127
Nasal · Retinal	127
Retinal (Scleral Plug, Retinal Tack)	128
Vitreoretinal (Spreading, Trocar, Foreign Body, Magnet)	129
Vitreoretinal 20 Gauge	130 – 131
Vitreoretinal 23 & 25 Gauge	132
Vitreoretinal · Interchangeable · Rotatable	133
Chalazion · Ptosis	134
Mosquito · Hemostat	135
Towel Clamps · Serrefines	136
11. NEEDLE HOLDERS	137
Needle Holders	139 – 144
Needle Holder with Scissors · Intraocular 25 Gauge	143
Blade Holders	145
12. MANIPULATORS · CHOPPERS · FIXATION RINGS	
HOOKS · SPATULAS · LENS LOOPS	147
Fixation Rings	149 – 150
Nucleus Choppers	151 – 154
Manipulators	153 – 154
Nucleus Manipulators	155
IOL Manipulators	155 – 157
Position Hooks	156 – 157
Iris Hooks	157 – 158

Nucleus Fragmentation Spatulas	158
Phaco Spatulas	159
Capsule Polishers	160
Lens Loops	161
Orbital Spatulas	162
Keratoplasty · Synechiae Spatulas	163
Cyclodialysis Spatulas	164
Iris Spatulas	165 – 166
Enucleation	167
LASIK Spatulas	167 – 168
DMEK Hooks · DSAEK Glide	168
Femtosecond Cataract	169
Femtosecond LASIK	169 – 172
Strabismus Muscle Tucker	172
Strabismus Hooks	173 – 174
Retinal Detachment Hook	174
Membrane Spatulas · Membrane Peelers	175
13. SCLERAL DEPRESSORS · LOCALIZERS	177
Scleral Depressors	179 – 180
Localizers	180
14. PROBES · DILATORS	181
Lacrimal Probes	183 – 186
Corneal Implant Probes	184
Trabecular Probes · U-Shaped Probe	186
15. CANNULAS · I/A INSTRUMENTS	187
Lacrimal	189 – 190
Lacrimal Probes	190
Retrobulbar · Eye Douche · A.C. Maintainers · Cystotome · Enzyme	191
Hydrodissection · Cyclodialysis	192
Irrigation Vectis · Healon · Trabeculotomy · Sub-Tenon	193
Irrigation	194
Irrigation/Aspiration	194 – 195
Capsule Polishers	195 – 196
LASIK	196 – 197
Irrigation Cannula Handpieces	197
I/A Handpieces · Bimanual	198
Bimanual Irrigation/Aspiration	199 – 201
Irrigation Choppers · Irrigation Hook · Irrigation Scraper	202
Coaxial Irrigation/Aspiration	203 – 205

CHAPTERS AND SUBCHAPTERS

DMEK Irrigation.....	203
Coaxial Irrigation/Aspiration.....	204
Retinal · Backflush Handpieces.....	206
Retinal · Backflush Cannulas	207
Retinal · Backflush Tips	208
Retinal · Infusion/Injection Cannulas	209
Retinal · Infusion · Aspiration · PFCL.....	210
Retinal.....	211
Miscellaneous.....	212
16. SINGLE USE INSTRUMENTS.....	213
Manual Trehines · Cutting Block · Non-Vacuum Punches.....	215
Vacuum Trehines & Punches	216
Surgical Knives	217
Scissors.....	218
Forceps.....	219
Endo Probes	220
Trocars.....	221
Fibre Optics	221 – 222
Vitreous Cutters · Infusion Tube	222
Speculum · DMEK Cartridge · Iris Retractors · Lacrimal Intubation	223
Intubation Sets · Cannulas (Lacrimal, Retrobulbar, Cystotome)	224
Cannulas (A.C. Maintainer, Rycroft, Hydrodissection)	225
Cannulas (Irrigation Vectis, Capsule Polisher, Simcoe, Viscocanalostomy, Irrigation)	226
Bimanual I/A Handpiece · Backflush Handpiece · Backflush Cannula · Silicone Brush	227
Cannulas (Backflush with Silicone Tip, PFCL, Silicone Oil, Infusion)	228
Cannulas (Injection, Viscous Fluids, Gas)	229
Oil Injection	230
Stopcock · Scalpel Blades · Absorbent Swabs · Marking Pen.....	231
Silicone Tires.....	232
Silicone Bands	233
Silicone Strips	234
Silicone Sleeves.....	234 – 235
Silicone Wedge.....	235
Silicone Sponge Implants	235 – 238
Bonn Injection Set	238
17. MISCELLANEOUS	239
Injectors (for MORCHER Capsular Tension Ring and Pupil Dilators)	241
Lacrimal Instruments · Rongeur · Nasal Speculum · Hammer · Elevatorium · Raspatorium	242
Bone Chisels · LASEK Funnel · Fascia Needle	243

Diamond Burr · Vitrectomy Lenses	244
Eye Shield · Eye Spheres.....	245
Battery Cautery · Intraocular Magnet.....	246
18. STERILIZATION CONTAINERS AND ADAPTERS.....	247
MICROSTOP®	249 – 251
Steel Cases	252 – 255
Trays & Racks for Steel Cases	256
Silicone Strips & Silicone Mats.....	257
Seals · Small Cases.....	258
STERISAFE® DURO	259 – 262
STERISAFE® DURO · Melag.....	263
Adapters	264 – 266
19. SETS	S1
ECCE – Extra Capsular Cataract Extraction	S2
Phaco / Cataract	S3
Phaco Emergency / MSICS – Manual Small Incision Cataract Surgery.....	S4
MICS – Micro Incision Cataract Surgery / FLACS – Femtosecond Laser-assisted Cataract Surgery	S5
Snyder-Osher IOL Explantation · Toric IOL · Hattenbach Iris suture	S6
LASEK · LASIK.....	S7
Femto LASIK · ReLex SMILE.....	S8
Keratoplasty	S9
DSEK / DSAEK · DALK.....	S10
DMEK.....	S11
Glaucoma.....	S12
Emergency	S13
Enucleation · Pterygium	S14
Office set (Optometrist) · Chalazion	S15
Lid / Ptosis	S16
DCR / Lacrimal.....	S17
Strabismus / Muscle.....	S18
Retinal Detachment	S19
Vitrectomy 20G	S20 – S21
Vitrectomy 23G	S22 – S23
Vitrectomy 25G	S24
20. INDEX	I1
Numerical Index.....	I2 – I9
Keyword Index	I10 – I17
Name Index	I18 – I20

© GEUDER AG

Although every effort has been made to ensure the accuracy of the content presented in this catalogue, GEUDER AG cannot be held liable for any errors that may occur.

Due to the constant technological progress all product specifications and illustrations are subject to change without notice.

No part of this catalogue may be reproduced, edited, copied or disseminated in any form (analog or digital) without prior written consent of GEUDER AG.

All rights reserved by GEUDER AG, Heidelberg, Germany.
Any violation of the copyright will be prosecuted.

Edition 06/2018

© GEUDER AG

Manufacturer of ophthalmic surgical instruments and systems


Some Geuder products may not be available for sale in certain countries where the products have not yet received investigational, clinical or marketing approval. Please contact your local distributor or GEUDER AG for more information.

GEUDER AG
Hertzstr. 4
69126 Heidelberg
Germany

Phone: +49 6221 306-6
Fax: +49 6221 303122
info@geuder.de
www.geuder.de

„NanOEdge®“, „UNO Colorline®“, „Densiron®“, „Megatron®“, „endoTRON®“ and „MACH2“ are registered trademarks of Geuder AG.
„EasyGas®“ is a registered trademark of Fluoron GmbH.
„Sterisafe®“ is a registered trademark of SAVUNA GmbH & Co. KG.
„microStop®“ is a registered trademark of Gebrüder Martin GmbH & Co. KG.
„ReLEX® Smile“ and „SMILE“ are registered trademarks of Carl Zeiss Meditec AG.
„IRIDEX®“ is a registered trademark of IRIDEX Corporation.
„ACCURUS®“ is a registered trademark of Alcon Laboratories, Inc.
„Pentasy®“ is a registered trademark of Fritz RUCK Ophthalmologische Systeme GmbH.
„Millennium®“ is a registered trademark of MBI Millennium Biomedical, Inc.

EYE SPECULA




WIRE	15
WIRE AND ADJUSTABLE	18
ADJUSTABLE / LOCKING	20
WITH ASPIRATION	27
ADJUSTABLE / BLADE ELEVATION	31

G-17020*BARRAQUER***EYE SPECULUM**

fenestrated

14 mm blades

overall length 38 mm



01

WIRE


G-17025*BARRAQUER***EYE SPECULUM**

for infants

fenestrated

10 mm blades

overall length 40 mm


**G-17023***BARRAQUER***EYE SPECULUM**

for premature infants

fenestrated


7 mm blades

overall length 34 mm

**G-17026***BARRAQUER***EYE SPECULUM**

15 mm solid blades

overall length 40 mm



EYE SPECULA


01

WIRE

G-32136

BARRAQUER
EYE SPECULUM


fenestrated
14 mm round blades
overall length 38 mm



G-32137

BARRAQUER
EYE SPECULUM


for infants
fenestrated
10 mm round blades
overall length 38 mm



G-32148

BARRAQUER
EYE SPECULUM


fenestrated
14.5 mm open blades
overall length 40 mm



G-32149

BARRAQUER FRANKFURT MODEL
EYE SPECULUM

fenestrated, large
14.5 mm blades
overall length 40 mm




G-32133*KRATZ-BARRAQUER***EYE SPECULUM**

fenestrated

15 mm round wing blades

overall length 40 mm

to lift and retract the lids away from the globe
 eliminating external pressure



01

WIRE

*OHNO***EYE SPECULUM**

fenestrated

14 mm round blades


overall length 45 mm

middle bar controls redundant lid tissue

G-32126

for the right eye

for the left eye


**G-17004***H.-R. KOCH***TEMPORAL APPROACH EYE SPECULUM**

for nasal placement

fenestrated

14 mm open blades

overall length 40 mm


**G-32129****TEMPORAL APPROACH EYE SPECULUM**

for nasal placement

fenestrated

15 mm round blades

overall length 40 mm



EYE SPECULA

01

WIRE · ADJUSTABLE


BARRAQUER-OOSTERHUIS

EYE SPECULUM

fenestrated, adjustable

14 mm blades, overall length 40 mm

10 mm blades, overall length 30 mm, for infants



G-17010

G-17011

BARRAQUER-OOSTERHUIS

EYE SPECULUM


solid blades, adjustable

G-17015

G-17016

15 mm blades, overall length 40 mm

11 mm blades, overall length 30 mm, for infants



BARRAQUER-OOSTERHUIS

EYE SPECULUM

fenestrated round blades


adjustable

G-32138

G-32139

14 mm blades, overall length 40 mm

11 mm blades, overall length 30 mm, for infants



G-17017


BARRAQUER-OOSTERHUIS-KRATZ

EYE SPECULUM

adjustable, fenestrated


15 mm open blades

overall length 40 mm



G-32119*KRATZ-BARRAQUER***EYE SPECULUM**

adjustable, fenestrated
 15 mm round wing blades
 overall length 40 mm



01

G-17027*THOMAS NEUHANN***EYE SPECULUM**


adjustable
 15 mm solid blades
 overall length 40 mm

**G-17030***SAUER***EYE SPECULUM**

12 mm solid blades
 overall length 30 mm

**G-32134***TAKEDA***EYE SPECULUM**


for infants
 9 mm fenestrated blades
 with finger grips and
 angled interlacing valves
 overall length 28 mm




WIRE · ADJUSTABLE

G-15945**MELLINGER
EYE SPECULUM**


large
adjustable, self retaining
16.5 mm fenestrated blades
overall length 65 mm

**G-32130****MELLINGER
EYE SPECULUM**


small, for infants
adjustable, self retaining
12.5 mm fenestrated blades
overall length 55 mm

**G-15946****MELLINGER-BEARD
EYE SPECULUM**

large
adjustable, self retaining
15.5 mm solid blades
overall length 65 mm

**G-32135****MELLINGER-BEARD
EYE SPECULUM**

small, for infants
adjustable, self retaining
11.5 mm solid blades
overall length 55 mm




*PRETORI***EYE SPECULUM**

13.5 mm fenestrated blades
adjustable, self retaining
overall length 63 mm

G-16050


G-16055



G-15950

*WEISS***EYE SPECULUM**


fenestrated
adjustable with locking screw
17.5 mm blades
overall length 85 mm



G-15951

*WEISS***EYE SPECULUM**


fenestrated
adjustable with locking screw
17.5 mm open blades
overall length 85 mm



G-15955

*WEISS***EYE SPECULUM**

small
adjustable with locking screw
15 mm fenestrated blades
overall length 55 mm



EYE SPECULA

G-15956*WEISS***EYE SPECULUM**

small

adjustable with locking screw

15 mm open blades

overall length 55 mm

**G-15970***WEISS***EYE SPECULUM**

adjustable with locking screw

16.5 mm solid blades

overall length 85 mm

**G-15975***WEISS***EYE SPECULUM**

small

adjustable with locking screw

14 mm solid blades

overall length 55 mm


**G-37130***BARTL***EYE SPECULUM**

adjustable

with silicone band to secure vitrectomy lens

15 mm solid blades


overall length 39 mm



G-16010

*BANGERTER***EYE SPECULUM**

adjustable, with locking screw
 15 mm solid blades
 overall length 52 mm



G-16020

*CASTROVIEJO***EYE SPECULUM**

fenestrated, adjustable
 overall length 96 mm
 14 mm blades, small
 16 mm blades, medium
 18 mm blades, large




G-16025

G-16030

G-15990

*COOK***EYE SPECULUM**

large
 16 mm solid blades
 adjustable, with locking screw
 overall length 50 mm



G-15995

*COOK***EYE SPECULUM**

medium
 11 mm solid blades
 adjustable, with locking screw
 overall length 40 mm



EYE SPECULA

G-16000

COOK
EYE SPECULUM

small
8 mm solid blades
adjustable, with locking screw
overall length 35 mm

**G-32128**


COOK
EYE SPECULUM

for infants
5 mm solid blades
adjustable, with locking screw
overall length 30 mm

**G-17000**


MURDOCK
EYE SPECULUM

left
14 mm fenestrated blades
adjustable, self retaining
overall length 50 mm

**G-17001**

MURDOCK
EYE SPECULUM

right
14 mm fenestrated blades
adjustable, self retaining
overall length 50 mm



MURDOCK

EYE SPECULUM

small


12 mm fenestrated blades
adjustable, self retaining
overall length 45 mm

G-17002


G-17003

left


right




G-15960

*LIEBERMANN-TENNANT***EYE SPECULUM**adjustable, fenestrated
for temporal placement
15 mm open blades
overall length 80 mm

G-15965

*KERSHNER***EYE SPECULUM**adjustable, reversible
for nasal and temporal placement
15 mm fenestrated blades
overall length 75 mm


G-15966

*KERSHNER***EYE SPECULUM**adjustable, reversible
for nasal and temporal placement
15 mm solid blades
overall length 75 mm

G-33965

DICK
FEMTOCATARACT EYE SPECULUM

adjustable
 18 mm flattened fenestrated blades with
 middle bar to control redundant lid tissue
 overall length 80 mm

**G-33971**

SOLOMON
FEMTOCATARACT EYE SPECULUM

14.5 mm flattened fenestrated blades with
 middle bar to control redundant lid tissue
 double knurled screw for fast adjustment and
 fine tuning
 overall length 90 mm

**G-33919**


SPARE TUBING SET
 for eye speculum with aspiration

1/2

**G-33920**

KNORZ
LASIK EYE SPECULUM

adjustable, with aspiration
 16 mm fenestrated blades
 to remove excess fluid from the surgical field
 supplied with silicone tubing and adaptor





NODA

EYE SPECULUM

adjustable, with aspiration
to remove excess fluid from the surgical field
supplied with silicone tubing and adaptor

**G-17019**

15 mm open blades

G-17021

12 mm open blades, for infants

**G-17005**

PAULMANN


EYE SPECULUM

adjustable, with aspiration
for PPV, with extra aspiration bar
16 mm fenestrated blades
to remove excess fluid from the surgical field
supplied with silicone tubing and adaptor

G-17006**CLAMP FOR INFUSION TUBE**

for PAULMANN eye speculum (G-17005)






G-15961

LIEBERMANN-TENNANT**EYE SPECULUM**adjustable, with aspiration
for temporal placement

15 mm open blades

overall length 80 mm

to remove excess fluid from the surgical field
supplied with silicone tubing and adaptor

G-17024

KOGA / NODA**EYE SPECULUM**

with 5 mm aspiration slits


adjustable

12 mm open blades

to remove excess fluid from the surgical field

supplied with silicone tubing and adaptor





G-17031

**NODA
EYE SPECULUM**


adjustable, with aspiration
open blades
2 x 10 aspiration ports
supplied with silicone tubing and adaptor




G-15962

**KAMMANN
EYE SPECULUM**

adjustable, with aspiration
self retaining
11 mm fenestrated blades
to remove excess fluid from the surgical field
supplied with silicone tubing and adaptor






G-33914

ZIRM

LASIK EYE SPECULUM

adjustable, with aspiration
16 mm fenestrated blades
to remove excess fluid from the surgical field
supplied with silicone tubing and adaptor




G-33925

MAUS

LASIK EYE SPECULUM

adjustable, with middle aspiration bar
16 mm fenestrated blades
to remove excess fluid from the surgical field
supplied with silicone tubing and adaptor

**EYE SPECULUM**


fenestrated blades
 adjustable mechanism regulates opening and blade elevation
 optimal for glaucoma surgery

- G-15920** large, 16 mm blades
G-15925 small, 12 mm blades

**SCHOTT****EYE SPECULUM**

solid blades
 adjustable mechanism regulates opening and blade elevation
 optimal for glaucoma surgery

- G-15921** large, 15 mm blades
G-15926 small, 11 mm blades



G-15928

**SCHOTT
EYE SPECULUM**

with aspiration
13 mm fenestrated blades
for infants
adjustable mechanism regulates
opening and blade elevation
optimal for glaucoma surgery

02

RETRACTORS





G-16080 *DESMARRES*
LID RETRCTOR
for infants
7 mm wide
overall length 133 mm



G-32120 *DESMARRES*
LID RETRCTOR
fig. 0
9 mm wide
overall length 138 mm



G-16060 *DESMARRES*
LID RETRCTOR
fig. 1
12 mm wide
overall length 133 mm



G-16065 *DESMARRES*
LID RETRCTOR
fig. 2
14 mm wide
overall length 135 mm

RETRACTORS

02

LID

G-16070

DESMARRES
LID RETRCTOR

fig. 3
16 mm wide
overall length 138 mm



G-16075

DESMARRES
LID RETRCTOR

fig. 4
17 mm wide
overall length 138 mm



G-32115

GIESSEN MODEL
LID RETRCTOR

14 mm wide, fenestrated
overall length 133 mm



G-17040

JAEGER
LID PLATE

stainless steel
20 mm and 23 mm concave ends
overall length 110 mm





- G-15542** AXENFELD
LACRIMAL SAC RETRCTOR
 4 prongs, blunt
 9 mm wide
 overall length 140 mm

02

LACRIMAL SAC



- AXENFELD
LACRIMAL SAC RETRCTOR
 3 prongs, 7 mm wide
 overall length 140 mm

- G-15540** blunt
G-15545 sharp



- AXENFELD
LACRIMAL SAC RETRCTOR
 2 prongs, 4 mm wide
 overall length 140 mm

- G-15550** blunt
G-15555 sharp



- G-15560** ROLLET
LACRIMAL SAC RETRCTOR
 4 prongs, sharp
 6 mm wide
 overall length 140 mm

**GUTHRIE****RETRACTOR**

double skin hook, sharp
overall length 125 mm
small
large


G-15350
G-15355



G-17085 *AMERICAN MODEL*
LACRIMAL SAC RETRATOR
 4 x 4 prongs, blunt
 cross-action
 adjustable, self retaining
 overall length 100 mm



G-17070 *MUELLER*
LACRIMAL SAC RETRATOR
 3 x 3 prongs
 adjustable
 overall length 70 mm



G-17080 *GOLDSTEIN*
LACRIMAL SAC RETRATOR
 3 x 3 prongs
 adjustable
 overall length 30 mm

CALIPERS, MARKERS
MEASURING INSTRUMENTS
TONOMETER, OCULOPRESSOR



CORNEAL MARKERS	41
CORNEAL MARKERS · MUSCLE TUCKER	42
CALIPERS · RULERS	43
RULERS · TRABECULECTOMY MARKER · STRABISMUS MARKER	44
LASIK MARKERS · KERATOMETER	45
DEGREE GAUGES · AXIS MARKERS	46
TORIC MARKERS	47
TONOMETER · OCULOPRESSOR	48

CALIPERS, MARKERS, MEASURING INSTRUMENTS,
TONOMETER, OCULOPRESSOR

03



**THOMAS NEUHANN
CORNEAL MARKER**
with cross hairs
inside diameter 3.0 mm
outside diameter 12.0 mm
overall length 105 mm

G-33725 4 blades
G-33726 6 blades
G-33727 8 blades




03

CORNEAL MARKERS



G-32222 **THOMAS NEUHANN
CORNEAL MARKER**
for astigmatic keratotomy
10° scale and cross-hairs
6.0 mm inside diameter
8.0 mm outside diameter
overall length 105 mm



**HOFFER
CORNEAL MARKER**
for optic zone
with cross hairs
overall length 103 mm
G-33740 inside diameter 7.0 mm
G-33741 inside diameter 7.5 mm
G-33742 inside diameter 8.0 mm



G-38600 **THOMAS NEUHANN
CORNEAL MARKER**
for DSAEK/DSEK/DMEK procedures
combination of 8.0 mm / 9.0 mm diameters
to produce a guide for Descemetorhexis
overall length 123 mm



CALIPERS, MARKERS, MEASURING INSTRUMENTS, TONOMETER, OCULOPRESSOR

03

CORNEAL MARKERS · MUSCLE TUCKER

LINDSTROM

CORNEAL MARKER

to mark arcuate incisions
for surgical correction of astigmatism
overall length 100 mm

G-33780

G-33781



DUNCKER

SUTURE MARKER

8-pointed star and cross hairs
to center the trephine and perform
double running suture
overall length 102 mm
for 7.0 mm transplant
for 8.0 mm transplant

G-32228

G-32229



G-32224

THOMAS NEUHANN

SUTURE MARKER

8 blades
for keratoplasty
fits with KRUMEICH GTS
overall length 102 mm



G-19137

CUEPPERS

STRABISMUS MUSCLE TUCKER

15 mm scale
serrated J-shaped hook to
prevent the muscle from slipping



CASTROVIEJO

CALIPER

scale on both sides
with 1 mm increments

G-19130 15 mm scale, overall length 78 mm
G-19135 20 mm scale, overall length 85 mm



03


CALIPERS · RULERS

G-19138

OSHER

INCISION MEASURING CALIPER

1 to 10 mm
scale on both sides
for internal measurement of cataract incisions
overall length 75 mm



G-37400

BRAUNSTEIN

FIXED CALIPER

for marking of sclerotomy site
at 3.5 mm posterior to the surgical limbus
overall length 45 mm



BRAUNSTEIN

FIXED CALIPER

for marking of sclerotomy site posterior
to the surgical limbus
overall length 80 mm

G-37410

3.0 mm / 3.5 mm

G-37412

3.0 mm / 4.0 mm

G-37414

3.5 mm / 4.0 mm

G-37418

4.0 mm / 5.0 mm



G-18650

RULER

150 mm
0.5 mm and 1.0 mm increments

CALIPERS, MARKERS, MEASURING INSTRUMENTS, TONOMETER, OCULOPRESSOR

03

RULERS · TRABECULECTOMY MARKER · STRABISMUS MARKER

G-18645 *GIESSEN MODEL*

RULER

vaulted

1 to 30 mm

overall length 60 mm



G-32225 *VOGT*
SCLERAL RULER

vaulted

1 to 20 mm

overall length 126 mm



G-32004 *ENGELS*
MEASURING SPATULA

1 to 15 mm

angled 45°

overall length 88 mm



G-32210 *FRONIMOPOLUS*
TRABECULECTOMY MARKER

5 x 5 mm

U-shaped marking pattern

overall length 70 mm



G-32209 *TOBIAS NEUHANN*
SCLERAL FLAP MARKER

5.0 mm and 4.0 mm

U-shaped marking patterns

swivel handle allows marking from temporal
or superior approach



G-15797 *FINK*
STRABISMUS MARKER

1 mm criss-cross marking tips

for accurate 9 mm recession of the
inferior oblique muscle


overall length 87 mm





G-33785

GRENE
VISUAL AXIS MARKER
0.5 mm diameter
concave tip
overall length 103 mm




03




G-33745

MENDEZ
LASIK CORNEAL MARKER
for lamellar flap identification
and realignment
overall length 105 mm




G-33746

CHAYET
LASIK CORNEAL MARKER
for lamellar flap identification
and realignment
overall length 105 mm



G-33747

LU
LASIK CORNEAL MARKER
for lamellar flap identification
and realignment
overall length 105 mm



G-32318

MALONEY
KERATOMETER
for qualitative measurement of astigmatism
overall length 20 mm



CALIPERS, MARKERS, MEASURING INSTRUMENTS, TONOMETER, OCULOPRESSOR

G-33760

MENDEZ


DEGREE GAUGE

measures from 0° to 180°

5° increments

12 mm inside diameter

overall length 110 mm



03

G-33762

DEGREE GAUGE

for toric IOLs

with 5° increments and

rotating alignment rings

12 mm inside diameter

overall length 110 mm



G-33763

GERTEN

HORIZONTAL AXIS MARKER

for precise marking of the 0° axis

during toric IOL implantation

or refractive laser surgery

overall length 105 mm



G-33774

SCOTT

CORNEAL PENDULUM MARKER

for temporary placement of

0°, 90°, and 180° pre-op toric marks

self-leveling pendulum assures precise

horizontal axis

overall length 110 mm



G-33777

TOMARK

ADJUSTABLE CORNEAL MARKER

to mark the axis of astigmatism during toric IOL placement

0° to 180° degree gauge

5° increments

overall length 32 mm



for use in combination with slit lamp or G-33776

03

G-33776

TOMARK PENDULUM INSTRUMENT

self-leveling handheld pendulum assures precise

horizontal axis

overall length 140 mm

to retain TOMARK corneal marker G-33777



G-33778


RETAINER FOR PENDULUM INSTRUMENT

for G-33776 TOMARK

overall length 180 mm

1/2





G-19230 SCHIOETZ
TONOMETER
with two weights of 7.5 g and 10 g

G-15392 HOMBURG MODEL
OCULOPRESSOR

G-15393 SILICONE BALL
with band, for Oculopressor
HOMBURG MODEL G-15392



04

CURETTES SPOONS





G-15900 *MEYHOFER*
CHALAZION CURETTE
1.0 mm diameter cup
overall length 129 mm

04



G-15870 *MEYHOFER*
CHALAZION CURETTE
2.0 mm diameter cup
overall length 129 mm

CHALAZION



G-15875 *MEYHOFER*
CHALAZION CURETTE
2.5 mm diameter cup
overall length 129 mm



G-15880 *MEYHOFER*
CHALAZION CURETTE
3.0 mm diameter cup
overall length 129 mm



G-15885 *MEYHOFER*
CHALAZION CURETTE
4.0 mm diameter cup
overall length 129 mm

G-15840

HEBRA
CHALAZION CURETTE
 small
 2.5 mm x 4.0 mm oval scoop
 overall length 133 mm

**04****G-15845**

HEBRA
CHALAZION CURETTE
 medium
 3 mm x 5 mm oval scoop
 overall length 133 mm

**G-15910**

AXENFELD
CHALAZION CURETTE
 double ended
 4.5 x 2.5 mm and 6.5 x 3.0 mm oval curettes
 overall length 144 mm

**G-31545**

BUNGE
EXENTERATION CURETTE
 medium
 10 mm spoon
 overall length 142 mm

**G-31555**

BUNGE
EXENTERATION CURETTE
 large
 12 mm spoon
 overall length 144 mm



FOREIGN BODY INSTRUMENTS





- G-14550 FOREIGN BODY INSTRUMENT**
double-ended with sliding handle
2.2 mm wide curved needle
1 mm wide *PLANGE* spud
overall length 95 mm



- G-14540 FOREIGN BODY INSTRUMENT**
double ended with sliding handle
2.2 mm wide curved needle
1.4 mm wide gouge
overall length 95 mm



- G-14530 FOREIGN BODY INSTRUMENT**
double ended with sliding handle
2.2 mm wide curved needle
1.4 mm wide spud
overall length 95 mm

05

**G-14490 FOREIGN BODY NEEDLE**

straight
4.5 mm long triangular tip
2.4 mm wide
overall length 120 mm

**G-14500 FOREIGN BODY NEEDLE**

curved
4.5 mm long triangular tip
2.4 mm wide
overall length 120 mm

**G-14570 PLANGE**

FOREIGN BODY SPUD
curved and rounded
1 mm wide spud
overall length 125 mm

**G-14480 FOREIGN BODY GOUGE**

straight
1.4 mm wide gouge
overall length 116 mm

KNIVES
DIAMOND KNIVES

PARACENTESIS · SIDEPORT · UNIVERSAL.....	59
TRIFACET · CLEAR CORNEA.....	60
CRESCENT · STEP TRIFACET · LRI.....	61
MICROMETER · VISCOCANALOSTOMY · RETINAL SURGERY.....	62
KERATOME · CHALAZION · SCLERAL.....	63
CORNEAL SCARIFIER · DMEK SCRAPER.....	64
SATO.....	65
CRVO.....	66



G-38475

DIACLEAN DIAMOND KNIFE

for paracentesis



1 mm wide spear

titanium handle

width 1.0 mm

length 4.0 mm

thickness 0.2 mm



06

G-38344

DIACLEAN DIAMOND KNIFE



1 mm wide, 45° diamond blade

titanium handle

width 1.0 mm

length 3.0 mm

thickness 0.2 mm



G-38341

DIACLEAN DIAMOND KNIFE



60° spear diamond blade

titanium handle

width 1.4 mm

length 3.0 mm

thickness 0.2 mm



G-38343

DIACLEAN DIAMOND KNIFE



45° spear diamond blade

titanium handle

width 1.4 mm

length 3.0 mm



thickness 0.2 mm



G-38342

DIACLEAN TRI-FACET DIAMOND KNIFE

1 mm wide tri-facet diamond blade
 titanium handle
 width 1.0 mm
 length 3.0 mm
 thickness 0.2mm





06

G-38473

*PHAM***DIACLEAN TRI-FACET DIAMOND KNIFE**



titanium handle
 width 1.5 mm
 length 3.0 mm
 thickness 0.2 mm



G-38448

**DIACLEAN CLEAR-CORNEA-INCISION
DIAMOND KNIFE, 1.5/2.0 mm**



MICS trapezoid blade, 35° angled
 4 cutting edges
 titanium handle
 width 1.5 / 2.0 mm
 length 5.0 mm
 thickness 0.2 mm



G-38444



**DIACLEAN CLEAR-CORNEA-INCISION
DIAMOND KNIFE, 2.2 MM**

4 cutting edges, 35° angled
 titanium handle
 width 2.2 mm
 length 5.0 mm
 thickness 0.2 mm






G-38453 DIACLEAN CLEAR-CORNEA-INCISION DIAMOND KNIFE, 2.85 MM



4 cutting edges, 35° angled
 titanium handle
 width 2.85 mm
 length 5.0 mm
 thickness 0.2 mm


G-38466 DIACLEAN DIAMOND CRESCENT KNIFE, 2.0 MM



round spatula shape, 40° angled
 titanium handle
 width 2.0 mm
 length 4.0 mm
 thickness 0.2 mm


G-31470 TRI-FACET DIAMOND STEP KNIFE

single footplate, titanium handle
 stepwise adjustable
 in sterilizing case
 width 1.0 mm
 length 3.0 mm
 thickness 0.2 mm


G-31468 LIMBAL RELAXING INCISION DIAMOND KNIFE


LRI knife with 0.6 mm precalibrated depth
 single footplate, titanium handle
 width 1.0 mm
 length 3.0 mm
 thickness 0.2 mm



G-31480

DIAMOND MICROMETER KNIFE


tri-facet for keratotomy
 vertical and 15° cutting edges with
 0.2 mm flat base cutting edge
 double footplate, titanium handle
 width 0.6 mm
 length 3.0 mm
 thickness 0.1 mm



G-19345

DIAMOND MICROMETER KNIFE


for keratotomy
 single 45° cutting edge
 double footplate, titanium handle
 width 1.0 mm
 length 3.0 mm
 thickness 0.2 mm



G-31487

DIAMOND KNIFE


for viscocanalostomy
 0.5 mm wide, 70° diamond blade
 width 0.5 mm
 length 4.0 mm
 thickness 0.2 mm



G-31477

RETINA DIAMOND KNIFE

45° cutting edge
 20 gauge / 0.9 mm
 titanium handle
 width 0.9 mm
 length 3.0 mm
 thickness 0.2 mm



**G-14420 HAMBURG MODEL****KERATOME**


angled at 14 mm from tip
 10 mm cutting surface
 1.1 mm wide
 overall length 120 mm

**G-14000 CHALAZION KNIFE**

straight, delicate
 10 mm cutting surface
 1.8 mm wide
 overall length 120 mm

**G-14105 MICRO ROUND KNIFE**

for viscocanalostomy
 0.8 mm diameter disc-shaped blade
 overall length 112 mm

**G-14180 BONN MODEL****SCLERAL KNIFE**

7 mm blade, rounded tip
 overall length 126 mm

**G-14185 BONN MODEL****SCLERAL KNIFE**

4 mm blade, rounded tip
 overall length 126 mm

*KUHNT***CORNEAL SCARIFIER**

rounded tip

overall length 123 mm

G-14010**G-14020**

fig. 1, 5.5 mm x 2.5 mm blade

fig. 2, 4.5 mm x 2.0 mm blade

**G-38621***SZURMAN***DMEK PREPARATION SCRAPER**

4.0 mm lower blade, blunt

1.5 mm rounded tip, sharp

for removal of Descemet's membrane from graft

overall length 117 mm

**G-14200***PAUFIQUE***KERATOPLASTY KNIFE**

2 mm long cutting edge

1.7 mm wide semi-circular tip

overall length 105 mm

**G-33670***LOHMANN***LASEK EPI PEELER**

"Hockey" knife with

5 mm bottom blade and 2 mm front blade

overall length 123 mm

**G-33022****STILETTO**

45° lancet

overall length 53 mm

for instruments with

G-33025

23 gauge / 0.6 mm

G-33026

20 gauge / 0.9 mm

19 gauge / 1.1 mm



**G-14165****SATO KNIFE**

1.5 mm cutting edge
pointed tip
overall length 123 mm

**G-16200****DARDENNE****SATO KNIFE**

1.2 mm cutting edge
pointed tip
overall length 123 mm



06

SATO

**G-37555****FOERSTER****SATO KNIFE**

conical shaft
front: 30 gauge / 0.3 mm
rear: 20 gauge / 0.9 mm
overall length 125 mm

**G-14162****HANSEN****SATO KNIFE**

double edge,
for arteriovenous dissection
20 gauge / 0.9 mm
overall length 130 mm





G-12200

CRVO KNIFE

for Radial Opticus Neurotomy (RON)
with Central Retinal Vein Occlusion (CRVO),
20 gauge / 0.9 mm
overall length 120 mm



Geuder




G-12220

DITHMAR

CRVO KNIFE

with preset incision depths at 1.5 mm,
for radial opticus neurotomy (RON)
with Central Retinal Vein Occlusion (CRVO),
20 gauge / 0.9 mm
overall length 129 mm



Geuder

06

CRVO

TREPHINES



BEST
TREPHINE

overall length 58 mm

G-31815	1.0 mm	G-31835	1.5 mm
G-31820	1.2 mm	G-31840	1.6 mm
G-31825	1.3 mm	G-31845	1.7 mm
G-31830	1.4 mm	G-31850	1.8 mm


SCHMIDT-MUMM
DUCT TREPHINE

overall length 38 mm



G-20610	1.0 mm	G-20613	1.3 mm
G-20611	1.1 mm	G-20614	1.4 mm
G-20612	1.2 mm	G-20615	1.5 mm

G-20616	1.6 mm
G-20617	1.7 mm
G-20618	1.8 mm

07

HAMBURG MODEL
KERATOPLASTY TREPHINE

overall length 40 mm



G-20050	5.0 mm	G-20060	6.0 mm
G-20054	5.4 mm	G-20065	6.5 mm
G-20055	5.5 mm	G-20070	7.0 mm

G-20080	8.0 mm
G-20090	9.0 mm
G-20100	10.0 mm

BEST
KERATOPLASTY TREPHINE

overall length 75 mm



G-20450	5.0 mm	G-20470	7.0 mm
G-20455	5.5 mm	G-20480	8.0 mm
G-20460	6.0 mm	G-20490	9.0 mm
G-20465	6.5 mm	G-20500	10.0 mm

FRANCESCHETTI
KERATOPLASTY TREPHINE

overall length 75 mm



G-17250	5.0 mm	G-17270	7.0 mm
G-17251	5.1 mm	G-17271	7.1 mm
G-17260	6.0 mm	G-17280	8.0 mm
G-17261	6.1 mm	G-17281	8.1 mm

G-17290	9.0 mm
G-17291	9.1 mm
G-17300	10.0 mm


MANY TREPHINES ALSO AVAILABLE IN MINUS, INTERMEDIATE AND PLUS SIZES

G-31870

THOMAS NEUHANN

CORNEAL PUNCH

for keratoplasty trephines with thread from
G-33870 to G-33899, with silicone block for
donor cornea

**KERATOPLASTY TREPHINE**

with thread, for *NEUHANN*
corneal punch G-31870
overall length 40 mm



G-33870	7.0 mm
G-33871	7.1 mm
G-33872	7.2 mm
G-33873	7.3 mm
G-33874	7.4 mm
G-33875	7.5 mm
G-33876	7.6 mm
G-33877	7.7 mm
G-33878	7.8 mm
G-33879	7.9 mm

G-33880	8.0 mm
G-33881	8.1 mm
G-33882	8.2 mm
G-33883	8.3 mm
G-33884	8.4 mm
G-33885	8.5 mm
G-33886	8.6 mm
G-33887	8.7 mm
G-33888	8.8 mm
G-33889	8.9 mm

G-33890	9.0 mm
G-33891	9.1 mm
G-33892	9.2 mm
G-33893	9.3 mm
G-33894	9.4 mm
G-33895	9.5 mm
G-33896	9.6 mm
G-33897	9.7 mm
G-33898	9.8 mm
G-33899	9.9 mm

G-32671

ÖRK LINZ MODEL
SILICONE BLOCK FOR DMEK
for donor cornea
white-colored

**G-32672**

SILICONE BLOCK
for donor cornea
blue-colored



07

SILICONE BLOCKS

G-32674

SILICONE BLOCK WITH CENTER POINT
for donor cornea

**G-38620**

SZURMAN
DMEK GRAFT RETAINER
for preparation of DMEK graft



G-21175

METAL STAND

with removable silicone insert G-21176
for aspiration of donor cornea



07


G-21176

SILICONE INSERT

for metal stand G-21175



GLAUCOMA PUNCHES RONGEURS




G-18210

CITELLI
BONE NIBBLING RONGEUR
 3 mm x 4.5 mm bite
 overall length 150 mm

**G-18220**

BEYER
BONE NIBBLING RONGEUR
 2 mm x 3 mm bite
 overall length 140 mm




08

RONGEURS

G-18225

KERRISON
BONE NIBBLING RONGEUR
 5 mm x 3.5 mm bite
 overall length 160 mm



GLAUCOMA PUNCHES, RONGEURS

G-18238

KELLY


DESCEMET'S MEMBRANE PUNCH

for Glaucoma

19 gauge / 1.0 mm

0.9 mm x 0.45 mm bite

overall length 122 mm



G-18239

LUNTZ-DODICK


TRABECULECTOMY PUNCH

for Glaucoma

with 360° rotatable cutting head

1 mm x 0.5 mm bite

overall length 115 mm



G-18235

GREHN

MICRO TRABECULECTOMY PUNCH

for Glaucoma

20 gauge / 0.9 mm

0.7 mm x 0.3 mm bite

overall length 90 mm




**G-18237**

GAASS

CORNEOSCLERAL PUNCH

for Glaucoma
 with 360° rotatable cutting head
 1.5 mm x 0.6 mm bite
 overall length 85 mm

**G-18232**

MOHR

GLAUCOMA PUNCH

1 mm circular bite
 overall length 78 mm



SCISSORS



CORNEAL	81
STITCH	82
CAPSULE	84
IRIS	86
IOL	87
TENOTOMY	89
STRABISMUS	90
MISCELLANEOUS · RING HANDLE	91
ENUCLEATION	94
VITREORETINAL	95

TROUTMAN**CORNEAL SCISSORS**

7 mm blades
 9.5 mm from tip to pivot
 overall length 100 mm
G-19710 right, curved blades, blunt tips
G-19715 left, curved blades, blunt tips

**CASTROVIEJO****CORNEAL SCISSORS**

lower blade 9.5 mm
 upper blade 9.0 mm
 overall length 105 mm
G-19690 right, curved blades, blunt tips
G-19695 left, curved blades, blunt tips

**KATZIN****KERATOPLASTY SCISSORS**

6 mm blades
 8 mm from tip to pivot
 overall length 100 mm
G-19720 right, curved blades, blunt tips
G-19725 left, curved blades, blunt tips

**G-19718****SAUTTER****CORNEAL SCISSORS**

light curve, fine blunt tips
 10 mm blades
 14 mm from tip to pivot
 overall length 100 mm

**G-19780****CASTROVIEJO****CORNEAL SCISSORS**

strong curve, blunt tips
 12 mm blades
 15 mm from tip to pivot
 overall length 103 mm

**G-19790****CASTROVIEJO****CORNEAL SCISSORS**

angled to side, blunt tips
 12 mm blades
 overall length 110 mm



G-19600
G-19605

CORNEAL SCISSORS

light curve
8 mm blades
overall length 105 mm
sharp tips
blunt tips



G-19607
G-19615

CORNEAL SCISSORS

strong curve, blunt tips
6 mm blades, overall length 105 mm
11 mm blades, overall length 110 mm



G-19727

MACKENSEN**KERATOPLASTY SCISSORS**

strong curve, blunt tips
overall length 85 mm



G-19585

CORNEAL SECTION SCISSORS

angled up, 8 mm blades
blunt tips
overall length 100 mm



G-19570

MICRO STITCH SCISSORS

straight, 8 mm blades
sharp tips
overall length 105 mm



G-19580

FINE STITCH SCISSORS

straight, 14 mm blades
sharp tips
overall length 110 mm



G-19590
CASTROVIEJO
ANTERIOR SYNECHIA SCISSORS

 angled up
 9 mm vaulted blades, blunt tips
 overall length 105 mm
**G-19740**
CASTROVIEJO
ANTERIOR SYNECHIA SCISSORS

 angled up
 15 mm vaulted blades, blunt tips
 overall length 115 mm
**TRABECULUM SCISSORS**
 sharp tips
 8.5 mm blades
 1.5 mm pivot width
 overall length 80 mm
**G-19745**
 straight
 curved

09

CAPSULOTOMY SCISSORS
 angled up, sharp tips
 2 mm pivot width
 overall length 80 mm
**G-19776**
 9 mm blades
 11 mm blades
G-32100**CAPSULOTOMY SCISSORS**
 angled, lightly curved tips
 13 mm blades
 2 mm pivot width
 overall length 80 mm
**G-32101****CAPSULOTOMY SCISSORS**
 angled, strongly curved tips
 13 mm blades
 2 mm pivot width
 overall length 80 mm
**G-32105****CAPSULOTOMY SCISSORS**
 lightly curved, sharp tips
 11 mm blades
 1.8 mm pivot width
 overall length 88 mm


	CAPSULOTOMY SCISSORS curved, sharp tips 7 mm blades 1.8 mm pivot width overall length 80 mm	
G-19620	<i>TUEBINGEN MODEL</i> CAPSULOTOMY SCISSORS straight, sharp tips 9 mm blades 1.8 mm pivot width overall length 88 mm	
G-19630	<i>TUEBINGEN MODEL</i> CAPSULOTOMY SCISSORS curved, sharp tips 9 mm blades 1.8 mm pivot width overall length 88 mm	
G-19760	VANNAS CAPSULOTOMY SCISSORS extra delicate, straight 6 mm blades, sharp tips 1.7 mm pivot width overall length 80 mm	
G-19765	VANNAS CAPSULOTOMY SCISSORS extra delicate, curved 6 mm blades, sharp tips 1.7 mm pivot width overall length 80 mm	
G-19770	VANNAS CAPSULOTOMY SCISSORS extra delicate, angled to side 7 mm blades, sharp tips overall length 80 mm	

G-19775

VANNAS

CAPSULOTOMY SCISSORS

extra delicate, angled up

7 mm blades, sharp tips

overall length 80 mm

**G-32940**

H.-R. KOCH


CAPSULE SCISSORS

horizontal, 45° angled

2.5 mm blades, 360° rotatable head

22 gauge / 0.7 mm straight shaft

overall length 107 mm

**G-32942**

BRAUWEILER

CAPSULE SCISSORS

vertical, 1.25 mm blades

22 gauge / 0.7 mm curved shaft

overall length 100 mm

**G-38246**

YAE

CAPSULE SCISSORS

interchangeable, titanium, gold-colored

3 mm straight micro blades

23 gauge / 0.6 mm curved shaft

for cutting synechiae, iridotomies,

and to rescue the rhelix during MICS

overall length 42 mm

for use with handle G-38240

**G-38240**

JAPAN MODEL


HANDLE

for capsule scissors and forceps




6 mm, titanium

overall length 90 mm



*OSHER***CAPSULE SCISSORS**

360° rotatable and
interchangeable tips
22 gauge / 0.7 mm
overall length 45 mm
to use with handpiece
G-38215



OSHER
CAPSULE SCISSORS
45° angeld
G-38232

OSHER
CAPSULE SCISSORS
horizontal, 45° angeld
G-38236

OSHER
CAPSULE SCISSORS
horizontal, 90° angeld
G-38238

*OSHER***HANPIPE FOR CAPSULE SCISSORS**

8 mm titanium handpiece, gold-colored
G-38215

**G-32948****HATTENBACH**
HYBRID SCISSORS

to cut sutures during iris and
dislocated IOL fixation
1.25 mm blades
25 gauge / 0.5 mm angled shaft



G-19540
G-19545
G-19550

NOYES
IRIS SCISSORS
straight
18 mm blades
overall length 120 mm
sharp tips
sharp-blunt
blunt tips



G-19640
G-19655

WECKER
IRIS SCISSORS
11 mm blades
overall length 110 mm
sharp tips
blunt tips



WECKER
IRIS SCISSORS

7 mm blades
 overall length 110 mm
 sharp tips
 blunt tips


G-19660
G-19675


G-19680
BARRAQUER
IRIS SCISSORS

7 mm blades
 blunt tips
 overall length 56 mm


G-32945
STEINMETZ
IOL SCISSORS

6.5 mm upper blade
 serrated lower tip with groove
 overall length 92 mm
 For cutting acrylic and silicone foldable IOLs




09

IRIS · IOL

G-32946
STEINMETZ
IOL SCISSORS


for small incision technique
 6.5 mm upper blade
 serrated lower tip with groove
 overall length 92 mm
 For cutting acrylic and silicone foldable IOLs



G-32998
G-32996

SNYDER-OSUMER**IOL SCISSORS**

For cutting acrylic and silicone foldable IOLs
lower notched blade to create IOL traction
18 gauge / 1.2 mm shaft
overall length 150 mm
short model, 3.5 mm blades
long model, 4.5 mm blades



G-32800

H.-R. KOCH

SILICONE LENS BI-SECTOR

for dividing implanted silicone lenses
overall length without wire loop 130 mm
5 wire loops included (G-32811)



G-32810
G-32811

WIRE LOOP

for *H.-R. KOCH* silicone lens bi-sector
(G-32800)
overall length 65 mm (loop 12 mm)
1 piece
5 pieces



G-17440
G-17441

COMBINED**NEEDLE HOLDER / SCISSORS**

curved 12 mm blades with
3 mm needle holder and tying platform
used to hold, place, and cut fine suture material
overall length 125 mm



G-19750 **WESTCOTT****SCISSORS**

blunt tips

curved, standard blades

overall length 120 mm

**G-19751** **WESTCOTT****SCISSORS**

pointed tips

curved, standard blades

overall length 120 mm

**G-32250** **WESTCOTT****SCISSORS**

blunt tips

curved, delicate blades

overall length 100 mm

**G-19470** **STEVENS****TENOTOMY SCISSORS**

short blades, blunt tips

22 mm from tip to pivot

overall length 105 mm

G-19475 straight

curved

**G-19480** **STEVENS****TENOTOMY SCISSORS**

long blades, blunt tips

28 mm from tip to pivot

overall length 110 mm

G-19485 straight

curved



G-32614*BONN MODEL***TENOTOMY SCISSORS**

curved blades, blunt tips

24 mm from tip to pivot

overall length 90 mm

**G-19430****G-19435***BONN MODEL***STRABISMUS SCISSORS**

blunt tips

24 mm from tip to pivot

overall length 90 mm

straight

curved

**G-19420****G-19425****STRABISMUS SCISSORS**

long blade

28 mm from tip to pivot

overall length 105 mm

straight

curved



STRABISMUS SCISSORS

short blade
 28 mm from tip to pivot
 overall length 110 mm

G-19460
G-19465

straight
 curved

**EYE SCISSORS**

straight, sharp tips
 overall length: tip to pivot:
G-19350 90 mm 27 mm
G-19355 100 mm 29 mm
G-19360 110 mm 30 mm
G-19365 120 mm 36 mm

**EYE SCISSORS**

curved, sharp tips
 overall length: tip to pivot:
G-19370 90 mm 25 mm
G-19375 100 mm 28 mm
G-19380 110 mm 29 mm
G-19385 120 mm 34 mm



G-19390
G-19395

EYE SCISSORS

delicate, sharp tips
31 mm from tip to pivot
overall length 105 mm
straight
curved

*BONN MODEL***EYE SCISSORS**

extra delicate
sharp tips
19 mm from tip to pivot
overall length 90 mm
curved
straight

G-19400
G-19405



G-19410

SPENCER
STITCH SCISSORS

20 mm straight blades
with hook to lift and
secure sutures before cutting
overall length 90 mm



G-19450**GRAEFE
EYE SCISSORS**

angled up
blunt tips
overall length 90 mm

**G-19505****EYE SCISSORS**

angled to side, probe tip
sharp-blunt
overall length 115 mm

**G-19510****EYE SCISSORS**

angled to side, probe tips
blunt tips
overall length 115 mm



G-31315**OPERATING SCISSORS**

straight blades

sharp-blunt

overall length 145 mm

**G-19520****G-19525****G-19530****ENUCLEATION SCISSORS**

curved blades, blunt

40 mm from tip to pivot

overall length 135 mm

strong curve


medium curve

light curve



20 GAUGE
VITREORETINAL SCISSORS


20 gauge / 0.9 mm
 titanium handle
 gold-colored
 overall length 155 mm




45° angle
 horizontal
G-36502 6 mm handle
G-36504 8 mm handle




curved
 with membrane peeler
G-36506 6 mm handle
G-36508 8 mm handle



curved
G-36512 6 mm handle
G-36514 8 mm handle

09



KROLL
 angled, with membrane peeler
G-36516 6 mm handle
G-36518 8 mm handle




angled, blunt tips
 for MIS
G-36522 6 mm handle
G-36524 8 mm handle




HEIMANN
 angled
G-36526 6 mm handle
G-36528 8 mm handle




HEIMANN
 straight
G-36532 6 mm handle
G-36534 8 mm handle



HEIMANN
 straight, blunt tips
G-36536 6 mm handle
G-36538 8 mm handle




60° angled
G-36546 6 mm handle
G-36548 8 mm handle

VITREORETINAL 20 GAUGE


curved
 horizontal
G-36552 6 mm handle
G-36554 8 mm handle



NAITO
 long blade 3.5 mm
 curved, tapered tip
G-36556 6 mm handle
G-36558 8 mm handle



ROIDER
 for RPE transplantation, vitrectomies
 and retinectomies
G-36572 6 mm handle
G-36574 8 mm handle

23 GAUGE


VITREORETINAL SCISSORS

23 gauge / 0.6 mm


titanium handle

gold-colored

overall length 155 mm




curved



HEIMANN


straight

G-36576 6 mm handle**G-36578** 8 mm handle

HEIMANN

for use with 23 gauge trocar system

angled, sharp tips

G-36540 6 mm handle**G-36541** 8 mm handle

HEIMANN

straight, blunt tips

G-36582 6 mm handle**G-36584** 8 mm handle


25 GAUGE

VITREORETINAL SCISSORS


25 gauge / 0.5 mm

titanium handle

gold-colored



curved



HEIMANN

straight

G-36562 6 mm handle**G-36564** 8 mm handle**G-36586** 6 mm handle**G-36588** 8 mm handle


ERGOLINE

ERGOLINE SCISSORS

360° rotatable and interchangeable tips
for Ergoline handles 6 and 8 mm diameter
overall length 55 mm

ERGOLINE HANDPIECE

for Ergoline scissors
gold-colored
overall length 105 mm
G-36900 6 mm titanium handpiece
G-36901 8 mm titanium handpiece




ERGOLINE SCISSORS

45° angled
horizontal
20 gauge / 0.9 mm
G-36902




KROLL

ERGOLINE SCISSORS
angled, with membrane peeler
20 gauge / 0.9 mm
G-36915



HEIMANN
ERGOLINE SCISSORS
angled
20 gauge / 0.9 mm
G-36925



HEIMANN ERGOLINE SCISSORS

straight
G-36930 20 gauge / 0.9 mm
G-36975 23 gauge / 0.6 mm
G-36985 25 gauge / 0.5 mm



NAITO

ERGOLINE SCISSORS
long blade 3.5 mm, curved
tip tapered to 0.75 mm
20 gauge / 0.9 mm
G-36955 20 gauge / 0.9 mm



CLEANING ADAPTER
for rotatable and
interchangeable tips

G-36604

**SINGLE-USE
INTERCHANGEABLE SCISSORS
AVAILABLE ON P. 218**

UNO
COLORLINE

SCISSORS

13

FORCEPS



CILIA	101
CONJUNCTIVA	102
SUTURING	103
CORNEAL	104
TYING	105
FIXATION	109
IRIS	110
CAPSULORHEXIS	113
IOL	121
MUSCLE	126
VITREORETINAL	129
CHALAZION	134
MOSQUITO	135
SERREFINES · TOWEL CLAMPS	136

- G-18720** *HENRY*
CILIA FORCEPS
 straight jaws
 smooth, blunt, squared tips
 overall length 100 mm



- G-18730** *GRADLE*
CILIA FORCEPS
 curved jaws
 smooth, blunt, squared tips
 overall length 98 mm



- G-18750** *BARRAQUER*
CILIA FORCEPS
 lightly angled
 6.5 mm platform
 overall length 109 mm



- G-18740** *BEER*
CILIA FORCEPS
 3 mm platform
 rounded tips 4.5 mm wide
 overall length 78 mm



- G-18850** *SHAAF*
CILIA AND FOREIGN BODY FORCEPS
 grooved tips
 curved and pointed
 finger grip handle
 overall length 96 mm



- G-31620** *BARRAQUER*
CILIA AND SUTURE FORCEPS
 grooved tips, 7.0 mm platform
 curved and pointed
 finger grip handle
 overall length 112 mm




- G-32740** *HAMBURG MODEL*
CILIA AND FOREIGN BODY FORCEPS
 tapered fine tips
 overall length 90 mm




G-31230

KREFELD MODEL
CONJUNCTIVA FORCEPS
 standard
 1 x 2 teeth
 overall length 100 mm




G-31235

FREIBURG MODEL
CONJUNCTIVA FORCEPS
 delicate
 1 x 2 teeth
 overall length 90 mm




G-18847

KAUFMANN
CONJUNCTIVA FORCEPS
 standard
 1 x 2 teeth
 5 mm tying platform
 overall length 75mm




G-18810

DRAEGER
CONJUNCTIVA FORCEPS
 1.2 mm wide
 2 atraumatic fixation knobs
 overall length 97 mm



G-18848



FOERSTER
CONJUNCTIVA FORCEPS
 delicate
 1 x 2 teeth
 finger grip handle
 overall length 95 mm






CASTROVIEJO
SUTURING FORCEPS
 1 x 2 teeth
 6 mm tying platform
 overall length 100 mm

G-19001 0.28 mm, teeth 0.12 mm
G-19005 0.3 mm
G-19010 0.6 mm
G-19015 0.9 mm
G-19020 1.2 mm




CASTROVIEJO
SUTURING FORCEPS
 new model
 1 x 2 teeth
 6 mm tying platform
 overall length 109 mm

G-19022 0.25 mm, teeth 0.12 mm
G-19023 0.3 mm
G-19024 0.6 mm
G-19025 0.9 mm
G-19026 1.2 mm




G-19008 *MUNICH MODEL*
SUTURING FORCEPS
 1 x 2 teeth, 0.35 mm
 5.5 mm tapered tying platform
 overall length 100 mm




ZUERICH MODEL
SUTURING FORCEPS

delicate
1 x 2 teeth, 0.12 mm
5 mm tying platform
overall length 75 mm


G-18965
G-18966




G-18950 **COLIBRI FORCEPS**
1 x 2 teeth, 0.12 mm
5 mm tying platform
overall length 75 mm



G-18960 **MICRO COLIBRI FORCEPS**
delicate
1 x 2 teeth, 0.12 mm
5 mm tying platform
overall length 77 mm



G-15134 **MICRO COLIBRI FORCEPS**
extra delicate
1 x 2 teeth, 0.12 mm
4 mm tying platform
overall length 76 mm



G-18975 **BARRAQUER**
SUTURE REMOVING FORCEPS
very fine pointed tips
4 mm tying platform
overall length 77 mm



G-18972 **THOMAS NEUHANN**
KERATOPLASTY FORCEPS
for corneal suturing and fixation
two prongs with 1 x 2 teeth
0.8 mm prong spread
overall length 80mm




**G-18970**

REMKY
KERATOPLASTY FORCEPS
 two prongs with 1 x 1 teeth, angled
 1.7 mm prong spread
 for corneal suturing and fixation
 overall length 80 mm

**G-38623**


REMKY
DMEK KERATOPLASTY FORCEPS
 modified by SZURMAN
 two prongs with 1 x 1 teeth, straight
 1.7 mm prong spread
 33 mm shafts
 to hold graft during preparation
 overall length 89 mm




10


G-18940
G-18945

PAUFIQUE
UNIVERSAL SUTURING FORCEPS
 1 x 2 teeth
 6 mm tying platform
 overall length 88 mm



G-19045
G-19046

PAUFIQUE
TYING FORCEPS
 with platform
 overall length 88 mm
 standard, 6 mm platform
 delicate, 5 mm platform




G-19030**TUEBINGEN MODEL****TYING FORCEPS**

straight, smooth jaws
6 mm tying platform
overall length 101 mm

**G-19040****TUEBINGEN MODEL****TYING FORCEPS**

curved, smooth jaws
6 mm tying platform
overall length 99 mm




10


TYING

G-19035**TUEBINGEN MODEL****TYING FORCEPS**

delicate, smooth jaws
4 mm tying platform
overall length 84mm


**G-19034****BONN MODEL****TYING FORCEPS**

extra delicate, smooth jaws
3 mm tying platform
overall length 73 mm




**G-19032****TYING FORCEPS**

angled, 7.5 mm
 6 mm tying platform
 overall length 99 mm

**G-31625***DARDENNE***TYING FORCEPS**

angled, 6.5 mm
 6 mm tying platform
 round handle with guide pin
 overall length 102 mm

**G-19038***BANGERTER***TYING FORCEPS**

9.5 mm platform
 longitudinal v-shape groove
 overall length 100 mm

**G-19048****TYING FORCEPS**

utility serrated jaws
 1 mm wide
 overall length 106 mm



G-32947*HATTENBACH***HYBRID TYING FORCEPS**


angled, for the anterior chamber
 intraocular 1.5 mm micro tying tips
 25 gauge / 0.5 mm
 overall length 98 mm




for use in combination with
 hybrid needle holder G-32949

G-31623*KELMAN-MCPHERSON***FRAGMENT FORCEPS**


angled, smooth jaws
 7.5 mm platform
 overall length 100 mm

**G-33941***KELMAN-MCPHERSON***FRAGMENT FORCEPS**

angled, fine grooved jaws,
 for femto-LASIK to retrieve lenticule
 during SMILE procedures
 7.5 mm platform
 overall length 100 mm

**G-31622***GASKIN***FRAGMENT FORCEPS**

angled
 with tying platform
 9 mm platform, overall length 107 mm
G-31624
 11 mm platform, overall length 100 mm





G-37700 JEWELERS TYPE FORCEPS
 standard, straight
 overall length 120 mm



G-37705 JEWELERS TYPE FORCEPS
 delicate, straight
 overall length 115 mm



G-37710 JEWELERS TYPE FORCEPS
 extra delicate, straight
 overall length 115 mm

10

G-19080 BARRAQUER FIXATION FORCEPS
 3 sets of 1 x 2 teeth
 11 mm fixation platform
 overall length 100 mm



G-19050 GRAEFE FIXATION FORCEPS
 4.5 mm wide
 8 interlocking teeth
 overall length 110 mm
 without lock
G-19055 with lock




G-19060 HEIDELBERG MODEL FIXATION FORCEPS
 3.5 mm wide
 6 interlocking teeth
 overall length 110 mm



G-18780
G-18781
G-18782
G-18783
G-18784

DRESSING FORCEPS


serrated jaws
overall length 100 mm



G-18790
G-18791
G-18792
G-18793
G-18794

TISSUE FORCEPS


1 x 2 teeth
overall length 100 mm



G-18855
G-18856


IRIS FORCEPS

serrated, delicate
overall length 73 mm
straight
curved




G-18908
G-18906
G-18911

BONN MODEL
IRIS FORCEPS
serrated, extra delicate
overall length 73 mm
straight
angled at 6 mm
angled at 8 mm




G-18900
G-18905
G-18910

BONN MODEL
IRIS FORCEPS
1 x 2 teeth, 0.12 mm
extra delicate
overall length 73mm
straight
angled at 6 mm
angled at 8 mm




*BONN MODEL***IRIS FORCEPS**


long, straight
1 x 2 teeth, 0.12 mm
overall length 102 mm

G-18845**G-18890** delicate, with 4.5 mm tying platform**G-18840***SWISS MODEL***IRIS FORCEPS**


curved, 1 x 2 teeth
overall length 100 mm

**G-18860***HEIDELBERG MODEL***IRIS FORCEPS**


light curve, 1 x 2 teeth
overall length 72 mm

**G-18870***HEIDELBERG MODEL***IRIS FORCEPS**


strong curve, 1 x 2 teeth
overall length 70 mm

**G-18857****IRIS FORCEPS**


straight, delicate
1 x 2 teeth
overall length 72 mm

**G-31953****IRIS FORCEPS**


straight, strong
1 x 2 teeth
overall length 72 mm




G-18846 *HAMBURG MODEL*
IRIS FORCEPS
 straight, delicate
 1 x 2 teeth
 overall length 103 mm




G-32240 *SWISS MODEL*
IRIS FORCEPS
 straight, extra delicate
 1 x 2 teeth
 overall length 70 mm




G-18880 *HESS-ENGELKING*
IRIS FORCEPS
 1 x 2 teeth
 overall length 72 mm



G-18897 *ELZE*
IRIS FORCEPS
 extra delicate
 1 x 2 teeth
 overall length 78 mm




G-18913 *THOMAS NEUHANN*
IRIS FORCEPS
 curved shaft
 1 x 2 teeth
 23 gauge / 0.6 mm
 overall length 95 mm



G-31307

CAPSULORHEXIS FORCEPS

strong, with heel stop to protect the iris
 interlocking micro serrated tips
 12 mm straight shafts
 overall length 103 mm



G-31308

CAPSULORHEXIS FORCEPS

round handle with guide pin
 extra delicate with
 interlocking micro serrated tips
 13 mm curved shafts
 overall length 105 mm



10

CAPSULORHEXIS

G-31300

MASKET**CAPSULORHEXIS FORCEPS**

delicate, with heel stop to protect the iris
 interlocking micro serrated tips
 12 mm curved shafts
 overall length 103 mm



G-31304

CAPSULORHEXIS FORCEPS


new model, with guide pin
 extra delicate with
 interlocking micro serrated tips
 13 mm curved shafts
 overall length 103 mm



G-31301

CAPSULORHEXIS FORCEPS

extra delicate
with smooth triangular grasping tips
13 mm curved shafts
overall length 98 mm



G-31305

CAPSULORHEXIS FORCEPS


extra delicate with
interlocking micro serrated tips
13 mm curved shafts
overall length 98 mm



G-31303

*H.-R. KOCH***CAPSULORHEXIS FORCEPS**


extra delicate with
interlocking micro serrated tips
13 mm curved shafts
1 mm marks from tip to bent
overall length 98 mm



G-31299

CAPSULORHEXIS FORCEPS

short handle, extra delicate
cross serrated grasping tips
12 mm straight shafts
overall length 83 mm




**G-32926**


AUFFARTH
MICS DOUBLE CROSS ACTION CAPSULORHEXIS FORCEPS
 interlocking micro serrated tips
 12 mm curved shafts, angled 130°
 with overlap system
 prevents inadvertent iris capture
 for 1.6 mm incisions
 overall length 100 mm

**G-32928**


INAMURA
MICS CAPSULORHEXIS FORCEPS
 cross action
 interlocking micro serrated tips
 10.5 mm curved shafts tip to pivot
 angled 130°
 overall length 103 mm

**10****G-31306**

KAMPIK
MICS CAPSULORHEXIS FORCEPS
 22 gauge / 0.7 mm curved shaft
 90° cystotome tips, reverse action
 overall length 93 mm


**CAPSULORHEXIS · MICRO INCISION****G-32935**

H.-R. KOCH
MICS CAPSULE FORCEPS
 short handle
 22 gauge / 0.7 mm curved shaft with
 micro grasping tips
 overall length 98 mm




**G-32938***H.-R. KOCH***MICS CAPSULE FORCEPS**


22 gauge / 0.7 mm curved shaft with
micro grasping tips
6 mm titanium handle, blue-colored
overall length 144 mm

**G-32936***TOBIAS NEUHANN***MICRO FORCEPS**

multipunction asymmetrical tips with
grasping platform
23 gauge / 0.6 mm straight shaft
6 mm titanium handle, blue-colored
overall length 155 mm

**G-32937***IKEDA***MICS CAPSULORHEXIS FORCEPS**

23 gauge / 0.6 mm curved shaft with
micro grasping tips
6 mm titanium handle, blue-colored
overall length 144 mm

**G-38244
G-38242***IKEDA***CAPSULORHEXIS FORCEPS**

interchangeable tips
23 gauge / 0.6 mm curved shaft
overall length 40 mm
textured, micrograsping tips **"BEAK"**
textured, micrograsping tips

**G-38244**
"beak"**G-38242**

for use with handle G-38240


JAPAN MODEL**HANDLE**

for capsule scissors and forceps
interchangeable tips
6 mm, titanium
overall length 90 mm




G-32944

HATTENBACH
CAPSULORHEXIS FORCEPS
 25 gauge / 0.5 mm angled shaft with
 micro grasping tips
 overall length 88 mm

**G-18805**


TOBIAS NEUHANN
DESCEMET MEMBRANE FORCEPS
 delicate blunt grasping tips
 with platform
 overall length 75 mm

**G-38624**

SZURMAN
DMEK PREPARATION FORCEPS
 for removal of Descemet layer
 with smooth round tips to prevent tearing
 the donor endothelium
 overall length 75 mm

**G-38628**

YÖRÜK
DMEK TRANSPLANT FORCEPS
 atraumatic grasping and removal of
 Descemet's membrane
 grasping area 6 mm
 overall length 77 mm

**G-38634**


GERTEN
DSAEK/DMEK RHESIS FORCEPS
 for Descemetoxyhexis
 reversed triangular tips to
 score and remove the DM from the host cornea
 23 gauge / 0.6 mm curved shaft
 overall length 142 mm



G-38603

DSAEK IMPLANTATION FORCEPS

13 mm curved shafts with limited closing
to avoid injury of endothelial transplant
1.5 mm serrated platform
overall length 107 mm



Geuder

G-38606

*THOMAS NEUHANN***DSAEK FORCEPS**

1 x 2 teeth, 0.12 mm
21 gauge / 0.8 mm curved shaft
for endothelial transplant procedures
overall length 97 mm



Geuder

G-33428

*PFAEFFL***FLAPORHEXIS FORCEPS**

2.5 mm x 1.5 mm cross serrated tips, oval
to lift femto-LASIK flaps
overall length 103 mm



Geuder

G-33429


*KNORZ***FEMTO-LASIK FLAP FORCEPS**

curved 3.25 mm T-bar
textured
to grasp and lift LASIK flaps
overall length 102 mm




**G-31899 LASIK CORNEAL FLAP FORCEPS**

angled shafts, 9 mm
highly polished round and smooth shafts
overall length 108 mm

**G-18914 BURATTO****LASIK CORNEAL FLAP FORCEPS**


curved shafts
cross serrated tips, 0.8 mm diameter
overall length 100 mm

**G-33960 RUPAL SHAH****LENTICULE FORCEPS**

3.5 mm serrated tips
20 gauge / 0.9 mm curved shaft
to retrieve lenticule during SMILE procedures
6 mm titanium handle, blue-colored
overall length 144 mm

**G-33961 LENTICULE FORCEPS**

7.5 mm shafts
with fine grooves
angled to side
to retrieve lenticule during SMILE procedures
overall length 99 mm




G-31518

NUCLEUS SPREADING FORCEPS

cross-action

3.0 mm textured paddle tips
for cracking the nucleus
overall length 112 mm



G-31519

AKAHOSHI**NUCLEUS SPREADING FORCEPS**

phaco prechopper

4.5 mm delicate sharp tips, 1 mm profile
2.0 mm cross-action
for cracking, dividing, and rotating the nucleus
overall length 105 mm



G-31720

AKAHOSHI**PRECHOPPER**


for micro-incisions, straight shafts, cross action
round handle
overall length 125 mm



G-31722

AKAHOSHI**HYBRID COMBO PRECHOPPER**

for micro-incisions, cross action, straight
round handle
length of tip 5.3 mm
overall length 120 mm



**G-31724***AKAHOSHI***UNIVERSAL PRECHOPPER**

for micro-incisions, cross action, straight
 round handle
 overall length 120 mm

**G-33440***VON BARSEWISCH***ANTERIOR CHAMBER FORCEPS**

extra delicate
 8.5 mm angled shafts
 1.5 mm platform
 overall length 72 mm

**G-32165***FECHNER***IMPLANTATION FORCEPS**

smooth curved shafts
 overall length 80 mm

**G-18935***DODEN***IMPLANTATION FORCEPS**

curved
 delicate lens holding shafts
 overall length 110 mm




**G-18930***SHEPARD***IMPLANTATION FORCEPS**

curved, 1.4 mm disc-shaped tips
delicate lens holding shafts
overall length 120 mm

**G-31962***T. KOHNEN***IMPLANTATION FORCEPS**

silicone tube coated
10 mm angled shafts
overall length 107 mm

**G-18918***DICK***FIXATION AND HOLDING FORCEPS**

for Artisan/Verisyse IOL
with guide mark for lens maximum hold
and stability
18 gauge / 1.2 mm curved shaft
overall length 143

**G-18919***DICK***IOL HAPTIC EXPANDER**

for enclavation of iris-fixed IOLs
21 gauge / 0.8 mm curved shaft
overall length 140 mm





**G-18944**

T. NEUHANN

IMPLANTATION FORCEPS

for retropupillary implantation of Iris-Claw-IOL
 0.9 mm wide holding surface
 20 gauge / 0.9 mm curved shaft
 overall length 143 mm

**G-18928**

SEKUNDO

IMPLANTATION FORCEPS

for retropupillary implantation of Iris-Claw-IOL
 20 gauge / 0.9 mm curved shaft with
 lens holding plate, three-point fixation
 and IOL centration guide mark
 6 mm titanium handle, blue-colored
 overall length 143 mm



10

**G-31902****FOLDER**


foldable IOL insertion forceps
 cross action
 6.5 mm concave jaws
 overall length 103 mm

**G-31903**

SPAЛЕCK

FOLDER

foldable IOL insertion forceps
 for silicone and acrylic IOLs
 angled shafts, concave jaws
 overall length 102 mm




G-31908**FOLDER**

IOL folding forceps
for easy pickup and folding acrylic IOLs
overall length 105 mm

**G-31994***H.-R. KOCH***PREFOLDER**

for foldable silicone, hydrogel or acrylic IOLs
overall length 135 mm




10

G-31905*H.-R. KOCH***FOLDER**

for foldable 3-piece IOLs
6.0 mm concave jaws
overall length 92 mm


**G-31901****HOLDING FORCEPS**

for foldable IOLs
highly polished round and smooth shafts
overall length 108 mm




**G-18941***STROBEL***HOLDING FORCEPS**


for Pharmacia CeeOn Edge IOL
overall length 110 mm

**G-32823****EXPLANTATION FORCEPS**


for soft IOLs
2.0 mm x 1.4 mm paddle tip with
interlocking serrations and guide pin
overall length 108 mm

**G-32997***SNYDER-OSHER***IOL FORCEPS**

to hold the IOL while cutting it inside
the capsular bag
2 mm long jaws with 0.5 mm serrated tips
22 gauge / 0.7 mm lightly curved shaft
overall length 99 mm


**G-18820***REMKY***SUPERIOR RECTUS FORCEPS**

angled
1 x 2 teeth
overall length 109 mm







G-18800 *ELSCHNIG*
SUPERIOR RECTUS FORCEPS
 standard, straight
 1 x 2 teeth
 overall length 110 mm




G-18795 *ELSCHNIG*
SUPERIOR RECTUS FORCEPS
 delicate, straight
 1 x 2 teeth
 overall length 100 mm



G-19105 *LEYDHECKER*
SUPERIOR RECTUS FORCEPS
 curved, tip width 0.8 mm
 1 x 2 teeth, 1.1 mm
 overall length 110 mm



G-32335 *CASTROVIEJO*
SUPERIOR RECTUS FORCEPS
 angled with tying platform
 1 x 2 teeth, 0.9 mm
 overall length 100 mm



G-19112 *BANGERTER*
MUSCLE FORCEPS
 10 teeth
 self retaining
 overall length 100 mm



*SATTLER***MUSCLE FORCEPS**

with slide lock

overall length 100 mm

G-19120

3 teeth, bulbous tip, right

G-19125**G-32350***TROELTSCH***NASAL FORCEPS**

serrated tips

for packing and dressing during ENT procedures

overall length 150 mm

**G-19107***ESSEN MODEL***CERCLAGE FORCEPS**

strongly curved

3 mm platform with dual grooves

overall length 104 mm



10


G-33752*WATZKE***SPREADING FORCEPS**

to stretch silicone sleeves

cross action handle and

serrated tips


overall length 115 mm

**G-33434***PFENNIGWERTH***HOLDING FORCEPS**

for vitreoretinal lenses

curved with concave lens holding tips


overall length 110 mm




MUSCLE · NASAL · RETINAL

G-33439**SCLERAL PLUG FORCEPS**

cross action handle, angled jaws
to securely hold plugs during insertion
and removal
overall length 105 mm

**G-31570****SCLERAL PLUG FORCEPS**


to securely hold plugs during insertion
and removal
overall length 100 mm




10



G-33432
SCLERAL PLUG
serrated,
20 gauge / 0.9 mm



G-33436
SCLERAL PLUG
to suture onto sclera
20 gauge / 0.9 mm




G-33435
SCLERAL PLUG
20 gauge / 0.9 mm




G-33441
SCLERAL PLUG
23 gauge / 0.6 mm

**G-33438****FORCEPS**

for HEIMANN retinal tacks
overall length 100 mm

**G-33437**

HEIMANN
RETINAL TACK
titanium
- for single use only -




G-37420*BARTZ-SCHMIDT*


TROCAR CLOSURE FORCEPS
 with function for bulbus fixation
 overall length 72 mm

**G-37430****TROCAR FORCEPS**

for trocar systems with groove
 overall length 110 mm

**G-33831***REIMANN*


VITREORETINAL FORCEPS
 for foreign body removal
 with retractable magnet and 3-prongs
 14 gauge / 2.1 mm straight shaft
 shaft length 35 mm
 overall length 138 mm




20 GAUGE

VITREORETINAL FORCEPS


20 gauge / 0.9 mm
6 mm and 8 mm diameter handles
titanium, blue-colored
overall length 155 mm




straight
serrated tips
G-36202 6 mm handle




straight
textured tips
G-36206 6 mm handle
G-36208 8 mm handle




straight
serrated "crocodile" tips
G-36256 6 mm handle
G-36258 8 mm handle




straight
silicone coated
G-36212 6 mm handle
G-36214 8 mm handle



TANO
straight, serrated tips
G-36356 6 mm handle
G-36358 8 mm handle




A. MUELLER
for retinal vessels / sheatotomy
G-36332 6 mm handle




straight
grip tips
G-36216 6 mm handle
G-36218 8 mm handle





90° angle
pointed tips
G-36276 6 mm handle
G-36278 8 mm handle




NAI/TO
ILM
G-36316 6 mm handle
G-36318 8 mm handle




straight
endgripping, extra delicate
0.2 mm x 0.05 mm
G-36292 6 mm handle
G-36294 8 mm handle




straight
endgripping, extra delicate
0.3 mm x 0.05 mm
G-36296 6 mm handle
G-36298 8 mm handle



straight
endgripping, delicate
0.2 mm x 0.1 mm
G-36302 6 mm handle
G-36304 8 mm handle





straight
endgripping
0.3 mm x 0.1 mm
G-36306 6 mm handle
G-36308 8 mm handle



20 GAUGE
VITREORETINAL FORCEPS


20 gauge / 0.9 mm
 6 mm and 8 mm diameter handles
 titanium, blue-colored
 overall length 155 mm




90° angle
 membrane peeler
G-36226 6 mm handle



45° angle
 membrane peeler
G-36232 6 mm handle



45° angle
G-36222 6 mm handle




membrane peeler tips
G-36252 6 mm handle
G-36254 8 mm handle



curved
 membrane dissector
G-36266 6 mm handle




45° angle
 membrane dissector
G-36272 6 mm handle




45° angle
 horizontal subretinal
G-36236 6 mm handle




45° angle, vertical tips
 subretinal
G-36246 6 mm handle
G-36248 8 mm handle




angled, extra delicate
 3.0 mm tips
G-36336 6 mm handle



angled, extra delicate
 4.5 mm tips
G-36326 6 mm handle
G-36328 8 mm handle



straight
 circular tips
G-36242 6 mm handle




foreign body forceps
G-36262 6 mm handle
G-36264 8 mm handle

23 GAUGE


VITREORETINAL FORCEPS

23 gauge / 0.6 mm
6 mm and 8 mm diameter handles
titanium, blue-colored
overall length 155 mm




straight
serrated "crocodile" tips

G-36342 6 mm handle
G-36344 8 mm handle




TANO
straight, serrated tips

G-36362 6 mm handle
G-36364 8 mm handle




straight
endgripping, extra delicate
0.2 mm x 0.05 mm

G-36282 6 mm handle
G-36284 8 mm handle



straight
grip tips


G-36286 6 mm handle
G-36288 8 mm handle



25 GAUGE


VITREORETINAL FORCEPS

25 gauge / 0.5 mm
6 mm and 8 mm diameter handles
titanium, blue-colored
overall length 155 mm




straight
serrated "crocodile" tips

G-36322 6 mm handle
G-36324 8 mm handle




TANO
straight, serrated tips

G-36366 6 mm handle
G-36368 8 mm handle



straight
endgripping, extra delicate
0.2 mm x 0.05 mm

G-36312 6 mm handle
G-36314 8 mm handle









ERGOLINE

ERGOLINE FORCEPS

360° rotatable and interchangeable tips
overall length 55 mm
for Ergoline handpieces 6 mm and 8 mm diameters

ERGOLINE HANDPIECE

for Ergoline forceps
blue-colored
G-36600 6 mm titanium handpiece
G-36601 8 mm titanium handpiece



ERGOLINE FORCEPS

straight
serrated "crocodile" tips
G-36655 20 gauge / 0.9 mm
G-36740 23 gauge / 0.6 mm
G-36720 25 gauge / 0.5 mm

TANO







ERGOLINE FORCEPS
straight, serrated tips
G-36755 20 gauge / 0.9 mm
G-36765 23 gauge / 0.6 mm

NAITO

ERGOLINE FORCEPS ILM
G-36715 20 gauge / 0.9 mm

ERGOLINE FORCEPS

endgripping
extra delicate 0.2 mm x 0.05 mm
G-36690 20 gauge / 0.9 mm
G-36680 23 gauge / 0.6 mm
G-36710 25 gauge / 0.5 mm



ERGOLINE FORCEPS

endgripping
extra delicate 0.3 mm x 0.05 mm
G-36695 20 gauge / 0.9 mm

ERGOLINE FORCEPS

endgripping
delicate 0.2 mm x 0.1 mm
G-36700 20 gauge / 0.9 mm

ERGOLINE FORCEPS

endgripping
delicate 0.3 mm x 0.1 mm
G-36705 20 gauge / 0.9 mm

**SINGLE-USE
INTERCHANGEABLE FORCEPS
AVAILABLE ON P. 219**

UNO
COLORLINE

G-18680*LAMBERT***CHALAZION FORCEPS**

round, 9 mm inside diameter
solid and fenestrated plates
overall length 92 mm

**G-18670***AYER***CHALAZION FORCEPS**

round, 8 mm inside diameter
interlacing tips
solid and fenestrated plates
overall length 90 mm

**G-18685***DESMARRES***CHALAZION FORCEPS**

oval
solid and fenestrated plates
20 mm, overall length 88 mm
G-18690 25 mm, overall length 90 mm
G-18695 31 mm, overall length 95 mm

**G-18690****G-18695****G-18700***CAUER***CHALAZION FORCEPS**

oval, 13 mm x 11 mm inside diameter
interlacing serrated tips
solid and fenestrated plates
overall length 99 mm

**G-18710***LORDAN***CHALAZION FORCEPS**

triangular, 13 mm x 8 mm inside diameter
solid and fenestrated plates
overall length 95 mm

**G-18715***BERKE***PTOSIS FORCEPS**

22 mm long serrated tips
slide lock mechanism
overall length 110 mm





HARTMANN

MOSQUITO FORCEPS

serrated
18 mm jaws
overall length 100 mm
G-18170
G-18171
straight
curved



HALSTEAD

MOSQUITO FORCEPS

G-18180
G-18181
serrated
25 mm jaws
overall length 125 mm
straight
curved



HALSTEAD

MOSQUITO FORCEPS

G-18185
G-18182
1 x 2 teeth
serrated
25 mm jaws
overall length 125 mm
straight
curved

G-18160

BACKHAUS
TOWEL FORCEPS
 overall length 90 mm

**G-33200**

SCHAEDEL
TOWEL FORCEPS
 overall length 90 mm



10

G-18190

DIEFFENBACH
SERREFINE
 straight, overall length 35 mm

**G-18195**

DIEFFENBACH
SERREFINE
 curved, overall length 35 mm

**G-31420**

DIEFFENBACH
SERREFINE
 straight, overall length 50 mm

**G-31415**

DIEFFENBACH
SERREFINE
 curved, overall length 50 mm

**G-18197**

SERREFINE
 extra delicate, overall length 30 mm



11

NEEDLE HOLDERS



Geuder®
Precision made in Germany

BARRAQUER**NEEDLE HOLDER**

10 mm delicate jaws
 0.8 x 0.55 mm tip
 overall length 120 mm
 curved, without lock
 curved, with lock

G-17500
G-17501

**BARRAQUER****NEEDLE HOLDER**

10 mm delicate jaws
 0.8 x 0.45 mm tip
 overall length 120 mm
 curved, without lock
 curved, with lock

G-17505
G-17506

**JACOBI****NEEDLE HOLDER**

5 mm extra delicate jaws
 1.0 x 0.5 mm tips
 overall length 115 mm
 straight, without lock
 straight, with lock

G-17451
G-17453



NEEDLE HOLDERS

11

JACOBI**NEEDLE HOLDER**

5 mm extra delicate jaws
 1.0 x 0.5 mm tips
 overall length 115 mm
 curved, without lock
 curved, with lock

G-17450
G-17452

**BARRAQUER****NEEDLE HOLDER**

7.5 mm extra delicate jaws
 0.6 x 0.3 mm tip
 overall length 105 mm
 straight, without lock
 straight, with lock

G-17461
G-17463

**BARRAQUER****NEEDLE HOLDER**

7.5 mm extra delicate jaws
 0.6 x 0.3 mm tip
 overall length 105 mm
 curved, without lock
 curved, with lock

G-17460
G-17462



G-32260
G-32261

BARRAQUER-TROUTMAN**NEEDLE HOLDER**

10 mm extra delicate jaws
0.8 x 0.6 mm tip
handle diameter 6 mm
overall length 100 mm
straight, without lock
straight, with lock



G-17535
G-17530

BARRAQUER-TROUTMAN**NEEDLE HOLDER**

10 mm extra delicate jaws
0.8 x 0.6 mm tip
handle diameter 6 mm
overall length 100 mm
curved, without lock
curved, with lock



G-17510
G-17511

BARRAQUER**NEEDLE HOLDER**

11 mm delicate jaws
0.9 x 0.8 mm tips
handle diameter 8 mm
overall length 140 mm
curved, without lock
curved, with lock



G-17515
G-17516

BARRAQUER**NEEDLE HOLDER**

11 mm extra delicate jaws
0.9 x 0.55 mm tips
handle diameter 8 mm
overall length 140 mm
curved, without lock
curved, with lock




**BARRAQUER****NEEDLE HOLDER**

11 mm delicate jaws
 1.1 x 0.75 mm tips
 handle diameter 10 mm
 overall length 140 mm

G-17520**G-17525****G-32315****BARRAQUER****NEEDLE HOLDER**

curved, without lock
 8 mm standard jaws
 1.2 x 0.8 mm tips
 original model with longitudinal groove
 handle diameter 10 mm
 overall length 130 mm

**11****G-32320**
G-32321**BARRAQUER****NEEDLE HOLDER**

12 mm standard jaws
 1.6 x 1.0 mm tips
 handle diameter 10 mm
 overall length 140 mm

CASTROVIEJO

NEEDLE HOLDER

short model

9 mm extra delicate jaws

0.9 x 0.55 mm tips

overall length 90 mm

G-17545**G-17540**

straight, without lock

straight, with lock



CASTROVIEJO

NEEDLE HOLDER

10.5 mm delicate jaws

1.5 x 0.8 mm tips

overall length 140 mm

G-17560**G-17565**

straight, without lock

straight, with lock



CASTROVIEJO

NEEDLE HOLDER

10.5 mm delicate jaws

1.5 x 0.8 mm tips

overall length 140 mm


G-17570**G-17575**

curved, without lock

curved, with lock

**COMBINED NEEDLE HOLDER / SCISSORS**

curved 12 mm blades with
3 mm needle holder and tying platform
used to hold, place and cut fine suture material
overall length 125 mm

**G-17440**


without lock

G-17441

with lock

G-32949**HATTENBACH****HYBRID NEEDLE HOLDER**

angled, for the anterior chamber
25 gauge / 0.5 mm
overall length 90 mm



for use in combination with
hybrid tying forceps G-32947

**KALT****NEEDLE HOLDER**

standard jaws
serrated, with longitudinal groove
thumb rest and catch
overall length 135 mm

G-17610

straight, cross serrated jaws

G-17615

curved, cross serrated jaws

**ARRUGA****NEEDLE HOLDER**


14 mm heavy jaws
thumb rest and catch
overall length 140 mm
G-17640
G-17645
straight, smooth jaws
curved, smooth jaws

**G-17660****HEIDELBERG MODEL****NEEDLE HOLDER**

straight
14 mm heavy jaws
with ratchet lock
overall length 125 mm

**G-17650****BOYNTON****NEEDLE HOLDER**

straight
with ratchet lock
overall length 120 mm

**G-17665****ROHRSCHNEIDER****NEEDLE HOLDER**

straight
13 mm heavy jaws
1.3 x 1.5 mm tips
overall length 100 mm

G-19270 **BARRAQUER**
RAZOR BLADE HOLDER
 flat handle
 overall length 110 mm



G-19320 **CASTROVIEJO**
RAZOR BLADE HOLDER
 6 mm round handle
 overall length 120 mm



G-19330 **TROUTMAN**
RAZOR BLADE HOLDER
 straight
 6 mm round handle
 overall length 95 mm



G-19300 **RAZOR BLADES**
 10 pcs. per box



G-18150 **HOLDER FOR ABSORBENT STRIPS**
 serrated
 to secure sterile absorbent swabs (G-18127)
 overall length 175 mm

MANIPULATORS
CHOPPERS, FIXATION RINGS
HOOKS, SPATULAS
LENS LOOPS



FIXATION RINGS	149
NUCLEUS CHOPPERS	151
LENS & NUCLEUS MANIPULATORS	155
IRIS HOOKS	157
NUCLEUS SPATULAS	158
CAPSULE POLISHERS	160
LENS LOOPS	161
ORBITAL SPATULAS	162
CYCLODIALYSIS SPATULAS	164
IRIS SPATULAS	165
FEMTOSECOND LASER MANIPULATORS	169
STRABISMUS HOOKS	173
MEMBRANE PEELERS AND SPATULAS	175

MANIPULATORS, CHOPPERS, FIXATION RINGS,
HOOKS, SPATULAS, LENS LOOPS

12



G-37440

**HATTENBACH
GLOBE FIXATION RING**
for transconjunctival PPV
12 mm inside diameter



G-33755

**THORNTON
FIXATION RING**
16 mm inside diameter
12 atraumatic teeth



G-31490

**FINE-THORNTON
FIXATION RING**
for clear cornea and refractive procedures
16 mm inside diameter, 10 mm opening
12 atraumatic teeth



G-31491

**FINE-THORNTON
FIXATION RING**
for clear cornea and refractive procedures
13 mm inside diameter, 10 mm opening
9 atraumatic teeth



MANIPULATORS, CHOPPERS, FIXATION RINGS, HOOKS, SPATULAS, LENS LOOPS

FIXATION RINGS

12



G-33757 **FINE-THORNTON
FIXATION RING**

for clear cornea and refractive procedures
swivel handle
16 mm inside diameter, 10 mm opening
12atraumatic teeth



G-33758 **FINE-THORNTON
FIXATION RING**

for clear cornea and refractive procedures
swivel handle
13 mm inside diameter, 10 mm opening
9 atraumatic teeth



G-33759 **FINE-THORNTON
FIXATION RING**

for clear cornea and refractive procedures
swivel handle
13 mm inside diameter, 10 mm opening
9 atraumatic teeth




G-37425 **FIXATION, HOLDING AND
MEASURING INSTRUMENT**
for trocar system, radius 3.5 mm



FLIERINGA SCLERAL FIXATION RING

G-15065	15 mm
G-15070	16 mm
G-15075	17 mm
G-15080	18 mm
G-15085	19 mm
G-15090	20 mm







G-32505 THOMAS NEUHANN
NUCLEUS DIVIDING HOOK
for vertical chopping techniques
0.5 mm wide fine tip
overall length 108 mm




MIYOSHI
NUCLEUS DIVIDING HOOK
1 mm tip, 0.6 mm inner edge
overall length 104 mm
universal
G-31500
G-31501 for use in left hand
G-31502 for use in right hand



G-31504 VÖLCKER
NUCLEUS DIVIDING HOOK
for narrow pupils
with outer notch to push iris aside
1.5 mm inner edge
overall length 104 mm



G-31503 ZIRM
NUCLEUS DIVIDING HOOK
1.5 mm inner edge
overall length 104 mm



MANIPULATORS, CHOPPERS, FIXATION RINGS, HOOKS, SPATULAS, LENS LOOPS

G-32500


TOBIAS NEUHANN

CHOPPER

paddle shape nucleus divider
0.5 mm wide tip
overall length 108 mm



side view



front view

G-32282

AGARWAL

KARATE PHACO CHOPPER

0.9 mm long tip
0.6 mm inner edge
overall length 104 mm



HAEFLIGER

NUCLEUS DIVIDING HOOK

to crack, chop and manipulate lenses
1.0 mm triangular tip
0.7 mm inner sharp hook
overall length 104 mm

G-31505

G-31506




G-31507

NAGAHARA

NUCLEUS DIVIDING HOOK

for horizontal chopping techniques
1.5 mm long tip
0.9 mm inner edge
overall length 98 mm





G-32501 *THOMAS NEUHANN*
NUCLEUS ROTATOR and
TOBIAS NEUHANN
CHOPPER
 combination of
 horizontal Y-shaped nucleus rotator and
 paddle shaped nucleus divider
 overall length 119 mm



like G-32160




like G-32500



G-32168 *THOMAS NEUHANN*
SPATULA and
NAGAHARA
NUCLEUS DIVIDING HOOK
 combination of
 0.5 mm wide phaco spatula and
 horizontal phaco chopper
 overall length 116 mm



like G-16189



like G-31507



G-32161 *THOMAS NEUHANN*
NUCLEUS ROTATOR AND DIVIDING HOOK
 combination of
 horizontal Y-shaped nucleus rotator and
 vertical phaco chopper
 overall length 117 mm



like G-32160



like G-32505



G-32164 *THOMAS NEUHANN*
SPATULA AND NUCLEUS DIVIDING HOOK
 combination of
 0.5 mm wide phaco spatula and
 vertical phaco chopper
 overall length 125 mm



like G-16189




like G-32505

MANIPULATORS, CHOPPERS, FIXATION RINGS, HOOKS, SPATULAS, LENS LOOPS

G-32162

**THOMAS NEUHANN
IRIS HOOK AND NUCLEUS DIVIDING HOOK**
combination of push-pull iris hook and
NEUHANN vertical chopper
overall length 120 mm



like G-16175




like G-32505

G-32169

**T. KOHNEN
IRIS HOOK AND NUCLEUS DIVIDING HOOK**
combination of push-pull iris hook and
NAGAHARA horizontal chopper
overall length 114 mm



like G-16175



like G-31507


**G-16181
G-16182**

**LIEBERMANN
MICROFINGER**
Nucleus Manipulator
1.5 mm curved microfinger
overall length 103 mm
for use in right hand
for use in left hand



G-31727

**AKAHOSHI
NUCLEUS SUSTAINER**
to manipulate and sustain nucleus
used in combination with prechopper
0.5 mm ball tip
1.8 mm angled
overall length 100 mm



G-16185

**DRYSDALE
NUCLEUS MANIPULATOR AND
ANTERIOR CAPSULE POLISHER**

to manipulate, rotate, and divide the nucleus
1.0 mm x 1.5 mm paddle shaped tip
overall length 100 mm



G-32160

**THOMAS NEUHANN
NUCLEUS ROTATOR**

horizontal Y-shaped rotator and manipulator
0.15 mm diameter tip
overall length 104 mm



G-32034

**HENNIG / FRIEDRICH
LENS EXTRACTION HOOK
(Hefri Hook)**

overall length 105 mm



G-32550

**MAINZ MODEL
LENS MANIPULATOR**
straight
0.15 mm diameter tip
overall length 114 mm



G-32556

**MAINZ MODEL
MANIPULATOR**
vaulted
0.15 mm diameter tip
overall length 109 mm



MANIPULATORS, CHOPPERS, FIXATION RINGS, HOOKS, SPATULAS, LENS LOOPS



JACOBI
LENS MANIPULATOR
vaulted
overall length 99 mm
G-16160
for 0.3 mm IOL hole
G-16161
for 0.2 mm IOL hole



SINSKEY
LENS MANIPULATION HOOK
for 0.2 mm IOL hole
10 mm angled shaft
overall length 103 mm
G-16166
G-16167
flat handle
round handle




G-32557
BREMEN MODEL
LENS MANIPULATOR
vertical Y-hook
0.15 mm diameter tip
overall length 116 mm



G-32555
RENTSCH
LENS MANIPULATOR AND HOOK
straight
0.18 mm diameter boat hook
overall length 137 mm



G-32530
FRANKFURT MODEL
MANIPULATOR
straight
0.18 mm diameter lens hook
overall length 137 mm



G-32275

ZIRM
LENS LOOP MANIPULATOR
vertical Y-hook
0.15 mm diameter tip
overall length 108 mm



G-32276

LESTER
LENS LOOP MANIPULATOR
vaulted
0.25 mm diameter tip
overall length 102 mm



G-16169

TANDOGAN
EXPLANTATION HOOK
for IOL
0.2 mm diameter tip
overall length 103 mm




G-16210

HIRSCHMANN
LENS SPATULA
0.6 mm wide spatula with
0.2 mm notches
overall length 104 mm



G-16175

DARDENNE
IRIS HOOK AND LENS MANIPULATOR
push-pull tip
overall length 105 mm



MANIPULATORS, CHOPPERS, FIXATION RINGS, HOOKS, SPATULAS, LENS LOOPS

G-32227

LUND

IRIS HOOK

0.8 mm collar button
overall length 105 mm



Geuder

G-16205

BONN MODEL

IRIS HOOK

extra delicate
0.15 mm diameter hook
overall length 111 mm



Geuder

G-31911

FRANKFURT MODEL

IRIS HOOK AND LENS MANIPULATOR

0.65 mm collar button
overall length 135 mm




Geuder

G-16191

STEINMETZ

SPATULA

for manual lens fragmentation
1 mm wide nucleus bisector
for use in left hand
overall length 110 mm




Geuder

G-16193

SPIRIG

SPATULA


for manual lens fragmentation
0.5 mm wide nucleus bisector
overall length 110 mm



Geuder



G-32728 *RENTSCH*
SPATULA
 12 mm long spatula
 tapered from 0.6 mm to 0.3 mm
 overall length 118 mm



Geuder



G-32005 *ENGELS*
PHACO SPOON
 0.6 mm wide round spatula
 overall length 86 mm



Geuder



G-32007 *D. D. KOCH*
PHACO SPATULA
 curved with notches
 duckbill 0.5mm wide tip
 overall length 105 mm



Geuder




G-32523 *AMON*
PHACO SPATULA
 paddle-shaped tip
 0.8 mm wide x 1.6 mm long
 0.2 mm thick
 overall length 103 mm



Geuder



SPATULA
 12 mm spatula
 angled, 45°
 overall length 108 mm
G-16189 0.5 mm wide x 0.25 mm thick
G-16190 1.0 mm wide x 0.25 mm thick




Geuder

MANIPULATORS, CHOPPERS, FIXATION RINGS,
HOOKS, SPATULAS, LENS LOOPS

G-32724

**SEKUNDO
ENCLAVATION SPATULA**
for retropupillary implantation of
Iris-Claw IOL
14 mm long spatula
tapered from 0.75 mm to 0.3 mm
overall length 118m



G-31378

**ANIS
BALL CAPSULE POLISHER**
1.5 mm diameter hemisphere tip
textured
overall length 113 mm



G-32563

**RENTSCH
CAPSULE POLISHING CURETTE**
angled 135°
1.5 mm diameter tip
anterior and posterior semi-sharp edges
overall length 108 mm



G-32564

**RENTSCH
CAPSULE POLISHING CURETTE**
angled 90°
1.5 mm diameter tip
anterior and posterior semi-sharp edges
overall length 110 mm





G-15620 *ARLT*
LENS LOOP
serrated, curved
3.5 mm wide
overall length 135 mm




G-15630 *SNELLEN*
LENS LOOP
smooth, curved
3.5 mm wide
overall length 135 mm



G-15640 *WEBER*
LENS LOOP
small
4 mm wide
overall length 130 mm



G-15645 *WEBER*
LENS LOOP
large
5 mm wide
overall length 135 mm



G-15635 *SPIR/G*
LENS LOOP
for manual lens fragmentation
serrated, curved
2.2 mm wide
overall length 135 mm

MANIPULATORS, CHOPPERS, FIXATION RINGS, HOOKS, SPATULAS, LENS LOOPS

ORBITAL SPATULAS

12



BANGERTER

SPATULA

lightly curved

overall length 139 mm

G-15720

6 mm wide

G-15725

10 mm wide

G-15730

12 mm wide



G-15731

BEHRENS-BAUMANN

ORBITAL SPATULA

with aspiration

10 mm wide, lightly curved

overall length 137 mm



SAUTTER

ORBITAL SPATULA

overall length 142 mm

G-15740

22 x 13 mm curved plate

G-15745

22 x 16 mm curved plate

G-15750

22 x 19 mm curved plate



G-32445

HAMBURG MODEL

ORBITAL SPATULA

18 x 12 mm curved plate

overall length 145 mm



G-32446

ORBITAL SPATULA

18 x 8 mm curved plate
overall length 145 mm



G-15758

KERATOPLASTY

SPATULA
perforated



G-32727

BINKHORST

SYNECHIAE SPATULA
10 mm long, semi-sharp
extra delicate
overall length 129 mm



G-15477

CASTROVIEJO

SYNECHIAE SPATULA
7 mm long, semi-sharp
0.5 mm wide
overall length 86 mm

MANIPULATORS, CHOPPERS, FIXATION RINGS,
HOOKS, SPATULAS, LENS LOOPS



G-15475 *CASTROVIEJO*
CYCLODIALYSIS SPATULA
double ended, blunt
10 and 15 mm long
overall length 134 mm



G-15470 *ELSCHNIG*
CYCLODIALYSIS SPATULA
10 mm long, blunt
1.2 mm wide
overall length 127 mm



G-15480 *CASTROVIEJO*
CYCLODIALYSIS SPATULA
10 mm long, blunt
0.7 mm wide
overall length 129 mm



G-15485 *CASTROVIEJO*
CYCLODIALYSIS SPATULA
15 mm long, blunt
0.7 mm wide
overall length 129 mm



G-15400 WECKER
IRIS SPATULA

slim
1.5 mm wide, lightly curved
overall length 129 mm



G-15405 WECKER
IRIS SPATULA

wide
2.5 mm wide, lightly curved
overall length 132 mm



G-31990 WECKER
IRIS SPATULA

conical
0.9 mm wide, lightly curved
overall length 134 mm



G-15390 FRANKFURT MODEL
IRIS SPATULA

0.9 mm wide, curved
overall length 125 mm



G-15391 FRANKFURT MODEL
IRIS SPATULA

extra delicate
0.25 mm wide, curved
overall length 125 mm

MANIPULATORS, CHOPPERS, FIXATION RINGS, HOOKS, SPATULAS, LENS LOOPS

G-32729

IRIS SPATULA

extra delicate, conical
7 mm long spatula, angled
tapers from 0.5 mm to 0.3 mm
overall length 114 mm



G-31975

G-31980

G-31985

BANGERTER

IRIS SPATULA

round handle
8 mm long spatula, angled
overall length 118 mm

1.0 mm wide

1.5 mm wide

2.0 mm wide



BANGERTER

IRIS SPATULA

flat handle
8 mm long spatula, angled
overall length 125 mm

1.0 mm wide

1.5 mm wide

2.0 mm wide



G-15360

G-15365

G-15370

G-15755

ESSEN MODEL

SPATULA

10 mm wide
25 mm angled, flat
overall length 149 mm





G-15757 *GIESSEN MODEL*
SPATULA AND RETRCTOR
9 mm wide
angled, concave
overall length 162 mm



G-15780 *WELLS*
ENUCLEATION SPOON
overall length 147 mm



G-33912 *ZIRM*
LASIK CORNEAL PROTECTOR
7.0 x 2.5 mm corneal shield
overall length 120 mm



G-33930 *THOMAS NEUHANN*
LASIK SPATULA
with beveled blunt tip
rounded anterior and
flat posterior surface
for lifting the corneal flap
overall length 106 mm

MANIPULATORS, CHOPPERS, FIXATION RINGS, HOOKS, SPATULAS, LENS LOOPS

G-33935

RAU

LASEK SPATULA

5.5 mm long spatula
for repositing and smoothing
the epithelial flap after laser ablation
overall length 116 mm



G-38607

DESCEMET INCISION HOOK

0.2 mm diameter blunt tip
for marking the area of Descemetorhexis
in recipient's eye
overall length 103 mm



G-38608

DESCEMET SCRAPER

1 mm semi-circular tip
for scraping the recipient Descemet's
membrane
overall length 103 mm



G-38622

SZURMAN

DMEK PREPARATION SPATULA

disc-shaped spatula
2.5 mm diameter, blunt
for removal of Descemet's membrane from graft
overall length 121 mm



G-38604

THOMAS NEUHANN

ENDOTHELIAL INSERTION SPATULA

with teflon insert, angled 145°, 4 x 2 mm
for easy insertion of donor lamellar button
with notched tip for "pull-through" technique
overall length 104 mm
(used with G-38606)





G-33966

DICK

FEMTOCATARACT SPATULA

double ended with round semi-sharp tips
for easy opening of femto incisions
spatula 2.5 mm flattened, 0.65 mm wide
for main and sideport incisions
knife 90° angled, 0.9 mm wide
for arcuate incisions
overall length 100 mm



G-33967

BREYER

FEMTOCATARACT SPATULA

angled
1 mm blunt round tip
for opening femto incisions
overall length 100 mm




G-33968

HOFFMANN

**FEMTOCATARACT NUCLEUS ROTATOR
AND CHOPPER**

horizontal Y-shaped rotator, 0.75 mm wide
and 1 mm "Quick Chop" tip
overall length 120 mm





G-33950

CHANSUE

FEMTO DOUBLE INSTRUMENT

blunt spoon-shaped lenticule spatula
1.25 mm wide, vaulted and
hook, 2 mm long, angled
overall length 115 mm
to open the incision, delineate the edge
of the lenticules and separate the interfaces
during SMILE procedures



MANIPULATORS, CHOPPERS, FIXATION RINGS, HOOKS, SPATULAS, LENS LOOPS


G-33952

CHANSUE

FEMTO DOUBLE INSTRUMENT

conical lenticule spatula, 9 mm long, vaulted hook, 2 mm long, angled overall length 115 mm

to open the incision, delineate the edge of the lenticules and separate the interfaces during SMILE procedures




G-33954

GUELL

FEMTO DOUBLE INSTRUMENT

1.25 mm wide spoon-shaped spatula, straight 2 mm long hook, angled overall length 113 mm

to open the incision, delineate the edge of the lenticules and separate the interfaces during SMILE procedures



G-33940

BLUM

FEMTO-LASIK ROUND SPOON

spoon-shaped spatula, 1.25 mm wide semi-sharp, bevel up, vaulted overall length 105 mm

to undermine, lift and reposit LASIK flaps created with a femtosecond laser




G-33942

BREYER

FEMTO-LASIK DOUBLE-ENDED INSTRUMENT

combination of
BLUM spoon-shaped spatula, semi-sharp and hook with 1 mm inner edge (one blunt and one slightly sharp edge)
overall length 115 mm






G-33427 PFAEFFL

FEMTO FLAP LIFTER


double ended, both ends blunt
10.5 mm conical v-shaped hook, curved and
1.5 mm modified SINSKEY hook
overall length 115 mm
for finding the flap edge, lifting, and retreating it



G-33943 PFAEFFL

FEMTO FLAP LIFTER



double ended
10.5 mm conical v-shaped hook, straight and
1.5 mm modified SINSKEY hook
overall length 115 mm
for finding the flap edge, lifting and retreating it



G-33958 PFAEFFL-BREYER

FEMTO FLAP LIFTER

double ended
10.5 mm conical v-shaped hook, curved
with double bevel tip and
1.5 mm modified SINSKEY hook
overall length 115 mm
for finding the flap edge, lifting and retreating it



G-33959 BREYER

FEMTO-LASIK HOOK INSTRUMENT

angled
1.3 mm wide scraper hook
to delineate and remove the lenticule
overall length 105 mm



MANIPULATORS, CHOPPERS, FIXATION RINGS,
HOOKS, SPATULAS, LENS LOOPS




G-33944

KOHNEN

FEMTO-LASIK DOUBLE-ENDED INSTRUMENT

11 mm long x 0.5 mm wide spatula and
2 mm modified SINSKEY hook, blunt tip
overall length 120 mm



G-37450

BONN MODEL

SUTURE TOOL

for scleral sutures



G-19137

CUEPPERS

STRABISMUS MUSCLE TUCKER

15 mm scale, for muscle stretching
serrated J-shaped tendon hook to
prevent the muscle from slipping




G-32726

MUENSTER MODEL

SUTURE HOOK

hook 5 mm

overall length 149 mm





G-15799 *GRAEFE*
STRABISMUS HOOK
short, 5 mm
overall length 130 mm



G-15800 *GRAEFE*
STRABISMUS HOOK
small, 8 mm
overall length 131 mm



G-15805 *GRAEFE*
STRABISMUS HOOK
medium, 10 mm
overall length 133 mm



G-15810 *GRAEFE*
STRABISMUS HOOK
large, 12 mm
overall length 137 mm



G-32915 *HAMBURG MODEL*
STRABISMUS HOOK
10 mm long
slim 1 mm diameter hook
overall length 135 mm

MANIPULATORS, CHOPPERS, FIXATION RINGS,
HOOKS, SPATULAS, LENS LOOPS



G-15830 *MALBRAN*
STRABISMUS HOOK
9.5 mm serrated hook with right-angle knob
to prevent muscle slippage
overall length 140 mm



G-15820 *BONN MODEL*
STRABISMUS HOOK
11 mm
overall length 143 mm



G-15821 *BONN MODEL*
STRABISMUS HOOK
11 mm retinal detachment hook
with 1.0 mm hole
overall length 143 mm




G-32905 *BONN MODEL*
STRABISMUS HOOK
curved, tenotomy hook
overall length 145 mm




G-32920 *JAMESON*
STRABISMUS HOOK
11 mm, with bulbous tip
overall length 125 mm




G-37435 *HATTENBACH
HORSESHOE-HOOK*
to pull the trocars off the lance
and to remove the trocars from the sclera
overall length 98 mm




G-37551 *BARTZ-SCHMIDT
ILM-SCRAPER*
20 gauge / 0.9 mm
diamond dusted tip



MEMBRANE SPATULAS



angled at 1.5 mm
G-37530 20 gauge / 0.9 mm




angled at 2 mm and 4 mm
G-37535 20 gauge / 0.9 mm



subretinal pick, angled at 4 mm
G-37545 20 gauge / 0.9 mm



LUCKE
curved
G-37515 20 gauge / 0.9 mm




2.5 mm angled
G-37525 20 gauge / 0.9 mm



angled at 1 mm, long rakes
G-37570 20 gauge / 0.9 mm

MEMBRANE PEELERS



1.0 mm angled
G-37500 20 gauge / 0.9 mm
G-37503 23 gauge / 0.6 mm



MORRIS-KROLL
1.5 angled
G-37505 20 gauge / 0.9 mm



0.65 mm angled
black chromium
G-37506 20 gauge / 0.9 mm



0.4 mm, recessed, angled
black chromium
G-37511 20 gauge / 0.9 mm
G-37512 23 gauge / 0.6 mm
G-37513 25 gauge / 0.5 mm

MANIPULATORS, CHOPPERS, FIXATION RINGS,
HOOKS, SPATULAS, LENS LOOPS

13

SCLERAL DEPRESSORS LOCALIZERS





G-32152 LAQUA

SCLERAL DEPRESSOR

with pocket clip

3 x 1.8 mm teardrop-shaped tip

overall length 128 mm



G-32150 LAQUA

SCLERAL DEPRESSOR

double-ended with pocket clip

3 x 1.8 mm teardrop-shaped tip

2 mm diameter scleral marker

overall length 134 mm



G-32153 LUCKE

SCLERAL DEPRESSOR

with pocket clip

3.8 x 2.0 mm teardrop-shaped tip

overall length 128 mm



G-32027 KREISSIG

SCLERAL DEPRESSOR

double-ended with pocket clip

conical tip with 2.4 mm diameter

and 6.0 mm bar

overall length 128 mm

SCLERAL DEPRESSORS, LOCALIZERS

G-15710 MEYER-SCHWICKERATH
LOCALIZER AND SCLERAL DEPRESSOR
for fixation and rotation




G-32715 SCHOCKET
SCLERAL DEPRESSOR
double-ended with pocket clip
3 x 2 mm teardrop tip and
6.0 mm bar
overall length 132 mm



G-32714 SCLERAL DEPRESSOR
double-ended
2 mm and 4 mm diameter
diamond dusted tips
overall length 120 mm



G-37445 HATTENBACH
TRIANGLE SCLERAL DEPRESSOR
for gentle buckling and firm grip without the risk of rupture
overall length 133 mm



PROBES DILATORS



G-150301 **WILDER****LACRIMAL PROBE**

delicate, long taper
0.45 mm tip diameter
overall length 104 mm

**G-150302** **WILDER****LACRIMAL PROBE**

medium, medium taper
0.50 mm tip diameter
overall length 104 mm

**G-150303** **WILDER****LACRIMAL PROBE**

strong, short taper
0.55 mm tip diameter
overall length 104 mm

**G-150401** **WILDER****LACRIMAL PROBE**

slim
delicate, long taper
0.30 mm tip diameter
overall length 95 mm

**G-150402** **WILDER****LACRIMAL PROBE**

slim
medium, medium taper
0.35 mm tip diameter
overall length 95 mm


**G-150403** **WILDER****LACRIMAL PROBE**

slim
strong, short taper
0.40 mm tip diameter
overall length 95 mm

**G-15020** **GIESSEN MODEL****LACRIMAL PROBE**

11 mm tapered shaft, angled
0.4 mm tip width
overall length 77 mm





*BOWMAN***LACRIMAL PROBE**

cylindrical, silver

overall length 130 mm

G-150001 size 0000/000**G-150002** size 00/0**G-150003** size 1/2**G-150004** size 3/4**G-150005** size 5/6*BOWMAN***LACRIMAL PROBE**

cylindrical, stainless steel

overall length 130 mm

G-150101 size 0000/000**G-150102** size 00/0**G-150103** size 1/2**G-150104** size 3/4**G-150105** size 5/6*BOWMAN***LACRIMAL PROBE**

bulbous tips, silver

overall length 130 mm

G-150051 size 0000/000**G-150052** size 00/0**G-150053** size 1/2**G-150054** size 3/4**G-150055** size 5/6*BOWMAN***LACRIMAL PROBE**

bulbous tips, stainless steel

overall length 130 mm

G-150151 size 0000/000**G-150152** size 00/0**G-150153** size 1/2**G-150154** size 3/4**G-150155** size 5/6

JUENEMANN**LACRIMAL PROBE**

with 72 mm probe and 18 mm grip plate
to insert suture material

provided with cleaning adaptor G-31022
G-33193 0.6 mm diameter
G-33194 0.8 mm diameter
G-33195 1.0 mm diameter

**G-15315****KRUMEICH****DALK PROBE**

for deep lamellar keratoplasty
12 mm tapered probe
overall length 73 mm

**G-15105**
G-15110**MACKENSEN****TRABECULAR PROBE**

7 mm curved probe
overall length 92 mm
G-15105 right
G-15110 left

**G-15125**
G-15130**HARMS****TRABECULAR PROBE**

dual 9.5 mm curved probes
3 mm spread
overall length 46 mm
G-15125 right
G-15130 left

**G-15136**
G-15138**THOMAS NEUHANN****U-SHAPED PROBE**

with holding plate
dual 9 mm curved probes
3 mm spread
overall length 12 mm
G-15136 0.30 mm probe diameter
G-15138 0.25 mm probe diameter

**G-15137****THOMAS NEUHANN****HOLDER**

with lock
for U-shaped probes G-15136 and G-15138
overall length 92 mm




CANNULAS
I/A INSTRUMENTS


LACRIMAL	189
ANTERIOR CHAMBER MAINTAINERS	191
HYDRODISSECTION	192
IRRIGATION VECTIS	193
VISCOELASTIC	194
CAPSULE POLISHERS	195
LASIK	196
BIMANUAL IRRIGATION/ASPIRATION	198
COAXIAL IRRIGATION/ASPIRATION	203
BACKFLUSH	206
INFUSION	209

LACRIMAL CANNULA


conical, straight
tapers from 19 to 25 gauge
25 mm shaft

G-33140silver
stainless steel**LACRIMAL CANNULA**


conical, strongly curved
tapers from 19 to 25 gauge
25 mm shaft

G-15140silver
stainless steel**LACRIMAL CANNULA**


conical, angled at 10 mm
tapers from 19 to 25 gauge
25 mm shaft

G-15150silver
stainless steel**LACRIMAL CANNULA**


conical, angled at 18 mm
tapers from 19 to 25 gauge
25 mm shaft

G-15160silver
stainless steel**G-15190****ANEL****LACRIMAL CANNULA**

straight
15 mm shaft
23 gauge / 0.64 mm

**G-15195****ANEL****LACRIMAL CANNULA**


strongly curved
16 mm shaft
23 gauge / 0.64 mm



LACRIMAL CANNULA	
cylindrical, straight	
G-33160	20 gauge / 0.9 x 35 mm
G-33170	no. 12, 0.72 x 32 mm
G-33175	no. 14, 0.67 x 30 mm
G-33180	no. 16, 0.60 x 22 mm
G-33185	no. 18, 0.45 x 22 mm
G-33190	no. 20, 0.42 x 22 mm
G-15170	<i>BONN MODEL</i>
LACRIMAL CANNULA	
with probe tip	
1.15 mm diameter olive tip	
BANGERTER	
LACRIMAL PROBE CANNULA	
straight, 0.6 mm side port	
G-15271	23 gauge / 0.64 x 60 mm
G-15272	21 gauge / 0.80 x 40 mm
G-15273	23 gauge / 0.64 x 40 mm
G-33197	<i>JUENEMANN</i>
LACRIMAL PROBE CANNULA	
to insert suture material	
21 gauge / 0.8 x 35 mm	
G-33198	<i>JUENEMANN</i>
LACRIMAL PROBE CANNULA	
with grip plate	
65 mm probe	
to insert suture material	
21 gauge / 0.8 x 65 mm	

RETROBULBAR INJECTION CANNULA


- G-15375 24 gauge / 0.6 x 60 mm
 G-32605 24 gauge / 0.6 x 35 mm
 G-32600 25 gauge / 0.5 x 35 mm
 G-15376 25 gauge / 0.5 x 50 mm
 G-15377 25 gauge / 0.5 x 40 mm





**G-15200****CRANZ
EYE DOUCHE**irrigation through five front ports
silver**G-33416***BOPP***ANTERIOR CHAMBER INFUSION CANNULA**anterior chamber maintainer with
threaded tip for better fixation
and silicone tube
20 gauge / 0.9 mm**G-33422***BONN MODEL***ANTERIOR CHAMBER INFUSION CANNULA**self-retaining textured tip
dual 0.35 mm side ports
silicone tube 15 cm
23 gauge / 0.6 mm**G-16105****IRRIGATION CYSTOTOME CANNULA**

26 gauge / 0.45 mm

**G-15220****ENZYME CANNULA**


24 gauge / 0.6 mm




G-15275	CYCLODIALYSIS SPATULA CANNULA 12 mm long spatula anterior 0.2 mm port 26 gauge / 0.45 mm			
G-31355	BUSSE CYCLODIALYSIS SPATULA CANNULA 12 mm long spatula dual 0.2 mm port 26 gauge / 0.45 mm			
G-15250 G-32166 G-32167	SAUTTER HYDRODISSECTION CANNULA front opening flattened tip vertical, 27 gauge / 0.4 mm vertical, 26 gauge / 0.45 mm horizontal, 27 gauge / 0.4 mm	 vertical G-15250 G-32166	 horizontal G-32167	
G-32155	TUEBINGEN MODEL HYDRODISSECTION CANNULA 30 gauge / 0.3 mm			
G-15185	FRANKFURT MODEL HYDRODISSECTION CANNULA 27 gauge / 0.4 mm			
G-31895	HYDRODISSECTION CANNULA vertically flattened tip front opening 24 gauge / 0.6 x 20 mm			
G-15252 G-15253	SPIRIG HYDRODISSECTION CANNULA for manual lens fragmentation right left			

G-16115 IRRIGATION VECTIS CANNULA


4 front ports
5 mm wide lens loop
curved

**G-32759 IRRIGATION VECTIS CANNULA**


small model
4 front ports
4 mm wide lens loop
curved

**G-31510 KLÖCKNER EXPRESSOR CANNULA**

8 mm wide lens loop
4 front ports
lightly curved

**G-31807 HEALON - VISCOELASTIC CANNULA**

flattened tip, angled at 6 mm
30 mm shaft
21 gauge / 0.8 mm

**G-22619 PH. JACOBI TRABECULAR ASPIRATION CANNULA**

angled at 13 mm
24 gauge / 0.6 mm

**BURK IRRIGATION TRABECULOTOMY CANNULA**





8 mm curved irrigation probe
30 gauge / 0.3 mm

G-15188
G-15189

**G-33199 NAGATA SUB-TENON INJECTION CANNULA**


foratraumatic penetration and
injection of anesthetic agent
24 gauge / 0.6 mm



G-15180	IRRIGATION CANNULA bulbous tip, 0.75 mm diameter tip extra delicate 26 gauge / 0.45 mm		
G-15230 G-15240	BONN MODEL IRRIGATION CANNULA flattened tip, front opening delicate, 0.5 x 1.5 mm medium, 0.6 x 1.9 mm		
G-16134 G-16135 G-16136	BINKHORST IRRIGATION CANNULA U-shaped 22 gauge / 0.7 mm straight angled right angled left		
G-32760 G-32765	MUELLER IRRIGATION CANNULA J-shaped with probe tip 26 gauge / 0.45 x 19 mm angled left angled right		
G-16138	RENTSCH IRRIGATION CANNULA J-shaped angled at 10 mm 26 gauge / 0.45 x 24 mm		
G-16106	SIMCOE IRRIGATION / ASPIRATION CANNULA dual shaft 0.3 mm diameter aspiration port front opening irrigation port 23 gauge / 0.64 mm with silicone tube and Luer-Lock hub		

G-16120 IRRIGATION / ASPIRATION CANNULA**COAXIAL**

straight
aspiration port 0.4 mm
infusion tube 18 gauge / 1.2 mm

**ASPIRATION CANNULA**


lightly curved, front opening
delicate, 19 gauge / 1.1 mm
standard, 17 gauge / 1.5 mm

**G-16125 KRATZ CAPSULE POLISHING CANNULA**


angled at 2 mm
half-textured tip
23 gauge / 0.64 mm

**G-16133 TEPING ASPIRATION POLISHING CANNULA**


0.3 mm posterior port
half-textured tip, curved
23 gauge / 0.64 mm

**G-16130 CAPSULE POLISHING CANNULA**

with half-textured ball tip
18 gauge / 1.2 mm


**G-16132 PAPST CAPSULE POLISHING CANNULA**

0.75 mm hemisphere textured tip
25 gauge / 0.5 mm

**KYU CAPSULE POLISHING CANNULA**


2 mm disc-shaped tip
textured

G-32529 22 gauge / 0.7 mm
G-32527 23 gauge / 0.64 mm




G-32528

GREITE
CAPSULE POLISHING CANNULA
 curette tip
 textured
 22 gauge / 0.7 mm

**G-32774**

BRAUWEILER
ASPIRATION CANNULA
 angled at 13 mm
 textured polisher tip
 0.3 mm beveled port
 21 gauge / 0.8 mm

**G-32769**

H.-R. KOCH
ASPIRATION CANNULA
 angled at 14 mm
 0.3 mm anterior port
 half-textured, spatulated tip
 21 gauge / 0.8 mm

**G-32776**


BRAUWEILER
IRRIGATION CANNULA
 angled at 12 mm
 front opening
 21 gauge / 0.8 mm

**G-32772**

BRAUWEILER
IRRIGATION CANNULA
 angled at 13 mm
 textured 0.8 mm iris retractor and
 capsule polisher
 21 gauge / 0.8 mm

**G-33909**

W. BAUER
LASIK CANNULA
 double-armed
 6 ports
 for irrigation of the flap and stromal bed
 25 gauge / 0.5 mm




G-33910

ZIRM

LASIK DISSECTOR CANNULA

flattened
front opening
25 gauge / 0.5 mm

**G-33911**

ZIRM

LASIK IRRIGATION CANNULA

flattened and textured tip
with 5 irrigation ports
20 gauge / 0.9 mm

**G-15310**

KRUMEICH

CORNEAL IRRIGATION CANNULA

for big bubble technique during DALK
semi-sharp dissector tip
0.2 mm posterior port
27 gauge / 0.4 mm



KRUMEICH

CORNEAL IRRIGATION CANNULA

for big bubble technique during DALK
angled at 12 mm

**G-15312**

0.4 mm anterior port

G-15314

0.4 mm posterior port

**G-22550****IRRIGATION CANNULA HANDPIECE**

overall length 115 mm

**G-32015****IRRIGATION CANNULA HANDPIECE**

male / female connectors
overall length 112 mm

- G-32771** **BRAUWEILER**
IRRIGATION CANNULA HANDPIECE
 male / female connectors
 78 mm



- G-32773** **BRAUWEILER**
ASPIRATION CANNULA HANDPIECE
 male / male connectors
 80 mm



- G-32778** **KEERL**
IRRIGATION CANNULA HANDPIECE
 with silicone tube



- G-32010** **IRRIGATION CANNULA HANDPIECE**
 with silicone tube



SCHMACK
 21 gauge / 0.8 mm




- G-22632** **IRRIGATION HANDPIECE**
 0.5 mm port
G-22633 **ASPIRATION HANDPIECE**
 0.3 mm port



D.D. KOCH - T. KOHNEN
ASPIRATION HANDPIECE
 0.3 mm aspiration port
 textured tip

- G-22100** 22 gauge / 0.7 mm
G-22102 23 gauge / 0.6 mm




G-22101

D.D. KOCH - T. KOHNEN
IRRIGATION HANDPIECE
 front opening
 22 gauge / 0.7 mm




G-22105
G-22106


ASPIRATION HANDPIECE
 0.25 mm aspiration port
 conical tip
 oval tube 0.95 x 0.55 mm




IRRIGATION HANDPIECE
 0.5 mm irrigation port
 conical tip
 oval tube 0.95 x 0.55 mm

**G-22637**


SCHMACK
ASPIRATION HANDPIECE
 4 x 0.15 mm posterior ports
 textured tip
 23 gauge / 0.6 mm




G-32011 **ASPIRATION HANDPIECE**
0.3 mm aspiration port
textured, conical tip
21 gauge / 0.8 mm



G-32012 **IRRIGATION HANDPIECE**
dual 0.5 mm irrigation ports
textured, conical tip
21 gauge / 0.8 mm



G-32018 **BRAUWEILER**
ASPIRATION HANDPIECE
angled at 13 mm
0.3 mm beveled port
textured tip
21 gauge / 0.8 mm



G-32019 **H.-R. KOCH**
ASPIRATION HANDPIECE
angled at 14 mm
0.3 mm anterior port
half-textured, spatulated tip
21 gauge / 0.8 mm



**G-32016**


BRAUWEILER
IRRIGATION HANDPIECE
 angled at 12 mm
 front opening
 21 gauge / 0.8 mm

**G-32017**

BRAUWEILER
IRRIGATION HANDPIECE
 angled at 13 mm
 textured 0.8 mm iris retractor and
 capsule polisher
 21 gauge / 0.8 mm

**G-32026**

R. KOCH
ASPIRATION HANDPIECE
 0.3 mm aspiration port
 textured tip
 19 gauge / 1.1 mm

**G-32028**

R. KOCH
IRRIGATION HANDPIECE
 angled
 dual 0.75 mm side ports
 19 gauge / 1.1 mm



G-32014*R. KOCH***IRRIGATION CHOPPER**

angled

dual 0.75 mm side ports

1 mm angled chopper tip

19 gauge / 1.1 mm

**G-32059***DICK***IRRIGATION CHOPPER**

angled

dual 0.75 mm side ports

1 mm curved chopper tip

19 gauge / 1.1 mm

**G-38601***ALTHAUS / CARTSBURG***IRRIGATION DESCemet HOOK**

hook diameter: 0.25 mm

hook length: 0.6 mm

for Descemtorhexis

rhombic 0.9 x 0.6 mm irrigation tube

overall length 98 mm

**G-38602***THOMAS NEUHANN***IRRIGATION DESCemet SCRAPER**

scraper width: 0.8 mm

scraper length: 0.6 mm

rhombic 0.9 x 0.6 mm irrigation tube


for peripheral scraping of the recipient bed

overall length 94 mm



G-38631


SZURMAN
DMEK IRRIGATION HANDPIECE
 curved, closed tip
 side port 0.3 mm
 tube 23 gauge / 0.6 mm

**G-38632**

SZURMAN
DMEK IRRIGATION HANDPIECE
 curved, front opening
 tube 23 gauge / 0.6 mm


G-24318
G-24418


PEP 2.2 MM I/A HANDPIECE
 Pure Efficiency Phaco
 for 2.2 mm incisions
 0.3 mm aspiration port
 coaxial I/A straight tip
 coaxial I/A 30° angled tip



with G-24311 PEP infusion sleeve


G-24311

PEP 2.2 MM INFUSION SLEEVE
 Pure Efficiency Phaco
 for 2.2 mm incisions
 purple translucent silicone



**PEP 2.4 MM I/A HANDPIECE**


Pure Efficiency Phaco
for 2.4 mm incisions
0.3 mm aspiration port
G-24332
G-24432
coaxial I/A straight tip
coaxial I/A 30° angled tip




with G-24314 PEP infusion sleeve

G-24314 PEP 2.4 MM INFUSION SLEEVE


Pure Efficiency Phaco
for 2.4 mm incisions
green translucent silicone

**G-24330 I/A HANDPIECE**






coaxial straight tip
0.3 mm aspiration port
with G-24100 infusion sleeve

**G-24430 I/A HANDPIECE**

coaxial 30° angled tip
G-24430
G-24435
G-24440
0.3 mm aspiration port
0.35 mm aspiration port
0.4 mm aspiration port



with G-24100 infusion sleeve

G-24100**INFUSION SLEEVE**blue silicone
for I/A handpieces**G-22525***PAPST***I/A HANDPIECE**coaxial
infusion tube 16 gauge / 1.6 mm
aspiration port 0.3 mm, straight
aspiration port 0.4 mm, straight
aspiration port 0.3 mm, angled 30°
aspiration port 0.4 mm, angled 30°
aspiration port 0.5 mm, angled 30°**G-22527****G-22540****G-22542****G-22545****G-22615***PAPST***I/A HANDPIECE**coaxial 30° angled tip
aspiration port 0.3 mm
infusion tube 18 gauge / 1.3 mm**G-22556***HUNOLD***I/A HANDPIECE**coaxial, curved
aspiration port 0.3 mm
infusion tube 18 gauge / 1.3 mm



for backflush tips G-37010 up to G-37080

SILICONE CHAMBER
for backflush handpieces
G-31805, G-37000, G-37002



G-31806 reusable
G-34806 single use, 10 pcs. per box, sterile

G-37005 **BACKFLUSH HANDPIECE**
with large silicone chamber
rotatable chuck



for backflush cannulas G-37010 up to G-37080

G-37006 **SILICONE CHAMBER**
large, for backflush handpiece G-37005




BACKFLUSH CANNULA

35 mm shaft

G-33050 20 gauge / 0.9 mm

G-33054 23 gauge / 0.6 mm

G-33053 25 gauge / 0.5 mm



G-33070


BACKFLUSH CANNULA

35 mm shaft

20 gauge / 0.9 mm shaft

25 gauge / 0.5 mm tip

3.5 mm tapered tip




G-33071

BACKFLUSH CANNULA

35 mm shaft

0.5 mm side port

20 gauge / 0.9 mm




G-33072

*LAQUA***BACKFLUSH CANNULA**

35 mm shaft

inside diameter: 21 gauge / 0.8 mm

outside diameter: 19 gauge / 1.1 mm




G-33073

KAMPIK**BACKFLUSH CANNULA**

35 mm shaft

inside diameter: 22 gauge / 0.7 mm

outside diameter: 20 gauge / 0.9 mm




G-33075

BACKFLUSH CANNULA

35 mm shaft

0.8 mm silicone brush

20 gauge / 0.9 mm



G-33077

BACKFLUSH CANNULA

20 gauge / 0.9 mm shaft


24 gauge / 0.55 mm tip

7 mm tapered tip, angled at 3.5 mm



BACKFLUSH TIP


35 mm shaft
G-37010 20 gauge / 0.9 mm
G-37011 23 gauge / 0.6 mm
G-37012 25 gauge / 0.5 mm

**G-37015 BACKFLUSH TIP**


35 mm shaft
 angled at 3 mm
 20 gauge / 0.9 mm

**G-37020 BACKFLUSH TIP**


35 mm shaft
 0.8 mm silicone brush
 20 gauge / 0.9 mm

**G-37025 BACKFLUSH TIP**


35 mm shaft
 angled at 6 mm
 0.8 mm silicone brush
 20 gauge / 0.9 mm

**G-37035 BACKFLUSH TIP**

35 mm shaft
 0.8 mm silicone tip
 20 gauge / 0.9 mm

**G-37046 KLEIN BACKFLUSH TIP**

20 gauge / 0.9 mm shaft
 26 gauge / 0.45 mm tip
 2 mm tapered tip
 overall length 35 mm


**G-37065 BACKFLUSH TIP**

35 mm shaft
 20 gauge / 0.9 mm shaft
 30 gauge / 0.3 mm tip
 tapered tip, angled at 5 mm
 oval aspiration port




G-37080**BACKFLUSH TIP**

35 mm shaft
 20 gauge / 0.9 mm shaft
 23 gauge / 0.6 mm tip
 2 mm tapered tip
 overall length 35 mm

**G-33420****INFUSION CANNULA**

silicone tube 40 cm
 Luer-Lock metal adapter
 20 gauge / 0.9 mm
 beveled, 4.0 mm (single-use G-33463)

**G-33425**



beveled, 5.0 mm

G-33430

beveled, 6.0 mm

INFUSION CANNULA

silicone tube 40 cm
 Luer-Lock plastic adapter
 20 gauge / 0.9 mm
 beveled, 4.0 mm

**G-33450**


beveled, 5.0 mm

G-33451

beveled, 6.0 mm

INJECTION CANNULA

for viscous fluids, without silicone tube
 for butterfly-19 cannula
 20 gauge / 0.9 mm

**G-33484**

beveled, 4.0 mm

G-33485


beveled, 5.0 mm

G-33486

beveled, 6.0 mm

INJECTION CANNULA

for viscous fluids, without silicone tube
 for butterfly-21 cannula
 20 gauge / 0.9 mm
 beveled, 4.0 mm

**G-33494**


beveled, 5.0 mm

G-33495

beveled, 6.0 mm

DUISBURG MODEL**INJECTION CANNULA**

for viscous fluids
 25 cm silicone tube with 2 metal sleeves
 and Luer-Lock adapter
 20 gauge / 0.9 mm
 beveled, 4.0 mm

**G-32180**

beveled, 5.0 mm

G-32181

beveled, 6.0 mm

G-32182

beveled, 6.0 mm

*HAMBURG MODEL***INJECTION CANNULA**


for viscous fluids

25 cm silicone tube with metal sleeve

and Luer-Lock adapter

3 spare silicone tubes

- G-33470** 20 gauge / 0.9 mm, beveled, 4 mm
G-33471 20 gauge / 0.9 mm, beveled, 5 mm
G-33472 20 gauge / 0.9 mm, beveled, 6 mm
G-33473 23 gauge / 0.6 mm, beveled, 4 mm
G-33474 23 gauge / 0.6 mm, beveled, 6 mm

**G-33475****SPARE SILICONE TUBE**

for injection cannulas

25 cm

3 pcs. per box




- G-33051**
G-33052

*WEBER***ASPIRATION CANNULA**

for viscous fluids


inside diameter: 22 gauge / 0.7 mm

outside diameter: 20 gauge / 0.9 mm

**G-37390****ASPIRATION CANNULA**

23 gauge / 0.6 mm


10 mm shaft

**G-33056***ROIDER***ASPIRATION CANNULA**

for viscous fluids

0.7 mm side port


19 gauge / 1.0 mm

**G-33055****PFCL CANNULA**

with dual shaft

24 gauge / 0.55 mm

25 gauge / 0.5 mm, lateral tube


**G-33057***CHANG***PFCL CANNULA**

for injection of heavy fluids

dual bore, coaxial

25 gauge / 0.5 mm tip


20 gauge / 0.9 mm shaft



HEIDELBERG MODEL**CANNULA**

for injection or aspiration of viscous fluids
and Densiron® 68, bevel 30°

- G-32699** 19 gauge / 1.1 mm
G-32698 18 gauge / 1.2 mm


**STOPPER**

for viscous fluid aspiration
with tube connection for single-use syringe

- G-33065** 10 ml syringe
G-33066 20 ml syringe

**PRESSURE TUBE**

for injection of viscous fluids
Luer-Lock female / male
reusable
single use, 10 pcs. per box, sterile



- G-28752** **SYRINGE HOLDER**
for manual injection of
viscous fluids in glass syringes,
with clamp and retraction mechanism

RETINAL

SYRINGE CLIP
customized labeling
G-19960
G-19961
for 2 ml syringes
for 5 ml syringes



SILICONE BULB
with Luer Cone
G-32665
G-32670
20 ml
60 ml



SINGLE USE INSTRUMENTS



NANOEDGE TREPHINES AND PUNCHES	215
NANOEDGE SURGICAL KNIVES	217
UNO COLORLINE (SCISSORS · FORCEPS · PROBES · TROCAR · FIBER OPTICS · VITREOUS CUTTERS)	218
SPECULUM	223
DMEK CARTRIDGE	223
IRIS RETRACTOR	223
LACRIMAL INTUBATION	223
CANNULAS	224
BIMANUAL HANDPIECE SET	227
BACKFLUSH	227
INFUSION / INJECTION CANNULAS	228
SCLERAL BUCKLING COMPONENTS	232

NANOEDGE TREPHINES

length 16 mm
sterile, single-use

G-520400	4.00 mm	G-520825	8.25 mm	G-521350	13.50 mm
G-520500	5.00 mm	G-520850	8.50 mm	G-521400	14.00 mm
G-520550	5.50 mm	G-520875	8.75 mm	G-521450	14.50 mm
G-520600	6.00 mm	G-520900	9.00 mm	G-521500	15.00 mm
G-520625	6.25 mm	G-520950	9.50 mm	G-521550	15.50 mm
G-520650	6.50 mm	G-521000	10.00 mm	G-521600	16.00 mm
G-520675	6.75 mm	G-521050	10.50 mm	G-521650	16.50 mm
G-520700	7.00 mm	G-521100	11.00 mm	G-521700	17.00 mm
G-520725	7.25 mm	G-521150	11.50 mm	G-521800	18.00 mm
G-520750	7.50 mm	G-521200	12.00 mm		
G-520775	7.75 mm	G-521250	12.50 mm		
G-520800	8.00 mm	G-521300	13.00 mm		



G-500001 NANOEDGE SINGLE-USE CORNEAL CUTTING BLOCK

Teflon, sterile



NANOEDGE TREPHINE PUNCHES

sterile, single-use


G-540600	6.00 mm	G-540725	7.25 mm	G-540850	8.50 mm
G-540625	6.25 mm	G-540750	7.50 mm	G-540875	8.75 mm
G-540650	6.50 mm	G-540775	7.75 mm	G-540900	9.00 mm
G-540675	6.75 mm	G-540800	8.00 mm	G-540925	9.25 mm
G-540700	7.00 mm	G-540825	8.25 mm	G-540950	9.50 mm



NANOEDGE VACUUM TREPHINE PUNCHES

sterile, single-use


G-560600	6.00 mm	G-560725	7.25 mm	G-560850	8.50 mm
G-560625	6.25 mm	G-560750	7.50 mm	G-560875	8.75 mm
G-560650	6.50 mm	G-560775	7.75 mm	G-560900	9.00 mm
G-560675	6.75 mm	G-560800	8.00 mm	G-560925	9.25 mm
G-560700	7.00 mm	G-560825	8.25 mm	G-560950	9.50 mm



NANOEDGE VACUUM TREPHINE radial


sterile, single-use

G-570600	6.00 mm	G-570750	7.50 mm	G-570850	8.50 mm
G-570650	6.50 mm	G-570775	7.75 mm	G-570875	8.75 mm
G-570700	7.00 mm	G-570800	8.00 mm	G-570900	9.00 mm
G-570725	7.25 mm	G-570825	8.25 mm		




SINGLE-USE KNIVES FOR OPHTHALMIC SURGERY

6 pcs. per box
sterile



protect, with sliding guard

NANOEDGE PARACENTESIS KNIVES




G-34026 15°
G-34064 15° protect
with sliding guard



G-34031 22,5°



G-34036 30°



NANOEDGE MVR STILETT straight
diameter insertable instruments

G-34181 20 gauge / 0.9 mm
G-34192 23 gauge / 0.6 mm
G-34198 25 gauge / 0.5 mm

incision / blade width

1.2 mm
0.8 mm
0.7 mm



NANOEDGE MVR STILETT angled
diameter insertable instruments

G-34191 20 gauge / 0.9 mm
G-34196 23 gauge / 0.6 mm
G-34197 25 gauge / 0.5 mm

incision / blade width

1.2 mm
0.8 mm
0.7 mm



NANOEDGE PHACO SLIT KNIFE

angled, double bevel


G-34050 1.80 mm
G-34051 2.20 mm
G-34052 2.40 mm
G-34053 2.50 mm
G-34054 2.75 mm
G-34055 2.85 mm
G-34060 2.20 mm, protect, with sliding guard
G-34061 2.40 mm, protect, with sliding guard
G-34062 2.85 mm, protect, with sliding guard



CLEAR CORNEA KNIFE

angled, double bevel

G-34122 2.75 mm
G-34126 3.00 mm



NANOEDGE
CRESCENT KNIFE
angled, bevel up
G-34081 2.0 mm



NANOEDGE
SPOON KNIFE
angled, double bevel
G-34086 2.2 mm



NANOEDGE
IMPLANT KNIFE
angled, bevel up
G-34146 6.0 mm

20G RED: 20 GAUGE / 0.9 MM

23G GREEN: 23 GAUGE / 0.6 MM


25G BLUE: 25 GAUGE / 0.5 MM

SCISSORS

SINGLE-USE VITREOUS SCISSORS

UNO COLORLINE

straight, sterile



G-42061 20G, 1 pcs. per box

G-42361 23G, 1 pcs. per box

G-42561 25G, 1 pcs. per box

SINGLE-USE VITREOUS SCISSORS

UNO COLORLINE

curved, sterile



G-42362 23G, 1 pcs. per box
G-46362 23G, 6 pcs. per box

G-42562 25G, 1 pcs. per box
G-46562 25G, 6 pcs. per box

HANDLES

HANDLE UNO COLORLINE 20G, 23G, 25G

titanium handle for UNO Colorline
single-use tips

G-40006 6 mm diameter

G-40008 8 mm diameter



FORCEPS

SINGLE-USE VITREOUS FORCEPS

ENDGRIPPING UNO COLORLINE

straight, extra delicate,
endgripping 0.3 mm x 0.1 mm, sterile



G-42041 20G, 1 pcs. per box
G-46041 20G, 6 pcs. per box

G-42341 23G, 1 pcs. per box
G-46341 23G, 6 pcs. per box

G-42541 25G, 1 pcs. per box
G-46541 25G, 6 pcs. per box

HATTENBACH

SINGLE-USE VITREOUS FORCEPS

WAVE DESIGN UNO COLORLINE

straight, wave design, 0.3 mm x 0.5 mm,
sterile



G-42042 20G, 1 pcs. per box
G-46042 20G, 6 pcs. per box

G-42342 23G, 1 pcs. per box
G-46342 23G, 6 pcs. per box

G-42542 25G, 1 pcs. per box
G-46542 25G, 6 pcs. per box

SINGLE-USE VITREOUS FORCEPS

CROCODILE DESIGN UNO COLORLINE

straight, crocodile design 0.3 mm x 2.5 mm,
sterile



G-42043 20G, 1 pcs. per box

G-42343 23G, 1 pcs. per box
G-46343 23G, 6 pcs. per box

SINGLE-USE TANO FORCEPS

UNO COLORLINE

straight, sterile



G-42044 20G, 1 pcs. per box

G-42344 23G, 1 pcs. per box

G-42544 25G, 1 pcs. per box

20G RED: 20 GAUGE / 0.9 MM

23G GREEN: 23 GAUGE / 0.6 MM

25G BLUE: 25 GAUGE / 0.5 MM

ENDO PROBES

SINGLE-USE ENDO-PROBE S

UNO COLORLINE

straight, plug for GEUDER endoTRON, sterile



G-61201 20G, 5 pcs. per box

G-61231 23G, 5 pcs. per box

G-61251 25G, 5 pcs. per box

SINGLE-USE ENDO-PROBE S

UNO COLORLINE

straight, plug for IRIDEX® Laser, sterile



G-61204 20G, 5 pcs. per box

G-61234 23G, 5 pcs. per box

G-61254 25G, 5 pcs. per box

SINGLE-USE ENDO-PROBE C

UNO COLORLINE

curved, plug for GEUDER endoTRON, sterile



G-61202 20G, 5 pcs. per box

G-61232 23G, 5 pcs. per box

G-61252 25G, 5 pcs. per box

SINGLE-USE ENDO-PROBE C

UNO COLORLINE

curved, plug for IRIDEX® Laser, sterile



G-61205 20G, 5 pcs. per box

G-61235 23G, 5 pcs. per box

G-61255 25G, 5 pcs. per box

20G RED: 20 GAUGE / 0.9 MM

23G GREEN: 23 GAUGE / 0.6 MM


25G BLUE: 25 GAUGE / 0.5 MM

TROCAR SYSTEMS

HATTENBACH / NIKOLIC

SINGLE-USE TROCAR SYSTEM ONE-STEP UNO COLORLINE

3 ports for one-step technique,
incl. infusion tubing, sterile



G-42011 20G, 1 set per box
G-46011 20G, 6 sets per box

G-42311 23G, 1 set per box
G-46311 23G, 6 sets per box

G-42511 25G, 1 set per box
G-46511 25G, 6 sets per box

SINGLE-USE TROCAR PLUGS UNO COLORLINE



G-42019 20G, 1 set (3 pcs.) per box

G-42319 23G, 1 set (3 pcs.) per box

FIBER OPTICS

SINGLE-USE FIBER OPTIC S UNO COLORLINE

spot, sterile



G-42021 20G, 1 pcs. per box
G-46021 20G, 6 pcs. per box

G-42321 23G, 1 pcs. per box
G-46321 23G, 6 pcs. per box

G-42521 25G, 1 pcs. per box
G-46521 25G, 6 pcs. per box

SINGLE-USE FIBER OPTIC W UNO COLORLINE

wide angle, sterile



G-42022 20G, 1 pcs. per box
G-46022 20G, 6 pcs. per box

G-42322 23G, 1 pcs. per box
G-46322 23G, 6 pcs. per box

G-42522 25G, 1 pcs. per box
G-46522 25G, 6 pcs. per box

20G RED: 20 GAUGE / 0.9 MM

23G GREEN: 23 GAUGE / 0.6 MM

25G BLUE: 25 GAUGE / 0.5 MM

SINGLE-USE FIBER OPTIC WS

UNO COLORLINE

wide angle, shielded, sterile



G-42023 20G, 1 pcs. per box
G-46023 20G, 6 pcs. per box


G-42323 23G, 1 pcs. per box
G-46323 23G, 6 pcs. per box

G-42523 25G, 1 pcs. per box
G-46523 25G, 6 pcs. per box

VITREOUS CUTTERS

SINGLE-USE HIGH SPEED VITREOUS CUTTER UNO COLORLINE

pneumatic, with GEUDER safety-connector,
sterile



G-42001 20G, 1 pcs. per box
G-46001 20G, 6 pcs. per box

G-42301 23G, 1 pcs. per box
G-46301 23G, 6 pcs. per box

G-42501 25G, 1 pcs. per box
G-46501 25G, 6 pcs. per box

SINGLE-USE HIGH SPEED VITREOUS CUTTER MACH 2 UNO COLORLINE

Doubleblade, pneumatic,
with GEUDER safety-connector, sterile



G-42002 20G, 1 pcs. per box
G-46002 20G, 6 pcs. per box

G-42302 23G, 1 pcs. per box
G-46302 23G, 6 pcs. per box


G-42502 25G, 1 pcs. per box
G-46502 25G, 6 pcs. per box

SINGLE-USE INFUSION TUBE 20G

UNO COLORLINE

for 20 gauge GEUDER
UNO Colorline vitreous cutters, sterile
G-42009 20G, 1 pcs. per box





- G-34340** **SINGLE-USE EYE SPECULUM**
 with aspiration, 2 x 4 aspiration ports
 fenestrated, open blades
 5 pcs. per box, sterile



- SZURMAN**
SINGLE-USE DMEK-CARTRIDGE
 for Descemet Membrane Endothelial Keratoplasty
 incl. Tube connection for loading
 the cartridge, sterile
G-38630 incision 3 mm, 14 gauge / 2.0 mm
G-38635 incision 2.4 to 2.75 mm, 16 gauge / 1.6 mm

- G-34523** **SINGLE-USE IRIS RETRACTOR**
 for mechanical dilation of the pupil, flexible
 4 iris retractors per set, 2 sets per box, sterile



- G-34810** **LACRIMAL INTUBATION TUBING**
 silicone tubing 35 cm
 OD 0.6 mm, ID 0.3 mm
 10 pcs. per box, sterile







- G-34811** **LACRIMAL INTUBATION SET**
 suture with silicone tubing
 0.6 mm
 10 pcs. per box, sterile



SINGLE USE INSTRUMENTS


INTUBATION SETS · CANNULAS (LACRIMAL, RETROBULBAR, CYSTOTOME)

G-34812	LACRIMAL INTUBATION SET silicone tubing 0.6 mm 2 probes 110 mm with olive tips 27 gauge / 0.4 mm 3 pcs. per box, sterile	
G-34814	DCR LACRIMAL INTUBATION SET silicone tubing 0.6 mm, 2 probes 45 mm 20 gauge / 0.9 mm 3 pcs. per box, sterile	
G-34290	SINGLE-USE LACRIMAL CANNULA curved 26 gauge / 0.45 x 30 mm 10 pcs. per box, sterile	
G-34200 G-34205	SINGLE-USE RETROBULBAR CANNULA 10 pcs. per box, sterile 25 gauge, 0.5 x 40 mm 26 gauge, 0.45 x 50 mm	
G-34210 G-34215	ATKINSON SINGLE-USE RETROBULBAR CANNULA 10 per box, sterile 25 gauge, 0.5 x 40 mm 23 gauge, 0.64 x 40 mm	
G-34220	SINGLE-USE PERIBULBAR CANNULA 25 gauge, 0.5 x 30 mm 10 pcs. per box, sterile	
G-34225 G-34230	SINGLE-USE CYSTOTOME CANNULA vaulted, 45° 25 gauge / 0.5 mm 10 pcs. per box, sterile 12.5 mm long, side cutting 12.5 mm long	

G-34261 LEWICKY

SINGLE-USE A. C. MAINTAINER

20 gauge / 0.9 mm, tip 5 mm
 3 pcs. per box, sterile



RYCROFT

SINGLE-USE A. C. CANNULA

angled at 5 mm, 45°
 10 pcs. per box, sterile
G-34235 20 gauge, 0.9 x 5 + 17 mm
G-34240 25 gauge, 0.5 x 25 mm



G-34245 SAUTTER

SINGLE-USE HYDRODISSECTION CANNULA

vertically flattened
 27 gauge / 0.4 x 22 mm
 10 pcs. per box, sterile



G-34270 TUEBINGEN MODEL

SINGLE-USE HYDRODISSECTION CANNULA


30 gauge / 0.3 x 22 mm
 10 pcs. per box, sterile








G-34271 FRANKFURT MODEL

SINGLE-USE HYDRODISSECTION CANNULA

27 gauge / 0.4 mm
 10 pcs. per box, sterile





SINGLE USE INSTRUMENTS

G-34265	SINGLE-USE IRRIGATION VECTIS CANNULA double port 25 gauge / 0.5 mm 10 pcs. per box, sterile		
G-34255	KRATZ SINGLE-USE CAPSULE POLISHING CANNULA angled at 9 mm, 45° closed tip with sideport, textured 23 gauge / 0.6 mm 10 pcs. per box, sterile		
G-34280	SIMCOE SINGLE-USE I/A CANNULA 0.4 mm diameter aspiration port front opening irrigation port 23 gauge / 0.64 mm 10 pcs. per box, sterile		
G-15197	SINGLE-USE IRRIGATION CANNULA for viscocanalostomy 26 gauge / 0.45 mm shaft 35 gauge / 0.21 mm angled tip 5 mm tapered tube 5 pcs. per box, sterile		
G-34335	SINGLE-USE IRRIGATION CANNULA for topical irrigation with bulbous tip, straight length 60 mm, 3 mm diameter tip 50 pcs. per box, sterile		
G-80975	SINGLE-USE INJECTION CANNULA for EasyGas® 30 gauge / 0.3 x 12 mm 100 pcs. per box, sterile		

**G-34400 SINGLE-USE BIMANUAL I/A HANDPIECE SET**

consisting of:


- Irrigation Handpiece, irrigation ports 2×0.6 mm
 - Aspiration Handpiece, aspiration port 0.4 mm
 - rough tip tube 21 gauge / 0.8 mm
- 5 pcs. per box, sterile





**G-34289 SINGLE-USE BACKFLUSH HANDPIECE**with silicone chamber, Luer-Lock
10 pcs. per box, sterile**G-34806 SINGLE-USE SILICONE CHAMBER**for backflush handpieces G-31805,
G-37000 and G-37002
10 pcs. per box, sterile**G-34275 SINGLE-USE BACKFLUSH CANNULA**inside diameter 25 gauge / 0.5 mm
outside diameter 21 gauge / 0.8 x 35 mm
10 pcs. per box, sterile**SINGLE-USE BACKFLUSH CANNULA**blunt tip
5 pcs. per box, sterile

- G-34291** 20 gauge / 0.9 mm
G-34296 23 gauge / 0.6 mm
G-34299 25 gauge / 0.5 mm

**SINGLE-USE BACKFLUSH CANNULA**with silicone brush
5 pcs. per box, sterile


- G-34293** 20 gauge / 0.9 mm
G-34294 23 gauge / 0.6 mm
G-34297 25 gauge / 0.5 mm




G-34292	SINGLE-USE BACKFLUSH CANNULA with silicone tip 5 pcs. per box, sterile 20 gauge / 0.9 mm			
G-34295	G-34298	23 gauge / 0.6 mm 25 gauge / 0.5 mm		
G-34285	SINGLE-USE PFCL CANNULA for injection of heavy fluids dual bore, coaxial 23 gauge / 0.64 mm 10 pcs. per box, sterile			
G-34497	G-34498	20 gauge / 0.9 mm x 8 mm 23 gauge / 0.6 mm x 8 mm		
G-33482	SINGLE-USE INJECTION CANNULA for viscous fluids, self-retaining, 4.0 mm tip with 25 cm PVC tube Luer-Lock plastic adapter 20 gauge / 0.9 mm 5 pcs. per box, sterile			
G-33463	SINGLE-USE INFUSION CANNULA self-retaining, 4.0 mm tip with 25 cm silicone tube Luer-Lock plastic adapter 20 gauge / 0.9 mm 5 pcs. per box, sterile			
G-33468	G-33469	SINGLE-USE INFUSION CANNULA with 25 cm silicone tube Luer-Lock plastic adapter 20 gauge / 0.9 mm 5 pcs. per box, sterile beveled tip, 4.0 mm beveled tip, 6.0 mm		
G-33462	SINGLE-USE INFUSION CANNULA 4.0 mm beveled tip with 25 cm PVC tube Luer-Lock plastic adapter 25 gauge / 0.5 mm 10 pcs. per box, sterile			

SINGLE-USE INJECTION CANNULA


for viscous fluids
 with 1 metal sleeve, Luer-Lock plastic adapter
 and 25 cm PVC tube
 20 gauge / 0.9 mm
 5 pcs. per box, sterile
G-33488
G-33489 beveled tip, 4.0 mm
 beveled tip, 6.0 mm

**G-34480****SINGLE-USE INJECTION CANNULA**

for viscous fluids
 4.0 mm beveled tip
 with 2 metal sleeves, Luer-Lock plastic adapter
 and 25 cm silicone tube
 19 gauge / 1.1. mm
 10 pcs. per box, sterile

**G-34492***KIRCHHOF***SINGLE-USE INJECTION CANNULA**

for gas / viscous fluids
 5.0 mm beveled tip with 4 infusion side ports
 2 metal sleeves, Luer-Lock plastic adapter and
 40 cm silicone tube
 20 gauge / 0.9 mm
 10 pcs. per box, sterile

**G-32696****SINGLE-USE PRESSURE TUBE**

for injection of viscous fluid, Luer-Lock female / male
 10 pcs. per box, sterile


**G-31891****SINGLE-USE SYRINGE**

Luer-Lock, 10 ml, sterile


**G-31520****SINGLE-USE MINI FILTER**

0.1 µm, sterile




**G-28766 SINGLE-USE OIL INJECTION SYSTEM**

to inject silicone oil pneumatically,
with protective cover for glass syringe,
pressure tube fits megaTRON® S3 / S4 HPS and
Pentasys® 2, sterile

**G-28767 SINGLE-USE OIL INJECTION SYSTEM**

to inject silicone oil pneumatically,
with protective cover for glass syringe,
pressure tube fits megaTRON® and Accurus®, sterile

**G-28768 SINGLE-USE OIL INJECTION SYSTEM**

to inject silicone oil pneumatically,
with protective cover for glass syringe,
pressure tube fits Millennium®, sterile

- G-23991** **SINGLE-USE THREE-WAY STOP COCK**
 2x Luer-Lock female, 1x Luer-Lock male
 plastic, sterile




- G-32700** **SCALPEL HANDLE**

- SCALPEL SURGICAL BLADES**
 100 pcs. per box, sterile
G-32705 #11, curved cutting edge
G-32710 #15, straight edge and pointed tip




- G-18127** **ABSORBENT SWABS**
 17 mm long
 box of 20 pouches
 10 swabs per pouch
 sterile




- G-33786** **MARKING PEN**
 gentian violet
 dual tip, fine and broad
 to ink corneal instruments
 20 pcs. per box, sterile




G-34602 **SILICONE TIRE (276 STYLE)**
 asymmetrical, 7.0 mm wide
 2.5 mm groove
 5 pcs. per box, sterile




G-34610 **SILICONE TIRE (287 STYLE)**
 convex, 7.0 mm wide
 2.5 mm groove
 5 pcs. per box, sterile



G-34616 **SILICONE TIRE (277 STYLE)**
 concave, 7.0 mm wide
 2.5 mm groove
 5 pcs. per box, sterile



G-34618 **SILICONE TIRE (279 STYLE)**
 concave, 9.0 mm wide
 2.5 mm groove
 5 pcs. per box, sterile





G-34620 SILICONE CIRLING BAND (40 STYLE)

0.75 mm x 2.0 mm

125 mm long



5 pcs. per box, sterile

**G-34622 SILICONE CIRLING BAND (240 STYLE)**

0.6 mm x 2.5 mm

125 mm long


5 pcs. per box, sterile

**G-34624 SILICONE STRIP (41 STYLE)**

0.75 mm x 3.5 mm

125 mm long



5 pcs. per box, sterile

**G-34626 SILICONE STRIP (42 STYLE)**

1.25 mm x 4.0 mm

125 mm long

5 pcs. per box, sterile





G-34632**GROOVED SILICONE STRIP (31 STYLE)**

2.0 mm x 4.5 mm

125 mm long, 2.0 mm groove



5 pcs. per box, sterile

**G-34634****GROOVED SILICONE STRIP (32 STYLE)**

2.5 mm x 5.0 mm

125 mm long, 2.0 mm groove



5 pcs. per box, sterile

**G-34636****GROOVED SILICONE STRIP (219 STYLE)**

1.25 mm x 4.5 mm

100 mm long, 2.5 mm groove

5 pcs. per box, sterile


**G-34644****SILICONE SLEEVE (72 STYLE)**

1.5 mm inside diameter

2.4 mm outside diameter

30 mm long



5 pcs. per box, sterile

**G-34642****SILICONE SLEEVE (70 STYLE)**

1.0 mm inside diameter

2.1 mm outside diameter

30 mm long



**G-34646****SILICONE SLEEVE (270 STYLE)**

0.76 mm inside diameter

1.65 mm outside diameter


30 mm long

5 pcs. per box, sterile




G-34645 OVAL SILICONE SLEEVE (3083 STYLE)


1.5 mm inside diameter
 2.5 mm outside diameter
 5 mm long
 5 pcs. per box, steril

**G-34647 OVAL SILICONE SLEEVE (3084 STYLE)**


1.5 mm inside diameter
 3.75 mm outside diameter
 5 mm long
 5 pcs. per box, steril

**G-34654 SILICONE WEDGE (137 STYLE)**


2.5 mm groove
 5 pcs. per box, sterile





**G-34650 MERIDIONAL SILICONE IMPLANT (106 STYLE)**

6.5 mm wide
 10.5 mm long
 5 pcs. per box, sterile


**G-34656 ROUND SILICONE SPONGE IMPLANT (503 STYLE)**

3.0 mm diameter
 80 mm long
 5 pcs. per box, sterile




G-34658	ROUND SILICONE SPONGE IMPLANT (504 STYLE)	4.0	
	4.0 mm diameter		
	80 mm long		
	5 pcs. per box, sterile		
G-34660	ROUND SILICONE SPONGE IMPLANT (505 STYLE)	5.0	
	5.0 mm diameter		
	80 mm long		
	5 pcs. per box, sterile		
G-34664	oval SILICONE SPONGE IMPLANT (501 STYLE)	4.0	
	2.5 mm x 4.0 mm		
	80 mm long		
	5 pcs. per box, sterile		
G-34666	oval SILICONE SPONGE IMPLANT (506 STYLE)	5.0	
	3.0 mm x 5.0 mm		
	100 mm long		
	5 pcs. per box, sterile		


G-34668 **OVAL SILICONE SPONGE IMPLANT (507 STYLE)**
 5.5 mm x 7.5 mm
 80 mm long
 5 pcs. per box, sterile




G-34672 **GROOVED SILICONE SPONGE IMPLANT (506G STYLE)**
 3.0 mm x 5.0 mm
 80 mm long, 2.5 mm groove
 5 pcs. per box, sterile






G-34674 **GROOVED SILICONE SPONGE IMPLANT (507G STYLE)**
 5.5 mm x 7.5 mm
 80 mm long, 2.5 mm groove
 5 pcs. per box, sterile



G-34684 **SEMI-CIRCULAR SILICONE SPONGE IMPLANT (510 STYLE)**
 2.5 mm x 5.0 mm
 80 mm long
 5 pcs. per box, sterile



G-34686	SEMI-OVAL SILICONE SPONGE IMPLANT (511 STYLE)	
G-34688	SEMI-OVAL SILICONE SPONGE IMPLANT (515 STYLE)	
G-34690	SEMI-OVAL SILICONE SPONGE IMPLANT (517 STYLE)	
G-34710	BONN INJECTION SET	 <p>single-use set for intravitreal injection consisting of:</p> <ul style="list-style-type: none"> - measuring and fixation instrument, 3.5 / 4.0 mm - eye speculum - irrigation cannula - cannula 0.9 x 70 mm - syringe 1 ml - syringe 10 ml - eye drape with pre-cut adhesive incise-film 140 x 140 cm - absorbent swab (3x) - gauze compress 10 x 10 (2x) - protective cover 60 x 60 cm <p>10 pcs. per box, sterile</p>


MISCELLANEOUS




CAPSULE TENSION RING INJECTORS	241
PUPIL DILATOR INJECTOR	241
LACRIMAL / DCR	242
LASEK	243
DIAMOND BURR	244
VITRECTOMY LENSES	244
EYE SHIELD	245
EYE SPHERES	245
CAUTERY	246

**G-32955****INJECTOR**


for capsular tension ring
 shaft diameter 1.2 mm
 for 1.9 mm incisions
 shaft length 30 mm
 overall length 156 mm

**G-32960****INJECTOR**


for *MORCHER* capsular tension ring
 shaft diameter 1.2 mm
 for 1.9 mm incisions
 shaft length 30 mm
 overall length 156 mm

**G-32965****INJECTOR**

for *MORCHER* capsular ring
 shaft diameter 1.6 mm
 for 2.4 - 2.5 mm incisions
 shaft length 30 mm
 overall length 156 mm

**G-32970****INJECTOR**

for *MORCHER* pupil dilator
 overall length 140 mm



G-31810 FRIEDMANN
RONGEUR
overall length 140 mm



G-32355 NASAL SPECULUM
overall length 150 mm



G-33125 MEAD
HAMMER
300 g, lead filled
overall length 190 mm



G-18630 HEIDELBERG MODEL
ELEVATORIUM
overall length 143 mm



G-18620 HEIDELBERG MODEL
RASPATORY
overall length 155 mm





G-18240 *WEST*
BONE CHISEL
 for lacrimal sac
 flat, 5.0 mm wide
 overall length 160 mm



G-18245 *WEST*
BONE GOUGE
 for lacrimal sac
 concave, 5.0 mm wide
 overall length 162 mm



ELSA
LASEK TREPHINE
 with crosshairs
 overall length 100 mm
G-33650
G-33652



LOHMANN
LASEK FUNNEL
 with crosshairs
 overall length 100 mm
G-33655
G-33657



G-32313 *WRIGHT*
FASCIA NEEDLE
 overall length 115 mm



DIAMOND BURR

top, inside, outside coat

G-18300 8.0 mm**G-18301** 9.0 mm**G-18302** 10.0 mm**G-36000***KILP***LENS**60 D flat, field of vision 20°, silicone ring
central posterior pole and central vitreous body**G-36010***LANDERS***LENS**93 D, biconcave with silicone ring
clear view of fundus in air-filled phakic eye**G-36015***MACHEMER***LENS**with silicone ring to magnify structures
in vitreous cavity or on retinal surface**G-36020***PEYMAN***LENS**wide field of vision 30°, with silicone ring
for central fundus and central vitreous body**G-36025***TOLENTINO***LENS**prism 20° with silicone ring
posterior peripheral fundus and
posterior peripheral vitreous body**G-36030***TOLENTINO***LENS**prism 30° with silicone ring
peripheral viewing beyond the equator**G-36045****SILICONE REPLACEMENT RING**for *KILP* lens G-36000

- G-36050 SILICONE REPLACEMENT RING**
 for lenses from G-36010 to G-36035



- G-36065 LENS HANDLE RING**
 for lenses without silicone ring
 to suture on to sclera
 for G-36000



- G-31105 EYESHIELD**
 plastic, transparent for use on both eyes
 100 pcs. per bag




- CIBIS*
EYE SPHERE
- | | |
|----------------|-------|
| G-20216 | 16 mm |
| G-20217 | 17 mm |
| G-20218 | 18 mm |
| G-20219 | 19 mm |
| G-20220 | 20 mm |



- HARMS*
EYE SPHERE
- | | |
|----------------|-------|
| G-20316 | 16 mm |
| G-20317 | 17 mm |
| G-20318 | 18 mm |
| G-20319 | 19 mm |
| G-20320 | 20 mm |



- WALSER*
EYE SPHERE
- | | |
|----------------|-------|
| G-20416 | 16 mm |
| G-20417 | 17 mm |
| G-20418 | 18 mm |
| G-20419 | 19 mm |
| G-20420 | 20 mm |





9/10



STERILIZATION CONTAINERS AND ADAPTERS



MICROSTOP®	249
STEEL CASES	252
TRAYS & RACKS FOR STEEL CASES	256
SILICONE STRIPS & MATS	257
SMALL CASES	258
STERISAFE® DURO	259
ADAPTERS	264

G-39745
**MICROSTOP®
CONTAINER LOW**

for one tray

outside dimensions: 310 x 189 x 70 mm

inside dimensions: 283 x 177 x 50 mm

**G-39750**
**MICROSTOP®
CONTAINER HIGH**

for two trays

outside dimensions: 310 x 189 x 90 mm

inside dimensions: 283 x 177 x 70 mm

**G-39783**
**MICROSTOP®
TRAY MINI**

29 mm, low

size: 277 x 172 x 29 mm

without lid

with lid

**G-39786**
**MICROSTOP®
TRAY MINI**

38 mm, high

size: 277 x 172 x 38 mm

without lid

with lid

**G-39785**

MICROSTOP®

G-39784

STERILIZATION CONTAINERS AND ADAPTERS


G-39787

MICROSTOP®

TRAY LID MINI

size: 277 x 172 mm

for G-39783 and G-39785



G-39782

MICROSTOP®

RETAINER ELEMENT

for MICROSTOP® tray mini

nap Style, 120 mm

separable into 6 single elements

29 mm, low

38 mm, high



G-39789

G-39788

MICROSTOP®

CLIP FOR RETAINER ELEMENTS

G-39782 and G-39789

10 pcs. per bag



G-39791

MICROSTOP®

ID LABEL

for tray, aluminium

customized labeling

silver-colored



G-39741

MICROSTOP®

LABEL

for lid, customized labeling



G-39746

MICROSTOP®

ID LABEL

for front, customized labeling




MICROSTOP®
LOGISTIC FRAME
 for lid to use with labels

- G-39753** color: red
G-39754 color: blue
G-39755 color: green



- G-39762** **MICROSTOP®**
SEAL
 with indicator
 1000 pcs. per box



- G-39764** **MICROSTOP®**
TAG
 with indicator
 1000 pcs. per box



- G-39680** **PROTOCOL CARDS**
 for ERGO-SAFE tray
 320 pcs. per box



G-19822	STERILIZATION CASE	for steam sterilization outside dimensions: 220 x 165 x 55 mm - textile filter in lid and base - silicone seal - perforated insert without edge - perforated silicone tray G-32785 - 2 labels with clamp	
G-19823	STERILIZATION CASE	for steam sterilization outside dimensions: 220 x 165 x 75 mm - textile filter in lid and base - silicone seal - perforated insert without edge - perforated silicone tray G-32785 - 2 labels with clamp	
G-28151	STERILIZATION CASE	for steam sterilization outside dimensions: 300 x 220 x 75 mm - textile filter in lid and base - silicone seal - perforated insert without edge - 3 silicone burl strips - 2 silicone paddings in lid - 2 labels with clamp	
G-28150	STERILIZATION CASE	for steam sterilization outside dimensions: 300 x 220 x 75 mm - textile filter in lid and base - silicone seal in the lid - perforated insert with edge, subdivided, with perforated silicone plate and perforated silicone tray - silicone insert, small, with perforated silicone tray - 2 labels with clamp	

G-19876**STERILIZATION CASE**

for steam sterilization
 outside dimensions: 300 x 220 x 55 mm
 - textile filter in lid and base
 - silicone seal
 - perforated insert without edge
 - 3 silicone burl strips
 - 2 silicone paddings in lid
 - 2 labels with clamp

**G-19856****STERILIZATION CASE**

for steam sterilization
 outside dimensions: 300 x 220 x 55 mm
 - paper filter with quick-clasp in lid and base
 - silicone seal in the lid
 - without insert and profile strips
 - 2 labels with clamp

**G-19843****STERILIZATION CASE**

for steam sterilization
 outside dimensions: 300 x 220 x 75 mm
 - textile filter in lid and base
 - silicone seal
 - perforated insert with edge
 - perforated silicone plate
 - perforated insert without edge
 - 3 silicone burl strips
 - 2 silicone paddings in lid
 - 2 labels with clamp



G-19842

STERILIZATION CASE

- for steam sterilization
- outside dimensions: 300 x 220 x 75 mm
- textile filter in lid and base
- silicone seal in the lid
- perforated insert with edge and perforated silicone plate
- perforated insert with edge perforated silicone plate and silicone instrument tray
- 2 labels with clamp



G-19877

STERILIZATION CASE

- for steam sterilization
- outside dimensions: 400 x 290 x 55 mm
- textile filter in lid and base
- silicone seal
- perforated insert without edge
- 4 silicone burl strips
- 3 silicone paddings in lid
- 2 labels with clamp



G-19892


STERILIZATION CASE

- for steam sterilization
- outside dimensions: 400 x 290 x 70 mm
- textile filter in lid and base
- silicone seal in the lid
- 2 perforated inserts without edge with 4 silicone nap strips each
- 3 silicone profile strips at the base of the upper perforated insert
- 3 silicone profile strips in the lid
- 2 labels with clamp



G-19894**STERILIZATION CASE**

for steam sterilization
 outside dimensions: 400 x 290 x 70 mm
 - textile filter in lid and base
 - silicone seal in the lid
 - perforated insert without edge
 with 4 silicone nap strips
 - perforated insert with edge and
 silicone nap plate
 - 3 silicone profile strips in the lid
 - 2 labels with clamp

**G-19890****STERILIZATION CASE**

for steam sterilization
 outside dimensions: 495 x 270 x 70 mm
 - textile filter in lid and base
 - silicone seal in the lid
 - perforated insert with edge
 4 silicone nap strips and lateral
 compartment (240 x 110 x 40 mm),
 - 3 silicone profile strips in the lid
 - 2 labels with clamp

**PERMANENT TEXTILE FILTER**

for steam sterilization cases

- G-31650** for G-19822, G-19823
- G-31651** for G-19841, G-19842, G-19843, G-19876, G-28150, G-28151
- G-31652** 6 perforations, for G-19877, G-19892, G-19894
- G-31656** for G-19877
- G-31661** for G-19890

**G-31649****INDICATOR**

for steam sterilization, slim
 1000 pcs. per box

Steri.-Datum / Date:	<input type="text"/>
Verfall / Valid until:	<input type="text"/>
Name:	<input type="text"/>

DAMPF = SCHWARZ
STEAM = BLACK
FREI = UND

STERILIZATION CONTAINERS AND ADAPTERS

G-32675 SILICONE TRAY

with bottom

outside dimensions: 160 x 85 x 20 mm



G-32660 STERILIZATION RACK


for 12 cannulas, silicone

outside dimensions: 85 x 56 x 8 mm




G-19870 STERILIZATION RACK

for 10 instruments with handle, silicone
outside dimensions: 160 x 85 x 20 mm



G-32785 SILICONE PERFORATED TRAY

outside dimensions: 194 x 135 x 24 mm



G-39650**STERILIZATION CASE**

for steam sterilization, plastic

for 6 vitreoretinal instruments

outside dimensions: 210 x 130 x 26 mm

**G-32676****SILICONE NAP STRIP BLUE**

340 x 30 mm

**G-32677****SILICONE STRIP BLUE**

350 mm long

**SILICONE BURL MAT**

perforated

G-19866


240 x 230 mm

G-19865

300 x 150 mm

G-19867

470 x 245 mm


**G-19899****SILICONE PLATE**

perforated, 270 x 185 mm


for hot air sterilization cases

G-19841, G-19842, G-19843,

G-19846, G-19847



STERILIZATION CONTAINERS AND ADAPTERS

G-31657	SAFETY SEAL for steam sterilization cases 100 pcs. per box	
G-31659	SEAL for hot-air sterilization cases 100 pcs. per box	
G-32695	SILICONE CLOSURE 17 cm long	
G-33080	STERILIZATION CASE for steam sterilization, plastic, for 5 instruments with handle outside dimensions: 160 x 85 x 20 mm	
G-33090	STERILIZATION CASE for one instrument or diamond knife, plastic outside dimensions: 160 x 60 x 20 mm	
G-33100	STERILIZATION CASE for one vitreoretinal instrument, plastic outside dimensions: 195 x 60 x 25 mm	

G-36090**STERILIZATION CASE**

plastic, for 6 vitreoretinal lenses

**G-39527****STERISAFE® DURO A4 PLUS**

container for steam, gas (EO, FO) and plasma sterilization, without trays and strips
 outside dimensions: 290 x 190 x 70 mm
 inside dimensions: 265 x 165 x 60 mm

**G-39526****STERISAFE® DURO A4 PLUS COMPLETE SET**

container for steam, gas (EO, FO) and plasma sterilization, including:
 1 x Sterisafe® wire basket A4-30
 5 x Sterisafe® silicone retainer strips 230 mm
 1 x Sterisafe® metal tray A4-30
 outside dimensions: 290 x 190 x 70 mm
 inside dimensions: 265 x 165 x 60 mm

**G-39525****STERISAFE® DURO A4**

container for steam, gas (EO, FO) and plasma sterilization
 outside dimensions: 290 x 190 x 60 mm
 inside dimensions: 265 x 165 x 49 mm

**G-39528****STERISAFE® DURO A6**

container for steam, gas (EO, FO) and plasma sterilization
 outside dimensions: 190 x 110 x 60 mm
 inside dimensions: 165 x 92 x 49 mm

STERILIZATION CONTAINERS AND ADAPTERS

G-39543 **STERISAFE® DURO WIRE BASKET A4-30**
outside dimensions: 260 x 163 x 30 mm
inside dimensions: 249 x 157 x 23 mm



G-39541 **STERISAFE® DURO METAL TRAY A4-15**
outside dimensions: 260 x 163 x 15 mm
inside dimensions: 249 x 157 x 11 mm



G-39542 **STERISAFE® DURO METAL TRAY A4-30**
outside dimensions: 260 x 163 x 30 mm
inside dimensions: 249 x 157 x 26 mm



G-39549 **STERISAFE® DURO WIRE BASKET A6-30**
outside dimensions: 162 x 90 x 30 mm
inside dimensions: 156 x 80 x 23 mm



G-39548 **STERISAFE® DURO METAL TRAY A6-15**
outside dimensions: 162 x 90 x 15 mm
inside dimensions: 156 x 80 x 11 mm



- G-39545** **STERISAFE® SILICONE RETAINER STRIP 75**
 for Sterisafe® DURO A4 and A6
 length: 75 mm



- G-39546** **STERISAFE® SILICONE RETAINER STRIP 145**
 for Sterisafe® DURO A4 and A6
 length: 145 mm



- G-39547** **STERISAFE® SILICONE RETAINER STRIP 230**
 for Sterisafe® DURO A4 and A6
 length: 230 mm

**STERISAFE® FILTER**

- G-39501** A4 10 pcs. per box
G-39511 A6 10 pcs. per box



- G-39530** **STERISAFE® FILTER TAG EXPIRATION DATE**
 10 pcs. per box



- G-39535** **SEAL / LABEL**
 seal with color indicator and blank label
 2000 pcs. per box


Steri.- Datum/Date	
Name	
	Dampf (braun) Steam (brown)

STERILIZATION CONTAINERS AND ADAPTERS


G-39540 SECURITY SEAL
100 pcs. per box



G-39502 STERISAFE® SILICONE BURL MAT
A4
G-39512 A6



G-39503 STERISAFE® SILICONE GRID INSERT
A4
G-39513 A6




G-39504 STERISAFE® METAL GRID INSERT
A4
G-39514 A6



G-39505

STERISAFE® METAL GRID INSERT

with 2 silicone burl strips




G-31670

MELAG METAL GRID INSERT

with 2 silicone nap strips

outside dimensions: 370 x 147 x 50 mm

for MELAG sterilizing case



STERILIZATION CONTAINERS AND ADAPTERS

G-31028

CLEANING ADAPTER

for ultrasonic handpieces (irrigation connector)
Luer female / male



G-31029

CLEANING ADAPTER

for ultrasonic handpieces (aspiration connector)
Luer female / female



G-31021

CLEANING ADAPTER

for ultrasonic tips



G-31013

CLEANING ADAPTER

for HATTENBACH 25 gauge hybrid instruments
(G-32944, G-32947 and G-32948)

G-22636

ADAPTER

for PAPST I/A handpiece
light model



G-31025

ASPIRATION ADAPTER

metal
for cleaning I/A handpieces
Luer-Lock female / female



G-31026

ASPIRATION ADAPTER

plastic
for cleaning I/A handpieces
Luer-Lock female / female



- G-31016** **CLEANING ADAPTER**
 for diamond knives, small



- G-31017** **CLEANING ADAPTER**
 for diamond knives, large



- G-31022** **CLEANING ADAPTER**
 for G-37010 - G-37080



- G-31031** **CLEANING ADAPTER**
 for aspiration handpiece G-37100 and G-37105
 Luer-Lock female / silicone tube
 2.1 x 1.0 mm



- G-31018** **CLEANING ADAPTER**
 for PFCL cannula G-33055
 Luer-Lock female / silicone tube



- G-31024** **CLEANING ADAPTER**
 for eye specula with aspiration



- G-31023** **CLEANING ADAPTER**
 for eye specula with aspiration
 G-33141, G-33142 and G-33143



- G-31027** **CLEANING ADAPTER**
 for JUENEMANN probes





G-31034

CLEANING ADAPTER

for foreign body instruments G-14530, G-14540, G-14550
Luer-Lock female

G-31014

CLEANING ADAPTER

for cannulas G-33057 and G-33059
Luer-Lock female



G-36604

CLEANING ADAPTER


for rotatable and interchangeable tips



PLEASE READ YOUR INSTRUMENT'S
CLEANING INSTRUCTIONS OR ASK
YOUR LOCAL SALES REPRESENTATIVE
FOR THE RIGHT CLEANING ADAPTER.


19

SETS



Geuder®
Precision made in Germany

ECCE - EXTRA CAPSULAR CATARACT EXTRACTION




page

1	G-17020	BARRAQUER EYE SPECULUM	fenestrated.....	15
2	G-15960	LIEBERMANN-TENNANT EYE SPECULUM	25
3	G-18180	HALSTEAD MOSQUITO FORCEPS	straight.....	135
4	G-18820	REMKY SUPERIOR RECTUS FORCEPS	125
5	G-17575	CASTROVIEJO NEEDLE HOLDER	curved, w. lock.....	142
6	G-19400	BONN EYE SCISSORS	curved, delicate.....	92
7	G-19750	WESTCOTT SCISSORS	blunt tips, curved, 12 cm.....	89
8	G-18791	DRESSING FORCEPS	1 x 2 teeth, 0.7 mm.....	110
9	G-18190	DIEFFENBACH SERREFINE	straight, 35 mm.....	136
10	G-15185	FRANKFURT HYDRODISSECTION CANNULA	27G.....	192
11	G-15250	SAUTTER HYDRODIS. CANNULA	27G, vert flat.....	192
12	G-16125	KRATZ CAPSULE POLISHING CANNULA	195
13	G-16115	IRRIGATION VECTIS CANNULA	193
14	G-16106	SIMCOE I/A CANNULA	23G, dual shaft.....	194
15	G-19710	TROUTMAN CORNEAL SCISSORS	right, curved.....	81
16	G-19715	TROUTMAN CORNEAL SCISSORS	left, curved.....	81
17	G-19776	CAPSULOTOMY SCISSORS	angled up.....	83
18	G-19765	VANNAS CAPSULOTOMY SCISSORS	curved.....	84
19	G-31305	CAPSULORHEXIS FORCEPS	extra delicate.....	114

page

20	G-18960	MICRO COLIBRI FORCEPS	with tying platform.....	104
21	G-19030	TUEBINGEN TYING FORCEPS	straight.....	106
22	G-19032	TYING FORCEPS	angled, 7.5 mm.....	107
23	G-16189	SPATULA	angled, 0.5 mm wide.....	159
24	G-15620	ARLT LENS LOOP	serrated.....	161
25	G-18910	BONN IRIS FORCEPS	angled, 8.0 mm, 1 x 2 teeth.....	110
26	G-19130	CASTROVIEJO CALIPER	15 mm.....	43
27	G-16175	DARDENNE IRIS HOOK	push-pull.....	157
28	G-18935	DODEN IMPLANTATION FORCEPS	121
29	G-16167	SINSKEY LENS MANIPULATION HOOK	0.2 mm.....	156
30	G-34026	NANOEDGE PARACENTESIS KNIFE	15° straight.....	217
31	G-34146	NANOEDGE IMPLANT KNIFE	angled, 6 mm.....	217
32	G-34523	SINGLE-USE IRIS RETRATOR	flexible, 2 x 4 pcs.....	223
33	G-34225	SINGLE-USE CYSTOTOME CANNULA	224
34	G-17500	BARRAQUER NEEDLE HOLDER	curved, w/o lock.....	139
35	G-39786	MICROSTOP® TRAY MINI	low, with lid.....	249
36	G-39783	MICROSTOP® TRAY MINI	low, without lid.....	249
37	G-39782	MICROSTOP® RETAINER ELEMENT	low.....	250
38	G-39750	MICROSTOP® CONTAINER HIGH	for two trays.....	249


PHACO / CATARACT



	page
1 G-15960 LIEBERMANN-TENNANT EYE SPECULUM	25
2 G-32148 BARRAQUER EYE SPECULUM fenestrated	16
3 G-34210 ATKINSON SINGLE-USE RETROBULBAR CANNULA	224
4 G-15250 SAUTTER HYDRODIS. CANNULA 27G, vert flat	192
5 G-15151 LACRIMAL CANNULA conical, short angle, steel	189
6 G-15200 CRANZ EYE DOUCHE silver	191
7 G-19023 CASTROVIEJO SUTURING FORCEPS new, 0.3 mm	103
8 G-33758 CLEAR-CORNEA-INCISION FIXATION RING	150
9 G-19485 STEVENS TENOTOMY SCISSORS curved, long	89
10 G-34026 NANOEDGE PARACENTESIS KNIFE 15° straight	217
11 G-34055 NANOEDGE PHACO SLIT KNIFE angled, 2.85 mm	217
12 G-34051 NANOEDGE PHACO SLIT KNIFE angled, 2.2 mm	217
13 G-34052 NANOEDGE PHACO SLIT KNIFE angled, 2.4 mm	217
14 G-18960 MICRO COLIBRI FORCEPS with tying platform	104
15 G-32162 THOMAS NEUHANN IRIS HOOK (Push-Pull) and NUCLEUS DIVIDING HOOK	154
16 G-19765 VANNAS CAPSULOTOMY SCISSORS curved	84
17 G-16189 SPATULA angled, 0.5 mm wide	159

	page
18 G-32926 AUFFARTH MICS DOUBLE CROSS ACTION CAPSULORHEXIS FORCEPS flat handle	115
19 G-31305 CAPSULORHEXIS FORCEPS extra delicate	114
20 G-31720 AKAHOSHI PRECHOPPER for micro-incisions	120
21 G-31507 NAGAHARA NUCLEUS DIVIDING HOOK	152
22 G-16167 SINSKEY LENS MANIPULATION HOOK 0.2 mm	156
23 G-19030 TUEBINGEN TYING FORCEPS straight	106
24 G-19032 TYING FORCEPS angled, 7.5 mm	107
25 G-17500 BARRAQUER NEEDLE HOLDER curved, w/o lock	139
26 G-19130 CASTROVIEJO CALIPER 15 mm	43
27 G-31901 HOLDING FORCEPS for foldable IOLs	124
28 G-19400 BONN EYE SCISSORS curved, delicate	92
29 G-22100 KOCH-KOHNNEN ASPIRATION HANDPIECE 22G	198
30 G-22101 KOCH-KOHNNEN IRRIGATION HANDPIECE 22G	199
31 G-39786 MICROSTOP® TRAY MINI low, with lid	249
32 G-39783 MICROSTOP® TRAY MINI low, without lid	249
33 G-39782 MICROSTOP® RETAINER ELEMENT low	250
34 G-39750 MICROSTOP® CONTAINER HIGH for two trays	249

PHACO EMERGENCY



	page
1 G-32996 SNYDER-OSHER IOL SCISSORS 18G/1.2mm, long	88
2 G-32997 SNYDER-OSHER IOL FORCEPS 22G / 0.7 mm	125
3 G-32945 STEINMETZ IOL SCISSORS	87
4 G-31994 H.-R. KOCH PREFOLDER for foldable IOLs	124
5 G-31905 H.-R. KOCH FOLDER for 3-piece IOLs	124
6 G-18935 DODEN IMPLANTATION FORCEPS	121
7 G-16167 SINSKEY LENS MANIPULATION HOOK 0.2 mm	156
8 G-32935 H.-R. KOCH CAPSULE FORCEPS shaft 22G	115
9 G-32940 H.-R. KOCH CAPSULE SCISSORS 45° angled	85
10 G-32944 HATTENBACH CAPSULORHEXIS FORCEPS angled	117


	page
11 G-32947 HATTENBACH HYBRID TYING FORCEPS 25G	108
12 G-32948 HATTENBACH HYBRID SCISSORS 25G	86
13 G-32949 HATTENBACH HYBRID NEEDLE HOLDER angled	143
14 G-34523 SINGLE-USE IRIS RETRACTOR flexible, 2 x 4 pcs.	223
15 G-32823 EXPLANTATION FORCEPS for soft IOLs	125
16 G-32800 H.-R. KOCH SILICONE LENS BI-SECTOR	88
17 G-32960 INJECTOR for MORCHER capsule tension ring	241
18 G-32970 INJECTOR for MORCHER pupil dilator	241
19 G-39526 STERISAFE DURO A4 PLUS COMPLETE SET	259

MSICS - MANUAL SMALL INCISION CATARACT SURGERY



	page
1 G-32136 BARRAQUER EYE SPECULUM fenest., round	16
2 G-15960 LIEBERMANN-TENNANT EYE SPECULUM	25
3 G-15620 ARLT LENS LOOP serrated, 3.5 mm wide	161
4 G-16167 SINSKEY LENS MANIPULATION HOOK 0.2 mm	156
5 G-17500 BARRAQUER NEEDLE HOLDER curved, w/o lock	139
6 G-18950 COLIBRI FORCEPS with tying platform	104
7 G-18935 DODEN IMPLANTATION FORCEPS	121

	page
8 G-31622 GASKIN FRAGMENT FORCEPS angled, 9.0 mm	108
9 G-15250 SAUTTER HYDRODIS. CANNULA 27G, vert flat	192
10 G-19765 VANNAS CAPSULOTOMY SCISSORS curved	84
11 G-16106 SIMCOE I/A CANNULA 23G, dual shaft	194
12 G-31307 CAPSULORHEXIS FORCEPS strong	113
13 G-19400 BONN EYE SCISSORS curved, delicate	92


MICS - MICRO INCISION CATARACT SURGERY


1 2 3 4 5 6 7 8 9

page

1	G-38448 DIACLEAN CCI DIAMOND KNIFE trapez.	1.5/2.0.....60
2	G-32926 AUFFARTH MICS CAPSULORHEXIS FORCEPS115
3	G-38240 JAPAN HANDLE for capsule forceps and scissors116
4	G-38244 IKEDA CAPSULORHEXIS FORCEPS 23G, micro-tips	116
5	G-38246 YAE CAPSULE SCISSORS 23G, interchangeable85

	page
6	G-32012 IRRIGATION HANDPIECE 2x0.5mm,con.,text.,21G.....200
7	G-32011 ASPIRATION HANDPIECE 0.3 mm port, 21G.....200
8	G-31507 NAGAHARA NUCLEUS DIVIDING HOOK152
9	G-32505 TH. NEUHANN NUCLEUS DIVIDING HOOK151

FLACS - FEMTOSECOND LASER-ASSISTED CATARACT SURGERY


1 2 3 4 5 6 7 8 9

page

1	G-33965 DICK FEMTOCATARACT EYE SPECULUM26
2	G-33966 DICK FEMTOCATARACT SPATULA double ended	169
3	G-33967 BREYER FEMTOCATARACT SPATULA blunt tip	169
4	G-33968 HOFFMANN FEMTOCATARACT NUCLEUS ROTATOR AND CHOPPER	169
5	G-32926 AUFFARTH MICS CAPSULORHEXIS FORCEPS115

	page
6	G-15250 SAUTTER HYDRODIS. CANNULA 27G, vert flat.....192
7	G-34245 SAUTTER SINGLE-USE HYDRODISSECTION CANNULA225
8	G-22100 KOCH-KOHNNEN ASPIRATION HANDPIECE 22G.....198
9	G-22101 KOCH-KOHNNEN IRRIGATION HANDPIECE 22G.....199
10	G-39526 STERISAFE DURO A4 PLUS COMPLETE SET259

SNYDER-OSUMER IOL EXPLANTATION



1 2

- page
1 G-32997 SNYDER-OSUMER IOL FORCEPS 22G / 0.7 mm.....125
2 G-32998 SNYDER-OSUMER IOL SCISSORS 18G/1.2mm, short....88

TORIC IOL



- page
1 G-33777 TOMARK ADJUSTABLE CORNEAL MARKER.....47
2 G-33776 TOMARK PENDULUM INSTRUMENT handheld.....47
3 G-33774 SCOTT CORNEAL PENDULUM MARKER.....46
4 G-33763 GERTEN HORIZONTAL AXIS MARKER.....46
5 G-16130 CAPSULE POLISHING CANNULA 18G / 1.2 mm....195
6 G-32564 RENTSCH CAPSULE POLISHING CURETTE 90°.....160
7 G-33786 MARKING PEN gentian violet, 20 pcs./box.....231


HATTENBACH IRIS SUTURE



1 2 3

- page
1 G-32949 HATTENBACH HYBRID NEEDLE HOLDER angled....143
2 G-32947 HATTENBACH HYBRID TYING FORCEPS 25G.....108
3 G-32948 HATTENBACH HYBRID SCISSORS 25G.....86

LASEK - LASER-ASSISTED EPITHELIAL KERATOMILEUSIS




	page		page
1 G-15961 LIEBERMANN-TENNANT ASPIRATION SPECULUM	28	5 G-33657 LOHMANN LASEK FUNNEL 9.5 mm, w crosshairs	243
2 G-33650 LASEK TREPHINE 8.0 mm, with crosshairs	243	6 G-33670 LOHMANN LASEK EPI PEELER	167
3 G-33652 LASEK TREPHINE 9.0 mm, with crosshairs	243	7 G-33930 TH. NEUHANN LASIK SPATULA	167
4 G-33655 LOHMANN LASEK FUNNEL 8.5 mm, w crosshairs	243	8 G-15180 IRRIGATION CANNULA bulbous tip, 26G	194

LASIK - LASER-ASSISTED IN SITU KERATOMILEUSIS




	page		page
1 G-33920 KNORZ LASIK EYE SPECULUM w aspiration	26	10 G-18914 BURATTO LASIK FLAP FORCEPS	119
2 G-15960 LIEBERMANN-TENNANT EYE SPECULUM	25	11 G-33943 PFAEFFL FLAP LIFTER V-shaped lifter and hook	171
3 G-33745 MENDEZ LASIK CORNEAL MARKER	46	12 G-33912 ZIRM LASIK CORNEAL PROTECTOR	167
4 G-33746 CHAYET LASIK CORNEAL MARKER	46	13 G-33940 BLUM FEMTO-LASIK ROUND SPOON bevel up	170
5 G-33747 LU LASIK CORNEAL MARKER	45	14 G-33759 FINE-THORNTON FIXATION RING 13 mm, blunt	150
6 G-33910 ZIRM LASIK DISSECTOR CANNULA 25G, flat	197	15 G-33786 MARKING PEN gentian violet, 20 pcs./box	231
7 G-33911 ZIRM LASIK IRRIGATION CANNULA 20G, flat	197	16 G-18150 HOLDER FOR ABSORBENT STRIPS	145
8 G-33909 BAUER LASIK CANNULA double-armed, 25G	196	17 G-18127 ABSORBENT SWABS 17 mm, sterile	231
9 G-33428 PFAEFFL FLAPORHEXIS FORCEPS oval tips	118	18 G-39526 STERISAFE DURO A4 PLUS COMPLETE SET	259

FEMTO LASIK



1 G-15960 LIEBERMANN-TENNANT EYE SPECULUM.....	25	7 G-33427 PFAEFFL FLAP LIFTER curved lifter and hook.....	171
2 G-15961 LIEBERMANN-TENNANT ASPIRATION SPECULUM.....	28	8 G-33943 PFAEFFL FLAP LIFTER V-shaped lifter and hook.....	171
3 G-18914 BURATTO LASIK FLAP FORCEPS.....	119	9 G-33758 CLEAR-CORNEA-INCISION FIXATION RING.....	150
4 G-33428 PFAEFFL FLAPORHEXIS FORCEPS oval tips.....	118	10 G-33930 TH. NEUHANN LASIK SPATULA.....	167
5 G-33429 KNORZ FEMTO-LASIK FLAP FORCEPS T-bar tips.....	118	11 G-18150 HOLDER FOR ABSORBENT STRIPS.....	145
6 G-33912 ZIRM LASIK CORNEAL PROTECTOR.....	167	12 G-18127 ABSORBENT SWABS 17 mm, sterile.....	231


RELEX®* SMILE - SMALL INCISION LENTICULE EXTRACTION



	page		page
1 G-15960 LIEBERMANN-TENNANT EYE SPECULUM.....	25	5 G-33960 RUPAL SHAH LENTICULE MICRO FORCEPS.....	119
2 G-15961 LIEBERMANN-TENNANT ASPIRATION SPECULUM.....	28	6 G-33961 LENTICULE FORCEPS angled to side, 7.5 mm.....	119
3 G-33952 CHANSUE FEMTO DOUBLE INSTRUMENT spatula/hook.....	170	7 G-33759 FINE-THORNTON FIXATION RING 13 mm, blunt.....	150
4 G-33950 CHANSUE FEMTO DOUBLE INSTRUMENT spoon-spatula/hook.....	169	8 G-33941 KELMAN-MCPHERSON FRAGMENT FORCEPS.....	108
		9 G-33959 BREYER FEMTO-LASIK HOOK.....	171
		10 G-33910 ZIRM LASIK DISSECTOR CANNULA 25G, flat.....	197

*ReLex® (Refractive Lenticule Extraction) is a registered trademark of Carl ZeissMeditec AG.

KERATOPLASTY




	page		page
1 G-15966 KERSHNER EYE SPECULUM solid blades.....	25	18 G-15758 KERATOPLASTY SPATULA perforated.....	163
2 G-18180 HALSTEAD MOSQUITO FORCEPS straight.....	135	19 G-14200 PAUFIQUE KERATOPLASTY KNIFE.....	64
3 G-18791 DRESSING FORCEPS 1 x 2 teeth, 0.7 mm	110	20 G-18972 THOMAS NEUHANN KERATOPLASTY FORCEPS	104
4 G-17450 JACOBI NEEDLE HOLDER curved, w/o lock.....	139	21 G-32229 DUNCKER SUTURE MARKER for 8.0 mm transplant.....	42
5 G-17500 BARRAQUER NEEDLE HOLDER curved, w/o lock.....	139	22 G-33727 THOMAS NEUHANN CORNEAL MARKER 8 blades.....	41
6 G-15065 FLIERINGA FIXATION RING 15 mm.....	150	23 G-18960 MICRO COLIBRI FORCEPS with tying platform.....	104
7 G-15070 FLIERINGA FIXATION RING 16 mm.....	150	24 G-19001 CASTROVIEJO SUTURING FORCEPS 0.12 mm.....	103
8 G-15075 FLIERINGA FIXATION RING 17 mm.....	150	25 G-19010 CASTROVIEJO SUTURING FORCEPS 0.6 mm.....	103
9 G-15080 FLIERINGA FIXATION RING 18 mm.....	150	26 G-19030 TUEBINGEN TYING FORCEPS straight.....	106
10 G-15085 FLIERINGA FIXATION RING 19 mm.....	150	27 G-19040 TUEBINGEN TYING FORCEPS curved.....	106
11 G-15090 FLIERINGA FIXATION RING 20 mm.....	150	28 G-15250 SAUTTER HYDRODIS. CANNULA 27G, vert flat.....	192
12 G-14020 KUHNT CORNEAL SCARIFIER fig. 2.....	64	29 G-32318 MALONEY KERATOMETER.....	45
13 G-19400 BONN EYE SCISSORS curved, delicate.....	92	30 G-570750 NANOEDGE VACUUM TREPHINE 7.5 mm.....	216
14 G-19750 WESTCOTT SCISSORS blunt tips, curved, 12 cm.....	89	31 G-570800 NANOEDGE VACUUM TREPHINE 8.0 mm.....	216
15 G-19720 KATZIN KERATOPLASTY SCISSORS right, 10 cm.....	81	32 G-560800 NANOEDGE VACUUM CORNEAL PUNCH 8.0 mm.....	216
16 G-19725 KATZIN KERATOPLASTY SCISSORS left, 10 cm.....	81	33 G-560850 NANOEDGE VACUUM CORNEAL PUNCH 8.5 mm.....	216
17 G-19765 VANNAS CAPSULOTOMY SCISSORS curved.....	84		

DSEK / DSAEK - DESCemet's STRIPPING (AUTOMATED) ENDOTHELIAL KERATOPLASTY




	page
1 G-15960 LIEBERMANN-TENNANT EYE SPECULUM	25
2 G-33416 BOPP A.C. INFUSION CANNULA silicone tube	191
3 G-38600 NEUHANN CORNEAL MARKER double 8 / 9 mm	41
4 G-38607 DESCemet INCISION HOOK 0.2 mm tip	168
5 G-38608 DESCemet SCRAPER	168
6 G-38601 ALTHAUS/CARTSBURG IRRIGATION DESCemet HOOK	202
	page
7 G-38602 NEUHANN IRRIGATION DESCemet SCRAPER	202
8 G-38634 GERTEN DSAEK/DMEK RHEXIS FORCEPS 23G	117
9 G-38603 DSAEK IMPLANTATION FORCEPS	118
10 G-38606 TH. NEUHANN DSAEK FORCEPS 21G	118
11 G-38604 TH. NEUHANN DSAEK INSERTION SPATULA	168

DALK - DEEP ANTERIOR LAMELLAR KERATOPLASTY




	page
1 G-15960 LIEBERMANN-TENNANT EYE SPECULUM	25
2 G-32728 RENTSCH SPATULA tapered 0.6 to 0.3 mm	159
3 G-19720 KATZIN KERATOPLASTY SCISSORS right, 10 cm	81
4 G-19725 KATZIN KERATOPLASTY SCISSORS left, 10 cm	81
5 G-19770 VANNAS CAPSULOTOMY SCISSORS angled to side	84
6 G-15310 KRUMEICH CORNEAL IRRIGATION CANNULA 27G	197
	page
7 G-15312 KRUMEICH CORNEAL IRRIGATION CANNULA 23G	197
8 G-19030 TUEBINGEN TYING FORCEPS straight	106
9 G-19040 TUEBINGEN TYING FORCEPS curved	106
10 G-17505 BARRAQUER NEEDLE HOLDER curved, w/o lock	139
11 G-570800 NANOEDGE VACUUM TREPHINE 8.0 mm	216

DMEK - DESCemet MEMBRANE ENDOTHELIAL KERATOPLASTY


	page		page
1 G-38600 NEUHANN CORNEAL MARKER double 8 / 9 mm.....	41	8 G-38607 DESCemet INCISION HOOK 0.2 mm tip	168
2 G-38620 SZURMAN DMEK GRAFT RETAINER.....	71	9 G-38608 DESCemet SCRAPER.....	168
3 G-38623 REMKY DMEK KERATOPLASTY FORCEPS.....	105	10 G-38601 ALTHAUS-CARTSBURG IRRIGATION DESCemet HOOK.....	202
4 G-38621 SZURMAN DMEK PREPARATION SCRAPER.....	64	11 G-38602 NEUHANN IRRIGATION DESCemet SCRAPER.....	202
5 G-38622 SZURMAN DMEK PREPARATION SPATULA.....	168	12 G-38634 GERTEL DSAEK/DMEK RHEXIS FORCEPS 23G.....	117
6 G-38624 SZURMAN DMEK PREPARATION FORCEPS.....	117	13 G-38635 SZURMAN SINGLE-USE DMEK CARTRIDGE 1.6 mm.....	223
7 G-32167 SAUTTER HYDRODIS. CANNULA 27G, horiz flat.....	192		


GLAUCOMA



	page
1 G-15960 LIEBERMANN-TENNANT EYE SPECULUM	25
2 G-17020 BARRAQUER EYE SPECULUM fenestrated	15
3 G-33200 SCHAEDEL TOWEL FORCEPS	136
4 G-34036 NANOEDGE PARACENTESIS KNIFE 30° straight	217
5 G-18180 HALSTEAD MOSQUITO FORCEPS straight	135
6 G-18820 REMKY SUPERIOR RECTUS FORCEPS	125
7 G-17565 CASTROVIEJO NEEDLE HOLDER straight, w. lock	142
8 G-17500 BARRAQUER NEEDLE HOLDER curved, w/o lock	139
9 G-19400 BONN EYE SCISSORS curved, delicate	92
10 G-19750 WESTCOTT SCISSORS blunt tips, curved, 12 cm	89
11 G-19630 TUEBINGEN CAPSULOTOMY SCISSORS curved	84
12 G-19680 BARRAQUER SCISSORS 7 mm blades	87
13 G-19130 CASTROVIEJO CALIPER 15 mm	43
14 G-15185 FRANKFURT HYDRODISSECTION CANNULA 27G	192
15 G-15250 SAUTTER HYDRODIS. CANNULA 27G, vert flat	192
16 G-22619 PH. JACOBI TRABECULAR ASPIRATION CANNULA	193
17 G-20612 SCHMIDT-MUMM DUCT TREPHINE 1.2 mm	69
18 G-15125 HARM'S TRABECULAR PROBE right	186
19 G-15130 HARM'S TRABECULAR PROBE left	186

	page
20 G-18235 GREHN GLAUCOMA PUNCH 20G	76
21 G-18238 KELLY DESCemet PUNCH for Glaucoma, 1.0 mm	76
22 G-18239 LUNTZ-DODICK TRABECULECTOMY PUNCH 360° rotatable	76
23 G-18791 DRESSING FORCEPS 1 x 2 teeth, 0.7 mm	110
24 G-18792 DRESSING FORCEPS 1 x 2 teeth, 1.0 mm	110
25 G-19030 TUEBINGEN TYING FORCEPS straight	106
26 G-18900 BONN IRIS FORCEPS straight, 1 x 2 teeth	110
27 G-16190 SPATULA angled, 1.0 mm wide	159
28 G-19230 SCHIOETZ TONOMETER	48
29 G-14105 MICRO ROUND KNIFE for viscocanalostomy	63
30 G-31487 DIAMOND KNIFE for viscocanalostomy	62
31 G-31470 TRI-FACET DIAMOND STEP KNIFE	61
32 G-15197 SINGLE-USE IRRIGATION CANNULA angled, 35G	226
33 G-32210 FRONIMOPOLUS MARKER U-shaped	44
34 G-33786 MARKING PEN gentian violet, 20 pcs./box	231
35 G-32209 TOBIAS NEUHANN SCLERAL FLAP MARKER	44
36 G-39526 STERISAFE DURO A4 PLUS COMPLETE SET	259


EMERGENCY



	page
1 G-15960 LIEBERMANN-TENNANT EYE SPECULUM	25
2 G-15995 COOK EYE SPECULUM medium 11 mm blades	23
3 G-17023 BARRAQUER EYE SPECULUM premat. infants	15
4 G-18181 HALSTEAD MOSQUITO FORCEPS curved	135
5 G-15185 FRANKFURT HYDRODISSECTION CANNULA 27G	192
6 G-15200 CRANZ EYE DOUCHE silver	191
7 G-15151 LACRIMAL CANNULA conical, short angle, steel	189
8 G-15141 LACRIMAL CANNULA conical, round curve, steel	189
9 G-34026 NANOEDGE PARACENTESIS KNIFE 15° straight	217
10 G-19750 WESTCOTT SCISSORS blunt tips, curved, 12 cm	89
11 G-19485 STEVENS TENOTOMY SCISSORS curved, long	89
12 G-19765 VANNAS CAPSULOTOMY SCISSORS curved	84
13 G-32920 JAMESON STRABISMUS HOOK 11 mm	174
14 G-15805 GRAEFE STRABISMUS HOOK medium, 10 mm	173
15 G-15799 GRAEFE STRABISMUS HOOK short, 5 mm	173
16 G-15821 BONN STRABISMUS HOOK 1.0 mm hole	174

	page
17 G-32120 DESMARRES LID RETRACTOR fig. 0	35
18 G-16060 DESMARRES LID RETRACTOR fig. 1	35
19 G-16065 DESMARRES LID RETRACTOR fig. 2	35
20 G-18782 DRESSING FORCEPS serrated 1.0 mm	110
21 G-19005 CASTROVIEJO SUTURING FORCEPS 0.3 mm	103
22 G-19010 CASTROVIEJO SUTURING FORCEPS 0.6 mm	103
23 G-19032 TYING FORCEPS angled, 7.5 mm	107
24 G-19030 TUEBINGEN TYING FORCEPS straight	106
25 G-18950 COLIBRI FORCEPS with tying platform	104
26 G-17500 BARRAQUER NEEDLE HOLDER curved, w/o lock	139
27 G-17565 CASTROVIEJO NEEDLE HOLDER straight, w. lock	142
28 G-19130 CASTROVIEJO CALIPER 15 mm	43
29 G-18850 SHAAF CILIA AND FOREIGN BODY FORCEPS	101
30 G-14550 FOREIGN BODY INSTRUMENT needle/spud	55


ENUCLEATION



	page
1 G-15990 COOK EYE SPECULUM large.....	23
2 G-18180 HALSTEAD MOSQUITO FORCEPS straight.....	135
3 G-18181 HALSTEAD MOSQUITO FORCEPS curved.....	135
4 G-19050 GRAEFE FIXATION FORCEPS without lock.....	109
5 G-18940 PAUFIQUE UNIVERSAL SUTURING FORCEPS.....	105
6 G-19375 EYE SCISSORS curved, sharp tips, 10 cm.....	91
7 G-19485 STEVENS TENOTOMY SCISSORS curved, long.....	89
8 G-15810 GRAEFE STRABISMUS HOOK large, 12 mm.....	173
9 G-19112 BANGERTER MUSCLE FORCEPS.....	126
10 G-15780 WELLS ENUCLEATION SPOON.....	167

	page
11 G-19525 ENUCLEATION SCISSORS medium curve, 13.5 cm.....	94
12 G-31545 BUNGE EXENTERATION CURETTE medium.....	52
13 G-31555 BUNGE EXENTERATION CURETTE large.....	52
14 G-17575 CASTROVIEJO NEEDLE HOLDER curved, w. lock.....	142
15 G-20416 WALSER EYE SPHERE 16 mm.....	245
16 G-20417 WALSER EYE SPHERE 17 mm.....	245
17 G-20418 WALSER EYE SPHERE 18 mm.....	245
18 G-20419 WALSER EYE SPHERE 19 mm.....	245
19 G-20420 WALSER EYE SPHERE 20 mm.....	245
20 G-39526 STERISAFE DURO A4 PLUS COMPLETE SET.....	259

PTERYGIUM



	page
1 G-15965 KERSHNER EYE SPECULUM fenestrated.....	25
2 G-32120 DESMARRES LID RETRCTOR fig. 0.....	35
3 G-19355 EYE SCISSORS straight, sharp tips, 10 cm.....	91
4 G-14010 KUHNT CORNEAL SCARIFIER fig. 1.....	64
5 G-17570 CASTROVIEJO NEEDLE HOLDER curved, w/o lock.....	142
6 G-19080 BARRAQUER FIXATION FORCEPS.....	109

	page
7 G-18940 PAUFIQUE UNIVERSAL SUTURING FORCEPS.....	105
8 G-18782 DRESSING FORCEPS serrated, 1.0 mm.....	110
9 G-18792 DRESSING FORCEPS 1 x 2 teeth, 1.0 mm.....	110
10 G-18300 DIAMOND BURR 8.0 mm.....	244
11 G-32700 SCALPEL HANDLE.....	231
12 G-32705 SCALPEL SURGICAL BLADES #11, curved, 100pcs.....	231


OFFICE SET (OPTOMETRIST)



	page
1 G-17020 BARRAQUER EYE SPECULUM fenestrated.....	15
2 G-17025 BARRAQUER EYE SPECULUM infants, fenestrated.....	15
3 G-18975 BARRAQUER SUTURE REMOVING FORCEPS	104
4 G-19765 VANNAS CAPSULOTOMY SCISSORS curved.....	84
5 G-18750 BARRAQUER CILIA FORCEPS	101
6 G-14550 FOREIGN BODY INSTRUMENT needle/spud.....	55
7 G-18850 SHAAF CILIA AND FOREIGN BODY FORCEPS	101
8 G-150401 WILDER LACRIMAL PROBE slim, long taper.....	183
9 G-32120 DESMARRES LID RETRACTOR fig. 0.....	35

	page
10 G-19460 STRABISMUS SCISSORS short, straight.....	91
11 G-32715 SCHOCKET SCLERAL DEPRESSOR.....	180
12 G-15870 MEYHOEFFER CHALAZION CURETTE 2 mm.....	51
13 G-18680 LAMBERT CHALAZION FORCEPS	134
14 G-31105 EYE SHIELD plastic, transparent, 100 pcs.....	245
15 G-15230 BONN IRRIGATION CANNULA delicate.....	194
16 G-32665 SILICONE BULB 20 ml.....	212
17 G-32920 JAMESON STRABISMUS HOOK 11 mm.....	174
18 G-39526 STERISAFE DURO A4 PLUS COMPLETE SET	259


CHALAZION



	page
1 G-17040 JAEGER LID PLATE stainless steel.....	36
2 G-19375 EYE SCISSORS curved, sharp tips, 10 cm.....	91
3 G-18791 DRESSING FORCEPS 1 x 2 teeth, 0.7 mm.....	110
4 G-18792 DRESSING FORCEPS 1 x 2 teeth, 1.0 mm.....	110
5 G-18685 DESMARRES CHALAZION FORCEPS 20 mm.....	134
6 G-18710 LORDAN CHALAZION FORCEPS triangular.....	134
7 G-15870 MEYHOEFFER CHALAZION CURETTE 2.0 mm.....	51

	page
8 G-15875 MEYHOEFFER CHALAZION CURETTE 2.5 mm.....	51
9 G-15880 MEYHOEFFER CHALAZION CURETTE 3.0 mm.....	51
10 G-15910 AXENFELD CHALAZION CURETTE	52
11 G-17575 CASTROVIEJO NEEDLE HOLDER curved, w. lock.....	142
12 G-32700 SCALPEL HANDLE	231
13 G-32705 SCALPEL SURGICAL BLADES #11, curved, 100 pcs.....	231


LID / PTOSIS



	page
1 G-17040 JAEGER LID PLATE stainless steel.....	36
2 G-18170 HARTMANN MOSQUITO FORCEPS straight.....	135
3 G-18171 HARTMANN MOSQUITO FORCEPS curved.....	135
4 G-19480 STEVENS TENOTOMY SCISSORS straight, long.....	89
5 G-19355 EYE SCISSORS straight, sharp tips, 10 cm.....	91
6 G-14010 KUHNT CORNEAL SCARIFIER fig. 1.....	64
7 G-19750 WESTCOTT SCISSORS blunt tips, curved, 12 cm.....	89
8 G-34026 NANOEDGE PARACENTESIS KNIFE 15° straight.....	217
9 G-19600 CORNEAL SCISSORS light curve, sharp tips.....	82
10 G-18715 BERKE PTOSIS FORCEPS.....	134
11 G-17575 CASTROVIEJO NEEDLE HOLDER curved, w. lock.....	142
12 G-19130 CASTROVIEJO CALIPER 15 mm.....	43
13 G-16060 DESMARRES LID RETRACTOR fig. 1.....	35
14 G-16065 DESMARRES LID RETRACTOR fig. 2.....	35
15 G-16070 DESMARRES LID RETRACTOR fig. 3.....	36
16 G-15800 GRAEFE STRABISMUS HOOK small, 8 mm.....	173
17 G-15805 GRAEFE STRABISMUS HOOK medium, 10 mm.....	173


	page
18 G-15355 GUTHRIE RETRACTOR sharp, large.....	36
19 G-15900 MEYHOEFER CHALAZION CURETTE 1.0 mm.....	51
20 G-15870 MEYHOEFER CHALAZION CURETTE 2.0 mm.....	51
21 G-15875 MEYHOEFER CHALAZION CURETTE 2.5 mm.....	51
22 G-19050 GRAEFE FIXATION FORCEPS without lock.....	109
23 G-18810 DRAEGER CONJUNCTIVA FORCEPS atraumatic.....	102
24 G-31953 IRIS FORCEPS straight, 1 x 2 teeth, strong.....	111
25 G-19005 CASTROVIEJO SUTURING FORCEPS 0.3 mm.....	103
26 G-19010 CASTROVIEJO SUTURING FORCEPS 0.6 mm.....	103
27 G-19040 TUEBINGEN TYING FORCEPS curved.....	106
28 G-18960 MICRO COLIBRI FORCEPS with tying platform.....	104
29 G-18710 LORDAN CHALAZION FORCEPS triangular.....	134
30 G-18680 LAMBERT CHALAZION FORCEPS.....	134
31 G-18685 DESMARRES CHALAZION FORCEPS 20 mm.....	134
32 G-15542 AXENFELD LACRIMAL SAC RETRACTOR.....	37
33 G-33786 MARKING PEN gentian violet, 20 pcs./box.....	231

DCR / LACRIMAL



page	page		
1 G-17070 MUELLER LACRIMAL SAC RETRATOR	38	18 G-150101 BOWMAN LACRIMAL PROBE size 0000/000	185
2 G-18180 HALSTEAD MOSQUITO FORCEPS straight	135	19 G-150102 BOWMAN LACRIMAL PROBE size 00/0	185
3 G-18181 HALSTEAD MOSQUITO FORCEPS curved	135	20 G-150103 BOWMAN LACRIMAL PROBE size 1/2	185
4 G-18225 KERRISON BONE NIBBLING RONGEUR	75	21 G-33198 JUENEMANN LACRIMAL PROBE CANNULA	190
5 G-32355 NASAL SPECULUM	242	22 G-33194 JUENEMANN LACRIMAL PROBE	186
6 G-31810 FRIEDMANN RONGEUR	242	0.8 mm, w/ grip plate	186
7 G-33125 MEAD HAMMER length 19 cm, 300 g	242	23 G-18782 DRESSING FORCEPS serrated, 1.0 mm	110
8 G-15020 GIESSEN LACRIMAL PROBE angled	183	24 G-18940 PAUFIQUE UNIVERSAL SUTURING FORCEPS	105
9 G-150401 WILDER LACRIMAL PROBE slim, long taper	183	25 G-18810 DRAEGER CONJUNCTIVA FORCEPS atraumatic	102
10 G-150302 WILDER LACRIMAL PROBE medium	183	26 G-34290 SINGLE-USE LACRIMAL CANNULA curved, 26G	224
11 G-18620 HEIDELBERG RASPATORY	242	27 G-15190 ANEL LACRIMAL CANNULA straight, 23G	189
12 G-18630 HEIDELBERG ELEVATORIUM	242	28 G-15195 ANEL LACRIMAL CANNULA round curve, 23G	189
13 G-18240 BONE CHISEL flat, 5.0 mm wide	243	29 G-19030 TUEBINGEN TYNG FORCEPS straight	106
14 G-18245 BONE CHISEL hollow, 5.0 mm wide	243	30 G-15542 AXENFELD LACRIMAL SAC RETRATOR	37
15 G-15095 PIGTAIL PROBE with suture holes	184	31 G-19485 STEVENS TENOTOMY SCISSORS curved, long	89
16 G-34812 LACRIMAL INTUBATION SET olive tips, 27G	224	32 G-19355 EYE SCISSORS straight, sharp tips, 10 cm	91
17 G-32726 MUENSTER SUTURE HOOK	172	33 G-17575 CASTROVIEJO NEEDLE HOLDER curved, w/ lock	142


STRABISMUS / MUSCLE



	page
1 G-15995 COOK EYE SPECULUM medium, 11 mm.....	23
2 G-16000 COOK EYE SPECULUM small, 8 mm.....	24
3 G-32128 COOK EYE SPECULUM infants, 5 mm.....	24
4 G-18180 HALSTEAD MOSQUITO FORCEPS straight.....	135
5 G-18181 HALSTEAD MOSQUITO FORCEPS curved.....	135
6 G-18190 DIEFFENBACH SERREFINE straight, 35 mm.....	136
7 G-19355 EYE SCISSORS straight, sharp tips, 10 cm.....	91
8 G-19475 STEVENS TENOTOMY SCISSORS curved, short.....	89
9 G-19600 CORNEAL SCISSORS light curve, sharp tips.....	82
10 G-19750 WESTCOTT SCISSORS blunt tips, curved, 12 cm.....	89
11 G-32225 VOGT SCLERAL RULER.....	44
12 G-16060 DESMARRES LID RETRACTOR fig. 1.....	35
13 G-16065 DESMARRES LID RETRACTOR fig. 2.....	35
14 G-15830 MALBRAN STRABISMUS HOOK.....	174
15 G-15800 GRAEFE STRABISMUS HOOK small, 8 mm.....	173


	page
16 G-15805 GRAEFE STRABISMUS HOOK medium, 10 mm.....	173
17 G-19135 CASTROVIEJO CALIPER 20 mm.....	43
18 G-15797 FINK MARKER for 9 mm recession.....	44
19 G-18782 DRESSING FORCEPS serrated, 1.0 mm.....	110
20 G-18810 DRAEGER CONJUNCTIVA FORCEPS atraumatic.....	102
21 G-19010 CASTROVIEJO SUTURING FORCEPS 0.6 mm.....	103
22 G-19005 CASTROVIEJO SUTURING FORCEPS 0.3 mm.....	103
23 G-19030 TUEBINGEN TYING FORCEPS straight.....	106
24 G-19040 TUEBINGEN TYING FORCEPS curved.....	106
25 G-19120 SATTLER MUSCLE FORCEPS w. slide lock, right.....	127
26 G-19125 SATTLER MUSCLE FORCEPS w. slide lock, left.....	127
27 G-19112 BANGERTER MUSCLE FORCEPS.....	126
28 G-19137 CUEPPERS STRABISMUS MUSCLE TUCKER.....	42
29 G-17575 CASTROVIEJO NEEDLE HOLDER curved, w. lock.....	142
30 G-17565 CASTROVIEJO NEEDLE HOLDER straight, w. lock.....	142

RETINAL DETACHMENT




	page		page
1 G-15966 KERSHNER EYE SPECULUM solid blades.....	25	18 G-15800 GRAEFE STRABISMUS HOOK small, 8 mm.....	173
2 G-18182 HALSTEAD MOSQUITO FORCEPS curved, 1x2 teeth.....	135	19 G-15805 GRAEFE STRABISMUS HOOK medium, 10 mm.....	173
3 G-33200 SCHAEDEL TOWEL FORCEPS.....	136	20 G-18950 COLIBRI FORCEPS with tying platform.....	104
4 G-18190 DIEFFENBACH SERREFINE straight, 35 mm.....	136	21 G-18782 DRESSING FORCEPS serrated, 1.0 mm.....	110
5 G-34026 NANOEDGE PARACENTESIS KNIFE 15° straight.....	217	22 G-18792 DRESSING FORCEPS 1 x 2 teeth, 1.0 mm.....	110
6 G-18820 REMKY SUPERIOR RECTUS FORCEPS.....	125	23 G-19005 CASTROVIEJO SUTURING FORCEPS 0.3 mm.....	103
7 G-17565 CASTROVIEJO NEEDLE HOLDER straight, w. lock.....	142	24 G-19120 SATTLER MUSCLE FORCEPS w slide lock, right.....	127
8 G-19750 WESTCOTT SCISSORS blunt tips, curved, 12 cm.....	89	25 G-19125 SATTLER MUSCLE FORCEPS w slide lock, left.....	127
9 G-19600 CORNEAL SCISSORS light curve, sharp tips.....	82	26 G-19032 TYING FORCEPS angled, 7.5 mm.....	107
10 G-19375 EYE SCISSORS curved, sharp tips, 10 cm.....	91	27 G-19030 TUEBINGEN TYING FORCEPS straight.....	106
11 G-19355 EYE SCISSORS straight, sharp tips, 10 cm.....	91	28 G-31953 IRIS FORCEPS straight, 1 x 2 teeth, strong.....	111
12 G-19475 STEVENS TENOTOMY SCISSORS curved, short.....	89	29 G-18855 IRIS FORCEPS straight, serrated.....	110
13 G-32445 HAMBURG ORBITAL SPATULA 18 x 12 mm.....	162	30 G-17500 BARRAQUER NEEDLE HOLDER curved, w/o lock.....	139
14 G-15720 BANGERTER SPATULA 6 mm wide.....	162	31 G-17501 BARRAQUER NEEDLE HOLDER curved, w. lock.....	139
15 G-15821 BONN STRABISMUS HOOK 1.0 mm hole.....	174	32 G-18650 RULER 150 mm, 0.5 mm increments.....	43
16 G-15830 MALBRAN STRABISMUS HOOK.....	174	33 G-19130 CASTROVIEJO CALIPER 15 mm.....	43
17 G-32150 LAQUA SCLERAL DEPRESSOR double-ended.....	179	34 G-33752 WATZKE SPREADING FORCEPS.....	127

VITRECTOMY 20G




	page		page
1 G-15966 KERSHNER EYE SPECULUM solid blades.....	25	29 G-33425 INFUSION CANNULA beveled, 5.0 mm.....	209
2 G-33200 SCHAEDEL TOWEL FORCEPS.....	136	30 G-32181 DUISBURG INJECTION CANNULA 5.0 mm.....	209
3 G-18190 DIEFFENBACH SERREFINE straight, 35 mm.....	136	31 G-37002 BACKFLUSH HANDPIECE with silicone chamber.....	206
4 G-18180 HALSTEAD MOSQUITO FORCEPS straight.....	135	32 G-34291 SINGLE-USE BACKFLUSH CANNULA 20G.....	227
5 G-18170 HARTMANN MOSQUITO FORCEPS straight.....	135	33 G-33050 BACKFLUSH CANNULA 20G.....	207
6 G-19375 EYE SCISSORS curved, sharp tips, 10 cm.....	91	34 G-34293 SINGLE-USE BACKFLUSH CANNULA 20G w brush.....	227
7 G-19485 STEVENS TENOTOMY SCISSORS curved, long.....	89	35 G-33075 BACKFLUSH CANNULA with silicone brush.....	207
8 G-19750 WESTCOTT SCISSORS blunt tips, curved, 12 cm.....	89	36 G-36000 KILP LENS 60 D flat, field of vision 20°.....	244
9 G-19760 VANNAS CAPSULOTOMY SCISSORS straight.....	84	37 G-33435 SCLERAL PLUG 20G / 0.9 mm.....	128
10 G-32446 ORBITAL SPATULA 18 x 8 mm.....	163	38 G-37551 BARTZ-SCHMIDT ILM-SCRAPER 20G.....	175
11 G-34181 NANOEDGE MVR STILETTO straight, 20G.....	217	39 G-37511 MEMBRANE PEELER 0.4 mm, recessed, angled, 20G.....	175
12 G-32920 JAMESON STRABISMUS HOOK 11 mm.....	174	40 G-33052 WEBER ASPIRATION CANNULA for viscous fluids.....	210
13 G-15821 BONN STRABISMUS HOOK 1.0 mm hole.....	174	41 G-42011 HATTENBACH / NIKOLIC UNO C TROCAR SYSTEM 20G.....	221
14 G-19130 CASTROVIEJO CALIPER 15 mm.....	43	42 G-40006 HANDLE UNO C 6 mm, titanium.....	218
15 G-37410 MARKING INSTRUMENT.....	43	43 G-46041 UNO C VITREOUS FORCEPS ENDGRIPPING 20G.....	219
16 G-18782 DRESSING FORCEPS serrated, 1.0 mm.....	110	44 G-42061 UNO C VITREOUS SCISSORS 20G, straight.....	218
17 G-18810 DRAEGER CONJUNCTIVA FORCEPS atraumatic.....	102	45 G-42043 UNO C VITREOUS FORCEPS CROCODILE DESIGN 20G.....	219
18 G-19010 CASTROVIEJO SUTURING FORCEPS 0.6 mm.....	103	46 G-42044 TANO UNO C FORCEPS 20G.....	219
19 G-18940 PAUFIQUE UNIVERSAL SUTURING FORCEPS.....	105	47 G-37425 FIXATION INSTRUMENT for trocar systems.....	150
20 G-19030 TUEBINGEN TYING FORCEPS straight.....	106	48 G-37430 TROCAR FORCEPS.....	129
21 G-18960 MICRO COLIBRI FORCEPS with tying platform.....	104	49 G-37435 HATTENBACH HORSESHOE-HOOK.....	175
22 G-17570 CASTROVIEJO NEEDLE HOLDER curved, w/o lock.....	142	50 G-39786 MICROSTOP® TRAY MINI low, with lid.....	249
23 G-36306 VITREORETINAL FORCEPS endgripping, 20G.....	130	51 G-39783 MICROSTOP® TRAY MINI low, without lid.....	249
24 G-36256 VITREORETINAL FORCEPS crocodile, 20G.....	130	52 G-39782 MICROSTOP® RETAINER ELEMENT low.....	250
25 G-36356 TANO VITREORETINAL FORCEPS serrated, 20G.....	130	53 G-39750 MICROSTOP® CONTAINER HIGH for two trays.....	249
26 G-36526 HEIMANN VITREORETINAL SCISSORS 20G, angled.....	95		
27 G-36532 HEIMANN VITREORETINAL SCISSORS 20G, straight.....	95		
28 G-31900 KUHN INTRAOCULAR MAGNET.....	246		

VITRECTOMY 23G



	page		page
1 G-15966 KERSHNER EYE SPECULUM solid blades.....	25	26 G-36540 HEIMANN VITREORETINAL SCISSORS	
2 G-33200 SCHAEDEL TOWEL FORCEPS.....	136	23G, angled, for trocars.....	96
3 G-18190 DIEFFENBACH SERREFINE straight, 35 mm.....	136	27 G-37002 BACKFLUSH HANDPIECE with silicone chamber.....	206
4 G-18180 HALSTEAD MOSQUITO FORCEPS straight.....	135	28 G-34296 SINGLE-USE BACKFLUSH CANNULA 23G.....	227
5 G-18170 HARTMANN MOSQUITO FORCEPS straight.....	135	29 G-34294 SINGLE-USE BACKFLUSH CANNULA 23 G w brush.....	227
6 G-19375 EYE SCISSORS curved, sharp tips, 10 cm.....	91	30 G-34285 SINGLE-USE PFCL CANNULA coaxial, 23G.....	228
7 G-19485 STEVENS TENOTOMY SCISSORS curved, long.....	89	31 G-36000 KILP LENS 60 D flat, field of vision 20°.....	244
8 G-19750 WESTCOTT SCISSORS blunt tips, curved, 12 cm.....	89	32 G-37512 MEMBRANE PEELER 0.4 mm, recessed, angled, 23G.....	175
9 G-19760 VANNAS CAPSULOTOMY SCISSORS straight.....	84	33 G-34498 SINGLE USE CANNULA for silicone oil, 23G.....	228
10 G-32446 ORBITAL SPATULA 18 x 8 mm.....	163	34 G-42311 HATTENBACH / NIKOLIC UNO C TROCAR SYSTEM	
11 G-32920 JAMESON STRABISMUS HOOK 11 mm.....	174	23G.....	221
12 G-15821 BONN STRABISMUS HOOK 1.0 mm hole.....	174	35 G-40006 HANDLE UNO C 6 mm, titanium.....	218
13 G-19130 CASTROVIEJO CALIPER 15 mm.....	43	36 G-42361 UNO C VITREOUS SCISSORS 23G, straight.....	218
14 G-37410 MARKING INSTRUMENT.....	43	37 G-46341 UNO C VITREOUS FORCEPS ENDGRIPPING 23G.....	219
15 G-18782 DRESSING FORCEPS serrated, 1.0 mm.....	110	38 G-46343 UNO C VITREOUS FORCEPS CROCODILE DESIGN	
16 G-18810 DRAEGER CONJUNCTIVA FORCEPS atraumatic.....	102	23G.....	219
17 G-19010 CASTROVIEJO SUTURING FORCEPS 0.6 mm.....	103	39 G-42344 TANO UNO C FORCEPS 23G.....	219
18 G-18940 PAUFIQUE UNIVERSAL SUTURING FORCEPS.....	105	40 G-37425 FIXATION INSTRUMENT for trocar systems.....	150
19 G-19030 TUEBINGEN TYING FORCEPS straight.....	106	41 G-37430 TROCAR FORCEPS.....	129
20 G-18960 MICRO COLIBRI FORCEPS with tying platform.....	104	42 G-37435 HATTENBACH HORSESHOE-HOOK.....	175
21 G-17570 CASTROVIEJO NEEDLE HOLDER curved, w/o lock.....	142	43 G-39786 MICROSTOP® TRAY MINI low, with lid.....	249
22 G-36282 VITREORETINAL FORCEPS endgripping, 23G.....	132	44 G-39783 MICROSTOP® TRAY MINI low, without lid.....	249
23 G-36342 VITREORETINAL FORCEPS crocodile, 23G.....	132	45 G-39782 MICROSTOP® RETAINER ELEMENT low.....	250
24 G-36362 TANO VITREORETINAL FORCEPS serrated, 23G.....	132	46 G-39750 MICROSTOP® CONTAINER HIGH for two trays.....	249
25 G-36576 HEIMANN VITREORETINAL SCISSORS 23G, straight.....	96		

VITRECTOMY 25G



	page
1 G-15966 KERSHNER EYE SPECULUM solid blades.....	25
2 G-33200 SCHAEDEL TOWEL FORCEPS.....	136
3 G-18190 DIEFFENBACH SERREFINE straight, 35 mm.....	136
4 G-18180 HALSTEAD MOSQUITO FORCEPS straight.....	135
5 G-18170 HARTMANN MOSQUITO FORCEPS straight.....	135
6 G-19375 EYE SCISSORS curved, sharp tips, 10 cm.....	91
7 G-19485 STEVENS TENOTOMY SCISSORS curved, long.....	89
8 G-19750 WESTCOTT SCISSORS blunt tips, curved, 12 cm.....	89
9 G-19760 VANNAS CAPSULOTOMY SCISSORS straight.....	84
10 G-32446 ORBITAL SPATULA 18 x 8 mm.....	163
11 G-32920 JAMESON STRABISMUS HOOK 11 mm.....	174
12 G-15821 BONN STRABISMUS HOOK 1.0 mm hole.....	174
13 G-19130 CASTROVIEJO CALIPER 15 mm.....	43
14 G-37410 MARKING INSTRUMENT.....	43
15 G-18782 DRESSING FORCEPS serrated, 1.0 mm.....	110
16 G-18810 DRAEGER CONJUNCTIVA FORCEPS atraumatic.....	102
17 G-19010 CASTROVIEJO SUTURING FORCEPS 0.6 mm.....	103
18 G-18940 PAUFIQUE UNIVERSAL SUTURING FORCEPS.....	105
19 G-19030 TÜEBINGEN TYNG FORCEPS straight.....	106
20 G-18960 MICRO COLIBRI FORCEPS with tying platform.....	104
21 G-17570 CASTROVIEJO NEEDLE HOLDER curved, w/o lock.....	142

	page
22 G-36312 VITREORETINAL FORCEPS endgripping, 25G.....	132
23 G-36322 VITREORETINAL FORCEPS crocodile, 25G.....	132
24 G-36366 TANO VITREORETINAL FORCEPS serrated, 25G.....	132
25 G-36586 HEIMANN VITREORETINAL SCISSORS straight, 25G.....	96
26 G-37002 BACKFLUSH HANDPIECE with silicone chamber.....	206
27 G-34299 SINGLE-USE BACKFLUSH CANNULA 25G.....	227
28 G-34297 SINGLE-USE BACKFLUSH CANNULA 25G w brush.....	227
29 G-36000 KILP LENS 60 D flat, field of vision 20°.....	244
30 G-37513 MEMBRANE PEELER 0.4 mm, recessed, angled, 25G.....	175
31 G-40006 HANDLE UNO C 6 mm, titanium.....	218
32 G-42561 UNO C VITREOUS SCISSORS 25G, straight.....	218
33 G-46541 UNO C VITREOUS FORCEPS ENDGRIPPING 25G.....	219
34 G-46542 HATTENBACH UNO C VITREOUS FORCEPS 25G wave design.....	219
35 G-42544 TANO UNO C FORCEPS 25G.....	219
36 G-37425 FIXATION INSTRUMENT for trocar systems.....	150
37 G-37435 HATTENBACH HORSESHOE-HOOK.....	175
38 G-39786 MICROSTOP® TRAY MINI low, with lid.....	249
39 G-39783 MICROSTOP® TRAY MINI low, without lid.....	249
40 G-39782 MICROSTOP® RETAINER ELEMENT low.....	250
41 G-39750 MICROSTOP® CONTAINER HIGH for two trays.....	249

IN

INDEX

Geuder®
Precision made in Germany

NUMERICAL INDEX

ITEM NO.	PAGE	ITEM NO.	PAGE	ITEM NO.	PAGE	ITEM NO.	PAGE
G-12...		15189	193	15750	162	16070	36, S16
12200	66	15190	189, S17	15755	166	16075	36
12220	66	15195	189, S17	15757	167	16080	35
		15197	, 226	15758	163, S9	16105	191
		15200	191, S3, S13	15780	167, S14	16106	194, S2, S4
G-14...		15220	191	15797	44, S18	16115	193, S2
		15225	195	15799	173, S13	16120	195
14000	63	15226	195	15800	173, S16, S18-S19	16125	195, S2
14010	64, S14, S16	15230	194	15805	173, S13, S16, S18-S19	16130	195, S6
14020	64, S9	15240	194	15810	173, S14	16132	195
14105	63, S12	15250	192, S3-S5, S9, S12	15820	174	16133	195
14162	65	15271	190	15821	174, S13, S19,S21, S23-S24	16134	194
14165	65	15272	190	15830	174, S18, S19	16135	194
14180	63	15273	190	15840	52	16136	194
14185	63	15275	192	15845	52	16138	194
14200	64, S9	15310	197, S10	15870	51, S15, S16	16161	156
14420	63	15312	197, S10	15875	51, S15, S16	16166	156
14480	56	15314	197	15880	51, S15	16167	156, S2-S4
14490	56	15315	186	15885	51	16169	157
14500	56	15350	38	15900	51, S16	16175	154, 157, S2
14530	55, 266	15355	38, S16	15910	52, S15	16181	154
14540	55, 266	15360	166	15920	31	16182	154
14550	55, 266, S13, S15	15365	166	15921	31	16185	155
14570	56	15370	166	15925	31	16189	153, 159, S2-S3
		15375	191	15926	31	16190	159, S12
G-15...		15376	191	15928	32	16191	158
		15377	191	15945	20	16193	158
15020	183, S17	15390	165	15946	20	16200	65
15050	184	15391	165	15950	21	16205	158
15055	184	15392	48	15951	21	16210	157
15065	150, S9	15393	48	15955	21	G-17...	
15070	150, S9	15400	165	15956	22		
15075	150, S9	15405	165	15960	25, S2-S4, S7-S8,S10, S12-S13	17000	24
15080	150, S9	15470	164	15961	28, S7, S8	17001	24
15085	150, S9	15475	164	15962	29	17002	25
15090	150, S9	15477	163	15965	25, S14	17003	25
15095	184, S17	15480	164	15966	25, S9, S19, S21,S23-S24	17004	17
15096	184	15485	164	15970	22	17005	27
15105	186	15540	37	15975	22	17006	27
15110	186	15542	37, S16-S17	15990	23, S14	17010	18
15125	186, S12	15545	37	15995	23, S13, S18	17011	18
15130	186, S12	15550	37			17015	18
15134	104	15555	37			17016	18
15136	186	15560	37	G-16...		17017	18
15137	186	15620	161, S2, S4	16000	24, S18	17019	27
15138	186	15630	161	16010	23	17020	15, S2, S12, S15
15140	189	15635	161	16020	23	17021	27
15141	189, S13	15640	161	16025	23	17023	15, S13
15150	189	15645	161	16030	23	17024	28
15151	189, S3, S13	15710	180	16050	21	17025	15
15160	189	15720	162, S19	16055	21	17026	15
15161	189	15725	162	16060	35, S13, S16, S18	17027	19
15170	190	15730	162	16065	35, S13, S16, S18	17030	19
15180	194, S7	15731	162			17031	29
15185	192, S2, S12-S13	15740	162				
15188	193	15745	162				

ITEM NO.	PAGE	ITEM NO.	PAGE	ITEM NO.	PAGE	ITEM NO.	PAGE
17040	36, S15-S16	18160	136		S23-S24	19026	103
17070	38, S17	18170	135, S16, S21, S23-S24	18820	125, S2, S12, S19	19030	106, S2-S3, S9-S10,
17080	38	18171	135, S16	18840	111		S12-S13, S17-S19,
17085	38	18180	135, S2, S9, S12, S14,	18845	111		S21, S23-S24
17250	69		S17-S18, S21, S23-S24	18846	112	19032	107, S2-S3, S13, S19
17251	69	18181	135, S13-S14, S17-S18	18847	102	19034	106
17260	69	18182	135, S19	18848	102	19035	106
17261	69	18185	135	18850	101, S13, S15	19038	107
17270	69	18190	136, S2, S18-S19,	18855	110	19040	106, S9-S10, S16, S18
17271	69		S21, S23-S24	18856	110	19045	105
17280	69	18195	136	18857	111	19046	105
17281	69	18197	136	18860	111	19048	107
17290	69	18210	75	18870	111	19050	109, S14, S16
17291	69	18220	75	18880	112	19055	109
17300	69	18225	75, S17	18890	111	19060	109
17440	88, 143	18232	77	18897	112	19080	109, S14
17441	88, 143	18235	76	18900	110, S12	19105	126
17450	139, S9	18237	77	18905	110	19107	127
17451	139	18238	76, S12	18906	110	19112	126, S14, S18
17452	139	18239	76, S12	18908	110	19120	127, S18-S19
17453	139	18240	243, S17	18910	110, S2	19125	127, S18-S19
17460	139	18245	243, S17	18911	110	19130	43, S2-S3, S12-S13,
17461	139	18300	244, S14	18913	112		S16, S19, S21, S23-S24
17462	139	18301	244	18914	119, S7-S8	19135	43, S18
17463	139	18302	244	18918	122	19137	42, 172, S18
17500	139, S2-S4, S9, S12-S13, S19	18620	242, S17	18919	122	19138	43
		18630	242, S17	18928	123	19230	48, S12
17501	139, S19	18645	44	18930	122	19270	145
17505	139, S10	18650	43, S19	18935	121, S2, S4	19300	145
17506	139	18670	134	18940	105, S14, S17,	19320	145
17510	140	18680	134, S15-S16		S21, S23-S24	19330	145
17511	140	18685	134, S15-S16	18941	125	19345	62
17515	140	18690	134	18944	123	19350	91
17516	140	18695	134	18945	105	19355	91, S14, S16-S19
17520	141	18700	134	18950	104, S4, S13, S19	19360	91
17525	141	18710	134, S15-S16	18960	104, S2-S3, S9, S16,	19365	91
17530	140	18715	134, S16		S21, S23-S24	19370	91
17535	140	18720	101	18965	104	19375	91, S14-S15, S19,
17540	142	18730	101	18966	104		S21, S23-S24
17545	142	18740	101	18970	105	19380	91
17560	142	18750	101	18972	104, S9	19385	91
17565	142, S12-S13, S18-S19	18780	110	18975	104	19390	92
17570	142, S14, S21, S23-S24	18781	110			19395	92
17575	142, S2, S14-S18	18782	110, S13-S14,	G-19...		19400	92, S2-S4, S9, S12
17610	143		S17-S19, S21, S23-S24			19405	92
17615	143	18783	110	19001	103, S9	19410	92
17640	144	18784	110	19005	103, S13, S16, S18-S19	19420	90
17645	144	18790	110	19008	103	19425	90
17650	144	18791	110, S2, S9, S12, S15	19010	103, S9, S13, S16,	19430	90
17660	144	18792	110, S12, S14, S15, S19		S18, S21, S23-S24	19435	90
17665	144	18793	110	19015	103	19450	93
		18794	110	19020	103	19460	91
G-18...		18795	126	19022	103	19465	91
		18800	126	19023	103, S3	19470	89
18127	145, 231, S7-S8	18805	117	19024	103	19475	89, S18-S19
18150	145, S7-S8	18810	102, S16-S18, S21,	19025	103	19480	89, S16

NUMERICAL INDEX

ITEM NO.	PAGE	ITEM NO.	PAGE	ITEM NO.	PAGE	ITEM NO.	PAGE
19485	89, S3, S13-S14, S17, S21, S23-S24	19856	253	20618	69	G-31...	
19505	93	19865	257	G-21...		31013	264
19510	93	19866	257	31014		31014	266
19520	94	19867	257	21175	72	31016	265
19525	94, S14	19870	256	21176	72	31017	265
19530	94	19876	253, 255	G-22...		31018	265
19540	86	19877	254-255	31021		31021	264
19545	86	19890	255	31022		31022	265
19550	86	19892	254-255	22100	198, S3, S5	31023	265
19570	82	19894	255	22101	199, S3, S5	31024	265
19580	82	19899	257	22102	198	31025	264
19585	82	19960	212	22105	199	31026	264
19590	83	19961	212	22106	199	31027	265
19600	82, S16, S18-S19	G-20...		22107	199	31028	264
19605	82	20050	69	22525	205	31029	264
19607	82	20054	69	22527	205	31031	265
19615	82	20055	69	22540	205	31034	266
19620	84	20060	69	22542	205	31105	245
19630	84, S12	20065	69	22545	205	31230	102
19640	86	20070	69	22550	197	31235	102
19655	86	20080	69	22556	205	31299	114
19660	87	20090	69	22615	205	31300	113
19675	87	20100	69	22619	193, S12	31301	114
19680	87, S12	20216	245	22632	198	31303	114
19690	81	20217	245	22633	198	31304	113
19695	81	20218	245	22636	264	31305	114, S2-S3
19710	81, S2	20219	245	22637	199	31306	115
19715	81, S2	20220	245	G-23...		31307	113, S4
19718	81	20316	245	31308		31308	113
19720	81, S9-S10	20317	245	23991	231	31315	94
19725	81, S9-S10	20318	245	G-24...		31355	192
19727	82	20319	245	31378		31378	160
19740	83	20320	245	31415		31415	136
19745	83	20416	245, S14	31420		31420	136
19747	83	20417	245, S14	24100	204-205	31468	61
19750	89, S2, S9, S12-S13, S16, S18-S19, S21, S23-S24	20418	245, S14	24311	203	31470	61, S12
19751	89	20419	245, S14	24314	204	31477	62
19760	84, S21, S23-S24	20420	245, S14	24318	203	31480	62
19765	84, S2-S4, S9, S13, S15	20450	69	24330	204	31487	62, S12
19770	84, S10	20455	69	24332	204	31490	149
19775	85	20460	69	24418	203	31491	149
19776	83, S2	20465	69	24430	204	31500	151
19777	83	20470	69	24432	204	31501	151
19780	81	20480	69	24435	204	31502	151
19790	81	20490	69	24440	204	31503	151
19820	246	20500	69	G-28...		31504	151
19822	252, 255	20610	69	31505		31505	152
19823	252, 255	20611	69	31506		31506	152
19825	246	20612	69, S12	28150	252, 255	31507	152-154, S3, S5
19841	255, 257	20613	69	28151	252, 255	31510	193
19842	254-255, 257	20614	69	28752	211	31518	120
19843	253, 255, 257	20615	69	28766	230	31519	120
19846	257	20616	69	28767	230	31520	229
19847	257	20617	69	28768	230	31545	52, S14
						31555	52, S14

ITEM NO.	PAGE	ITEM NO.	PAGE	ITEM NO.	PAGE	ITEM NO.	PAGE
31570	128	32005	159	32224	42	32714	180
31620	101	32007	159	32225	44, S18	32715	180
31622	108, S4	32010	198	32227	158	32724	160
31623	108	32011	200, S5	32228	42	32726	172, S17
31624	108	32012	200, S5	32229	42, S9	32727	163
31625	107	32014	202	32240	112	32728	159, S10
31649	255	32015	197	32250	89	32729	166
31650	255	32016	201	32260	140	32740	101
31651	255	32017	201	32261	140	32759	193
31652	255	32018	200	32275	157	32760	194
31656	255	32019	200	32276	157	32765	194
31657	258	32023	84	32282	152	32769	196
31659	258	32026	201	32313	243	32771	198
31661	255	32027	179	32315	141	32772	196
31670	263	32028	201	32318	45, S9	32773	198
31720	120, S3	32034	155	32320	141	32774	196
31722	120	32059	202	32321	141	32776	196
31724	121	32100	83	32335	126	32778	198
31727	154	32101	83	32350	127	32785	256
31805	206	32105	83	32355	242, S17	32800	88, S4
31806	206	32115	36	32445	162, S19	32810	88
31807	193	32119	19	32446	163, S21, S23-S24	32811	88
31810	242, S17	32120	35, S13-S15	32500	152-153	32823	125, S4
31815	69	32126	17	32501	153	32905	174
31820	69	32127	17	32505	151, 153-154, S5	32915	173
31825	69	32128	24, S18	32523	159	32920	174, S13, S15,
31830	69	32129	17	32527	195		S21, S23-S24
31835	69	32130	20	32528	196	32926	115, S3, S5
31840	69	32133	17	32529	195	32928	115
31845	69	32134	19	32530	156	32935	115, S4
31850	69	32135	20	32550	155	32936	116
31855	69	32136	16, S4	32555	156	32937	116
31860	69	32137	16	32556	155	32938	116
31865	69	32138	18	32557	156	32940	85, S4
31870	70	32139	18	32563	160	32942	85
31891	229	32148	16, S3	32564	160, S6	32944	117, 264, S4
31895	192	32149	16	32600	191	32945	87, S4
31899	119	32150	179, S19	32605	191	32946	87
31900	246, S21	32152	179	32614	90	32947	108, 264, S4, S6
31901	124, S3	32153	179	32660	256	32948	86, 264, S4, S6
31902	123	32155	192	32665	212	32949	108, 143, S4, S6
31903	123	32160	153, 155	32670	212	32955	241
31905	124, S4	32161	153	32671	71	32960	241, S4
31908	124	32162	154, S3	32672	71	32965	241
31911	158	32164	153	32674	71	32970	241, S4
31953	111, S16, S19	32165	121	32675	256	32996	88, S4
31962	122	32166	192	32676	257	32997	125, S4, S6
31975	166	32167	192, S11	32677	257	32998	88, S6
31980	166	32168	153	32695	258		
31985	166	32169	154	32696	211, 229	G-33...	
31990	165	32180	209	32697	211		
31994	124, S4	32181	209, S21	32698	211	33022	64
		32182	209	32699	211	33025	64
G-32...		32209	44, S12	32700	231, S14-S15	33026	64
		32210	44, S12	32705	231, S14-S15	33045	206
32004	44	32222	41	32710	231	33050	207, S21

NUMERICAL INDEX

ITEM NO.	PAGE	ITEM NO.	PAGE	ITEM NO.	PAGE	ITEM NO.	PAGE
33051	210	33450	209	33873	70	G-34...	
33052	210, S21	33451	209	33874	70	34026	217, S2-S3, S13, S16
33053	207	33452	209	33875	70	34031	217
33054	207	33462	228	33876	70	34036	217, S12
33055	210, 265	33463	228	33877	70	34050	217
33056	210	33468	228	33878	70	34051	217, S3
33057	210, 266	33469	228	33879	70	34052	217, S3
33059	266	33470	210	33880	70	34053	217
33065	211	33471	210	33881	70	34054	217
33066	211	33472	210	33882	70	34055	217, S3
33070	207	33473	210	33883	70	34060	217
33071	207	33474	210	33884	70	34061	217
33072	207	33475	210	33885	70	34062	217
33073	207	33482	228	33886	70	34064	217
33075	207, S21	33484	209	33887	70	34081	217
33077	207	33485	209	33888	70	34086	217
33080	258	33486	209	33889	70	34122	217
33090	258	33488	229	33890	70	34126	217
33100	258	33489	229	33891	70	34146	217, S2
33125	242, S17	33494	209	33892	70	34181	217, S21
33135	189	33495	209	33893	70	34191	217
33140	189	33496	209	33894	70	34192	217
33141	265	33650	243, S7	33895	70	34196	217
33142	265	33652	243, S7	33896	70	34197	217
33143	265	33655	243, S7	33897	70	34198	217
33160	190	33657	243, S7	33898	70	34200	224
33170	190	33670	64, S7	33899	70	34205	224
33175	190	33725	41	33909	196, S7	34210	224, S3
33180	190	33726	41	33910	197, S7-S8	34215	224
33185	190	33727	41, S9	33911	197, S7	34220	224
33190	190	33740	41	33912	167, S7-S8	34225	224, S2
33193	186	33741	41	33914	30	34230	224
33194	186, S17	33742	41	33919	26	34235	225
33195	186	33745	45, S7	33920	26, S7	34240	225
33197	190	33746	45, S7	33925	30	34245	225, S5
33198	190, S17	33747	45, S7	33930	167, S8	34255	226
33199	193	33752	127, S19	33935	168	34261	225
33200	136, S12, S19, S21, S23-S24	33755	149	33940	170, S7	34265	226
		33757	150	33941	108, S8	34270	225
33416	191, S10	33758	150, S3, S8	33942	170	34271	225
33420	209	33759	150, S7-S8	33943	171, S7-S8	34275	227
33422	191	33760	46	33944	172	34280	226
33425	209, S21	33762	46	33950	169, S8	34285	228, S23
33427	171, S8	33763	46, S6	33952	170, S8	34289	227
33428	118, S7-S8	33774	46, S6	33954	170	34290	224, S17
33429	118, S8	33776	47, S6	33958	171	34291	227, S21
33430	209	33777	47, S6	33959	171, S8	34292	228
33432	128	33778	47	33960	119, S8	34293	227, S21
33434	127	33780	42	33961	119, S8	34294	227, S23
33435	128, S21	33781	42	33965	26, S5	34295	228
33436	128	33785	45	33966	169, S5	34296	227, S23
33437T	128	33786	231, S6-S7, S12, S16	33967	169, S5	34297	227, S24
33438	128	33831	129	33968	169, S5	34298	228
33439	128	33870	70	33971	26	34299	227, S24
33440	121	33871	70			34335	226
33441	128	33872	70				

ITEM NO.	PAGE	ITEM NO.	PAGE	ITEM NO.	PAGE	ITEM NO.	PAGE
34340	223	36090	259	36504	95	36955	97
34400	227	36202	130	36506	95	36975	97
34480	229	36206	130	36508	95	36985	97
34492	229	36208	130	36512	95		
34497	228	36212	130	36514	95		
34498	228, S23	36214	130	36516	95		
34523	223, S2, S4	36216	130	36518	95	37000	206
34602	232	36218	130	36522	95	37002	206, S21, S23-S24
34610	232	36222	131	36524	95	37005	206
34616	232	36226	131	36526	95, S21	37006	206
34618	232	36232	131	36528	95	37010	206, 208, 265
34620	233	36236	131	36532	95, S21	37011	208
34622	233	36242	131	36534	95	37012	208
34624	233	36246	131	36536	95	37015	208
34626	233	36248	131	36538	95	37020	208
34632	234	36252	131	36540	96, S23	37025	208
34634	234	36254	131	36541	96	37035	208
34636	234	36256	130, S21	36546	95	37046	208
34642	234	36258	130	36548	95	37065	208
34644	234	36262	131	36552	95	37080	206, 209, 265
34645	235	36264	131	36554	95	37100	265
34646	234	36266	131	36556	95	37105	265
34647	235	36272	131	36558	95	37130	22
34650	235	36276	130	36562	96	37390	210
34654	235	36278	130	36564	96	37400	43
34656	235	36282	132, S23	36566	96	37410	43, S21, S23-S24
34658	236	36284	132	36568	96	37412	43
34660	236	36286	132	36572	95	37414	43
34664	236	36288	132	36574	95	37418	43
34666	236	36292	130	36576	96, S23	37420	129
34668	237	36294	130	36578	96	37425	150, S21, S23-S24
34672	237	36296	130	36582	96	37430	129, S21, S23
34674	237	36298	130	36584	96	37435	175, S21, S23-S24
34684	237	36302	130	36586	96, S24	37440	149
34686	238	36304	130	36588	96	37445	180
34688	238	36306	130, S21	36600	133	37450	172
34690	238	36308	130	36601	133	37500	175
34710	238	36312	132, S24	36604	97, 266	37503	175
34806	206, 227	36314	132	36655	133	37505	175
34810	223	36316	130	36680	133	37506	175
34811	223	36318	130	36690	133	37511	175, S21
34812	224, S17	36322	132, S24	36695	133	37512	175, S23
34814	224	36324	132	36700	133	37513	175, S24
		36326	131	36705	133	37515	175
G-36...		36328	131	36710	133	37525	175
		36332	130	36715	133	37530	175
36000	244-245, S21, S23-S24	36336	131	36720	133	37535	175
36010	244, 245	36342	132, S23	36740	133	37545	175
36015	244	36344	132	36755	133	37551	175, S21
36020	244	36356	130, S21	36765	133	37555	65
36025	244	36358	130	36900	97	37570	175
36030	244	36362	132, S23	36901	97	37602	184
36035	245	36364	132	36902	97	37604	184
36045	244	36366	132, S24	36915	97	37700	109
36050	245	36368	132	36925	97	37705	109
36065	245	36502	95	36930	97	37710	109

NUMERICAL INDEX

ITEM NO.	PAGE	ITEM NO.	PAGE	ITEM NO.	PAGE	ITEM NO.	PAGE
G-38...		39530	261	42301	222	G-61...	
		39535	261	42302	222		
38215	86	39540	262	42311	221, S23	61201	220
38232	86	39541	260	42319	221	61202	220
38236	86	39542	260	42321	221	61204	220
38238	86	39543	260	42322	221	61205	220
38240	85, 116, S5	39545	261	42323	222	61231	220
38242	116	39546	261	42341	219	61232	220
38244	116, S5	39547	261	42342	219	61234	220
38246	85, S5	39548	260	42343	219	61235	220
38341	59	39549	260	42344	219, S23	61251	220
38342	60	39650	257	42361	218, S23	61252	220
38343	59	39680	251	42362	218	61254	220
38344	59	39741	250	42501	222	61255	220
38444	60	39745	249	42502	222		
38448	60, S5	39746	250	42511	221	G-80...	
38453	61	39750	249, S2-S3, S21,	42521	221		
38466	61		S23-S24	42522	221	80975	226
38473	60	39753	251	42523	222	G-150...	
38475	59	39754	251	42541	219		
38600	41, S10-S11	39755	251	42542	219	150001	185
38601	202, S10-S11	39762	251	42544	219, S24	150002	185
38602	202, S10-S11	39764	251	42561	218, S24	150003	185
38603	118, S10	39782	250, S2-S3, S21,	42562	218	150004	185
38604	168, S10		S23-S24			G-46...	
38606	118, 168, S10	39783	249, 250, S2-S3, S21,	150005	185		
38607	168, S10-S11		S23, S24	150051	185		
38608	168, S10-S11	39784	249	46001	222		
38620	71, S11	39785	249, 250	46002	222		
38621	64, S11	39786	249, S2-S3, S21,	46011	221		
38622	168, S11		S23-S24	46021	221		
38623	105, S11	39787	250	46022	221		
38624	117, S11	39788	250	46023	222		
38628	117	39789	250	46041	219, S21		
38630	223	39791	250	46042	219		
38631	203			46301	222		
38632	203	G-40...		46302	222		
38634	117, S10-S11			46311	221		
38635	223, S11	40006	218, S21, S23	46321	221		
		40008	218	46322	221		
				46323	222		
G-39...		G-42...		46341	219, S23		
				46342	219		
39501	261			46343	219, S23		
39502	262	42001	222	46362	218		
39503	262	42002	222	46501	222		
39504	262	42009	222	46502	222		
39505	263	42011	221, S21	46511	221		
39511	261	42019	221	46521	221		
39512	262	42021	221	46522	221		
39513	262	42022	221	46523	222		
39514	262	42023	222	46541	219, S24		
39525	259	42041	219	46542	219, S24		
39526	259, S4-S5, S7,S12, S14-S15	42042	219	46562	218		
39527	259	42043	219, S21				
39528	259	42044	219, S21				
		42061	218, S21				

ITEM NO. PAGE ITEM NO. PAGE

520550	215	560625	216
520600	215	560650	216
520625	215	560675	216
520650	215	560700	216
520675	215	560725	216
520700	215	560750	216
520725	215	560775	216
520750	215	560800	216, S9
520775	215	560825	216
520800	215	560850	216, S9
520825	215	560875	216
520850	215	560900	216
520875	215	560925	216
520900	215	560950	216
520950	215		

G-570 ...**G-521 ...**

521000	215	570600	216
521050	215	570650	216
521100	215	570700	216
521150	215	570725	216
521150	215	570750	216, S9
521200	215	570775	216
521250	215	570800	216, S9
521300	215	570825	216
521350	215	570850	216
521400	215	570875	216
521450	215	570900	216
521500	215		
521550	215		
521600	215		
521650	215		
521700	215		
521800	215		

G-540 ...

540600	215
540625	215
540650	215
540675	215
540700	215
540725	215
540750	215
540775	215
540800	215
540825	215
540850	215
540875	215
540900	215
540925	215
540950	215

G-560 ...

560600	216
--------------	-----

KEYWORD	PAGE	KEYWORD	PAGE	KEYWORD	PAGE
A					
Absorbent Swabs.....	231, S8, 231	Bone Gouge, West.....	243	Chopper, Tobias Neumann.....	152
A.C. Infusion Cannula.....	S10	Bone Nibbling Rongeur.....	75, S17	Cilia And Foreign Body Forceps..	101, S13, S15
A.C. Infusion Cannula, Bopp.....	S10	Bone Nibbling Rongeur, Beyer.....	75	Cilia And Foreign Body Forceps, Hamburg	
Adapter.....	264	Bone Nibbling Rongeur, Citelli.....	75	Model.....	101
Anterior Chamber Forceps.....	121	Bone Nibbling Rongeur, Kerrison	75, S17	Cilia and Foreign Body Forceps, Shaaf.....	S13
Anterior Chamber Forceps, von Barsewisch	121	Bonn Injection Set.....	238	Cilia And Foreign Body Forceps, Shaaf.....	101
Anterior Chamber Infusion Cannula.....	191	C		Cilia And Suture Forceps.....	101
Anterior Chamber Infusion Cannula, Bonn Model	191	Caliper	43, S2-S3, S12-S13, S16,	Cilia And Suture Forceps, Barraquer.....	101
Anterior Chamber Infusion Cannula, Bopp.....	191S18-S19, S21, S23-S24		Cilia Forceps.....	101, S15
Anterior Synechia Scissors.....	83	Caliper, Castroviejo	43, S2-S3, S12-S13,	Cilia Forceps, Barraquer.....	101
Anterior Synechia Scissors, Castroviejo	83S16, S18-S19, S21, S23-S24		Cilia Forceps, Beer.....	101
Aspiration Adapter.....	264	Cannula	211	Cilia Forceps, Gradle.....	101
Aspiration Cannula.....	193-196, 198, 210,	Cannula, Heidelberg Model.....	211	Cilia Forceps, Henry.....	101
.....S12, S21		Capsule Forceps	S4	Clamp For Infusion Tube	27
Aspiration Cannula, Brauweiler.....	196, 198	Capsule Forceps, H.-R. Koch.....	S4	Cleaning Adapter.....	97, 264-266
Aspiration Cannula Handpiece	198	Capsule Polishing Cannula	195-196, 226,	Clear-Cornea-Incision Fixation Ring.....	S3, S8
Aspiration Cannula Handpiece, Brauweiler	198S2, S6		Clear Cornea Knife.....	217
Aspiration Cannula, H.-R. Koch.....	196	Capsule Polishing Cannula, Greite	196	Colibri Forceps.....	104, S2-S4, S9, S13,
Aspiration Cannula, Roider.....	210	Capsule Polishing Cannula, Kratz.....	195, S2S16, S19, S21, S23-S24	
Aspiration Cannula, Weber	210, S21	Capsule Polishing Cannula, Kyu.....	195	Combined Needle Holder / Scissors.....	88, 143
Aspiration Handpiece	198-201, 227, S3, S5	Capsule Polishing Cannula, Papst.....	195	Conjunctiva Forceps.....	102, S16-S18,
Aspiration Handpiece, Brauweiler.....	200	Capsule Polishing Curette	160, S6S21, S23-S24	
Aspiration Handpiece, D.D. Koch - T. Kohnen	198	Capsule Polishing Curette, Rentsch	160, S6	Conjunctiva Forceps, Draeger....	102, S16-S18,
Aspiration Handpiece, H.-R. Koch	200	Capsule Scissors	85-86, S4-S5S21, S23-S24	
Aspiration Handpiece, Koch-Kohnen.....	S3, S5	Capsule Scissors, Brauweiler.....	85	Conjunctiva Forceps, Foerster.....	102
Aspiration Handpiece, R. Koch	201	Capsule Scissors, H.-R. Koch.....	85, S4	Conjunctiva Forceps, Freiburg Model	102
Aspiration Handpiece, Schmack.....	198-199	Capsule Scissors, Osher.....	86	Conjunctiva Forceps, Kaufmann	102
Aspiration Polishing Cannula	195	Capsule Scissors, Yae.....	85, S5	Conjunctiva Forceps, Krefeld Model	102
Aspiration Polishing Cannula, Teping	195	Capsulorhexis Forceps	113-117, S2-S5	Corneal Irrigation Cannula	197, S10
Aspiration Speculum.....	S7-S8	Capsulorhexis Forceps, Hattenbach	117, S4	Corneal Irrigation Cannula, Krumeich	197, S10
Aspiration Speculum, Liebermann-Tennant	S7-S8	Capsulorhexis Forceps, H.-R. Koch	114	Corneal Marker	41-42, 45, 47, S7, S9-S11
Aspiration Handpiece, Koch-Kohnen.....	S5	Capsulorhexis Forceps, Ikeda	1166, S5	Corneal Marker, Hoffer	41
B		Capsulotomy Scissors	83-85, S2-S4,	Corneal Marker, Lindstrom	42
Backflush Cannula.....	207, 227-228,S9-S10, S12-S13, S15, S21, S23-S24		Corneal Marker, Thomas Neumann	41, S9-S11
.....S21, S23-S24		Capsulotomy Scissors, Vannas.....	84-85,	Corneal Pendulum Marker	46, S6
Backflush Cannula, Kampik	207S2-S4, S9-S10, S13, S15, S21, S23-S24		Corneal Pendulum Marker, Scott	46, S6
Backflush Cannula, Laqua	207	Cerclage Forceps	127	Corneal Punch	70, S9
Backflush Handpiece.....	206, 227, S21, S23	Cerclage Forceps, Essen Model.....	127	Corneal Punch, Thomas Neumann	70
Backflush Tip	208, 209	Chalazion Curette	51-52, S15-S16	Corneal Scarifier	64, 57, S9, S16
Backflush Tip, Klein.....	208	Chalazion Curette, Axenfeld	52, S15	Corneal Scarifier, Kuhnt	64, S9, S14, S16
Ball Capsule Polisher	160	Chalazion Curette, Hebra	52	Corneal Scissors	81-82, S2, S16, S18-S19
Ball Capsule Polisher, Anis	160	Chalazion Curette, Meyhoefer	51, S15-S16	Corneal Scissors, Castroviejo	81
Battery Cautery	246	Chalazion Curette, Meyhoeffer	S15	Corneal Scissors, Katzin	81, S9-S10
Bone Chisel	243, S17	Chalazion Forceps	134, S15	Corneal Scissors, Sautter	81, 162, 192, 225,
Bone Chisel, West	243	Chalazion Forceps, Ayer	134S3-S5, S9, S12	
Bone Gouge	243	Chalazion Forceps, Cauer	134	Corneal Scissors, Troutman	81, S2
		Chalazion Forceps, Desmarres	134, S16	Corneal Section Scissors	82
		Chalazion Forceps, Lambert	134, S16	Corneoscleral Punch	77
		Chalazion Forceps, Lordan.....	134, S15-S16	Corneoscleral Punch, Gass	77
		Chalazion Knife	63	Crescent Knife	61
		Chopper	152	Crvo Knife	66

KEYWORD	PAGE	KEYWORD	PAGE	KEYWORD	PAGE
Cyclodialysis Spatula, Elschnig.....	164	Duct Trehine.....	69, S12	Eye Speculum, Mellinger-Beard.....	20
D		Duct Trehine, Schmidt-Mumm.....	69, S12	Eye Speculum, Murdock.....	24-25
Dalk Probe.....	186	Elevatorium.....	242, S17	Eye Speculum, Noda.....	27-29
Dalk Probe, Krumeich.....	186	Elevatorium, Heidelberg Model.....	242, S17	Eye Speculum, Ohno.....	17
DCR Lacrimal Intubation Set.....	224	Enclavation Spatula.....	160	Eye Speculum, Paulmann.....	27
Degree Gauge.....	46	Enclavation Spatula, Sekundo.....	160	Eye Speculum, Pretor.....	21
Degree Gauge, Mendez.....	46	Endothelial Insertion Spatula.....	168	Eye Speculum, Sauer.....	19
Descemet Incision Hook.....	168, S10-S11	Endothelial Insertion Spatula,		Eye Speculum, Schott.....	31-32
Descemet Membrane Forceps.....	117	Thomas Neuhann.....	168	Eye Speculum, Takeda.....	19
Descemet Membrane Forceps,		Enucleation Scissors.....	94, S14	Eye Speculum, Thomas Neuhann.....	19
Tobias Neuhann.....	117	Enucleation Spoon.....	167, S14	Eye Speculum, Weiss.....	21-22
Descemet Punch.....	S12	Enucleation Spoon, Wells.....	167, S14	Eye Sphere.....	245, S14
Descemet Punch, Kelly.....	S12	Enzyme Cannula.....	191	Eye Sphere, Cibis.....	245
Descemet Scraper.....	168, S10-S11	Ergoline Forceps.....	133	Eye Sphere, Harms.....	245
Descemet's Membrane Punch.....	76	Ergoline Forceps ILM.....	133	Eye Sphere, Walser.....	245, S14
Descemet's Membrane Punch, Kelly.....	76	Ergoline Forceps ILM, Naito.....	133	F	
Diaclean CCI Diamond Knife.....	S5	Ergoline Forceps, Tano.....	133	Fascia Needle.....	243
Diaclean Clear-Cornea-Incision.....		Ergoline Scissors.....	97	Fascia Needle, Wright.....	243
Diamond Knife, 1.5/2.0 mm.....	60	Ergoline Scissors, Heimann.....	97	Femtocataract Eye Speculum.....	26, S5
Diaclean Clear-Cornea-Incision.....		Ergoline Scissors, Kroll.....	97	Femtocataract Eye Speculum, Dick.....	26, S5
Diamond Knife, 2.2 mm.....	60	Ergoline Scissors, Naito.....	97	Femtocataract Eye Speculum, Solomon.....	26
Diaclean Clear-Cornea-Incision.....		Exenteration Curette.....	52, S14	Femtocataract Nucleus Rotator	
Diamond Knife, 2.85 mm.....	61	Exenteration Curette, Bunge.....	52, S14	and Chopper.....	169, S5
Diaclean Diamond Crescent Knife.....	61	Explantation Forceps.....	125, S4	Femtocataract Nucleus Rotator	
Diaclean Diamond Knife.....	59	Explantation Hook.....	157	and Chopper, Hoffmann.....	169, S5
Diaclean Tri-Facet Diamond Knife.....	60	Explantation Hook, Tandogan.....	157	Femtocataract Spatula.....	169, S5
Diaclean Tri-Facet Diamond Knife, Pham.....	60	Expressor Cannula.....	193	Femtocataract Spatula, Breyer.....	169, S5
Diamond Burr.....	244, S14	Expressor Cannula, Klöckner.....	193	Femtocataract Spatula, Dick.....	169, S5
Diamond Knife.....	59- 62, S12	Eye Douche.....	191, S3, S13	Femto Double Instrument.....	169, 170, S8
Diamond Micrometer Knife.....	62	Eye Douche, Cranz.....	191, S3, S13	Femto Double Instrument,	
DMEK Graft Retainer.....	71, S11	Eye Scissors.....	91-93, S2-S4, S9, S12,	Chansue.....	169, 170, S8
DMEK Graft Retainer, Szurman.....	71, S11S14-S19, S21, S23-S24		Femto Double Instrument, Guell.....	170
DMEK Irrigation Handpiece.....	203	Eye Scissors, Bonn Model.....	92, S2, -S4, S9, S12	Femto Flap Lifter.....	171
DMEK Irrigation Handpiece, Szurman.....	203	Eye Scissors, Graefe.....	93	Femto Flap Lifter, Pfaeffl.....	171
DMEK Keratoplasty Forceps.....	105, S11	Eyeshield.....	245	Femto-LASIK Double-Ended	
DMEK Keratoplasty Forceps, Remky... <td>105, S11</td> <td>Eye Speculum.....</td> <td>15-32, 223, S2-S5, S7-S10,</td> <td>Instrument.....</td> <td>170, 172</td>	105, S11	Eye Speculum.....	15-32, 223, S2-S5, S7-S10,	Instrument.....	170, 172
DMEK Preparation Forceps.....	117, S11S1-S15, S18-S19, S21, S23-S24		FemtolASIK Double-Ended Instrument,	
DMEK Preparation Forceps, Szurman	117, S11	Eye Speculum, Bangerter.....	23	Breyer.....	170
DMEK Preparation Scraper.....	64, S11	Eye Speculum, Barraquer.....	15-17, 19,	FemtolASIK Double-Ended Instrument,	
DMEK Preparation Scraper, Szurman...	64, S11S2-S4, S12-S13, S15		Kohnen.....	172
DMEK Preparation Spatula.....	168, S11	Eye Speculum, Barraquer Frankfurt Model.....	16	Femto-LASIK Flap Forceps.....	118, S8
DMEK Preparation Spatula, Szurman	168, S11	Eye Speculum, Barraquer-Oosterhuis.....	18	Femto-LASIK Flap Forceps, Knorz.....	118, S8
DMEK Transplant Forceps.....	117	Eye Speculum, Bartl.....	22	Femto-LASIK Hook.....	S8
DMEK Transplant Forceps, Yörük.....	117	Eye Speculum, Castroviejo.....	23	Femto-LASIK Hook, Breyer.....	S8
Dressing Forceps.....	110, S2, S9, S12-S15,	Eye Speculum, Cook.....	23-24, S13-S14, S18	Femto-LASIK Hook Instrument.....	171
.....S17-S19, S21, S23-S24		Eye Speculum, Kammann.....	29	Femto-LASIK Round Spoon.....	170, S7
DSAEK/DMEK Rhexis Forceps.....	117, S10-S11	Eye Speculum, Kershner.....	25, S9, S14, S19,	Femto-LASIK Round Spoon, Blum.....	S7
DSAEK/DMEK Rhexis Forceps,	S21, S23-S24		FemtoLASIK Round Spoon, Blum.....	170
Gerten.....	117, S10-S11	Eye Speculum, Koga / Noda.....	28	Fixation And Holding Forceps.....	122
DSAEK Forceps.....	118, S10	Eye Speculum, Kratz-Barraquer.....	17, 19	Fixation And Holding Forceps, Dick.....	122
DSAEK Forceps, Thomas Neuhann....	118, S10	Eye Speculum, Liebermann-Tennant.....	25, 28,	Fixation Forceps.....	109, S14, S16
DSAEK Implantation Forceps.....	118, S10S2-S4, S7-S8, S10, S12-S13		Fixation Forceps, Barraquer.....	109, S14
DSAEK Insertion Spatula.....	S10	Eye Speculum, Mellinger.....	20	Fixation Forceps, Graefe.....	109, S14, S16
DSAEK Insertion Spatula, Thomas Neuhann.....S10					

KEYWORD	PAGE	KEYWORD	PAGE	KEYWORD	PAGE
Fixation Forceps, Heidelberg Model.....	109	Holder, Thomas Neuhann.....	186	IOL Forceps.....	125, S4, S6
Fixation, Holding and Measuring Instrument.....	150, S23-S24	Holding Forceps.....	122, 124, 125, 127, S3	IOL Forceps, Snyder-Osher.....	125, S4, S6
Fixation Instrument.....	S21	Holding Forceps, Pfennigwerth.....	127	IOL Haptic Expander.....	122
Fixation Ring.....	149-150, S3, S8-S9	Holding Forceps, Strobel.....	125	IOL Haptic Expander, Dick.....	122
Fixation Ring, Fine-Thornton.....	149-150	Horizontal Axis Marker.....	46, S6	IOL Scissors.....	88
Fixation RIng, Fine-Thornton.....	S7-S8	Horizontal Axis Marker, Gerten.....	46, S6	IOL Scissors.....	87-88, S4, S6
Fixation Ring, Thornton.....	149	Horseshoe-Hook.....	175, S21, S23-24	IOL Scissors, Snyder-Osher.....	88, S4, S6
Fixed Caliper.....	43	Horseshoe-Hook,		IOL Scissors, Steinmetz.....	87, S4
Fixed Caliper, Braunstein.....	43	Hattenbach.....	175, S21, S23-S24	Iris Forceps.....	110-112, S2, S12, S16, S19
Flap Forceps.....	S8	Hybrid Combo Prechopper.....	120	Iris Forceps, Bonn Model.....	110-111, S2, S12
Flap Lifter.....	S7-S8	Hybrid Combo Prechopper, Akahoshi.....	120	Iris Forceps, Elze.....	112
Flap Lifter, Pfaeffl.....	S7, S8	Hybrid Needle Holder.....	143, S4, S6	Iris Forceps, Hamburg Model.....	112
Flaporhexis Forceps.....	118, S7, S8	Hybrid Needle Holder,		Iris Forceps, Heidelberg Model.....	111
Flaporhexis Forceps, Pfaeffl.....	118, S7, S8	Hattenbach.....	143, S4, S6	Iris Forceps, Swiss Model.....	111-112
Folder.....	123, 124, S4	Hybrid Scissors.....	86, S4, S6	Iris Forceps, Thomas Neuhann.....	112
Folder, H.-R. Koch.....	124, S4	Hybrid Scissors, Hattenbach.....	86, S4, S6	Iris Hook.....	154, 157-158, S2-S3
Folder, Spaleck.....	123	Hybrid Tying Forceps.....	108, S4, S6	Iris Hook and Lens Manipulator.....	157-158
Forceps.....	128	Hybrid Tying Forceps, Hattenbach..	108, S4, S6	Iris Hook and Lens Manipulator, Dardenne.....	157
Foreign Body Gouge.....	56	Hydrodissection Cannula.....	192, 225, 217, S2-S5, S9, S12	Iris Hook and Lens Manipulator, Frankfurt Model.....	158
Foreign Body Instrument.....	55, S13, S15	Hydrodissection Cannula,		Iris Hook and Nucleus Dividing Hook.....	154
Foreign Body Needle.....	56	Frankfurt Model.....	192, S2, S12-S13	Iris Hook and Nucleus Dividing Hook, Thomas Neuhann.....	154, S3
Foreign Body Spud.....	56	Hydrodissection Cannula, Sautter.....	192, S2-S5, S9, S12	Iris Hook and Nucleus Dividing Hook, T. Kohnen.....	154
Foreign Body Spud, Plange.....	56	Hydrodissection Cannula, Spirig.....	192	Iris Hook, Bonn Model.....	158
Fragment Forceps.....	108, S4, S8	Hydrodissection Cannula,		Iris Hook, Dardenne.....	S2
Fragment Forceps, Gaskin.....	108, S4	Tuebingen Model.....	192	Iris Hook, Lund.....	158
Fragment Forceps, Kelman-McPherson	108, S8	I		Iris Scissors.....	86-87
G		I/A Cannula, Simcoe.....	S2, S4	Iris Scissors, Barraquer.....	87
Glaucoma Punch.....	77, S12	I/A Handpiece.....	203-205, 227	Iris Scissors, Noyes.....	86
Glaucoma Punch, Grehn.....	S12	I/A Handpiece, Hunold.....	205	Iris Scissors, Wecker.....	86-87
Glaucoma Punch, Mohr.....	77	I/A Handpiece, Papst.....	205	Iris Spatula.....	165-166
Globe Fixation Ring.....	149	ILM-Scraper.....	175, S21	Iris Spatula, Bangerter.....	166
Globe Fixation Ring, Hattenbach.....	149	ILM-Scraper, Bartz-Schmidt.....	175, S21	Iris Spatula, Frankfurt Model.....	165
Grooved Silicone Sponge Implant (506G Style).....	237	Implantation Forceps	118, 121-123, S2, S4	Iris Spatula, Wecker	165
Grooved Silicone Sponge Implant (507G Style).....	237	Implantation Forceps, Doden.....	121, S2, S4	Irrigating Chopper, Dick.....	202
Grooved Silicone Strip (31 Style).....	234	Implantation Forceps, Fechner.....	121	Irrigation / Aspiration Cannula.....	194-195
Grooved Silicone Strip (32 Style).....	234	Implantation Forceps, Sekundo.....	123	Irrigation / Aspiration Cannula Coaxial.....	195
Grooved Silicone Strip (219 Style).....	234	Implantation Forceps, Shepard.....	122	Irrigation / Aspiration Cannula, Simcoe	194
H		Implantation Forceps, T. Kohnen.....	122	Irrigation Cannula.....	194, 196-198, 226, S7, S10, S12
Hammer.....	242, S17	Implant Knife.....	S2	Irrigation Cannula, Binkhorst.....	194
Hammer, Mead.....	242, S17	Incision Measuring Caliper.....	43	Irrigation Cannula, Bonn Model.....	194, S15
Handle.....	85, 89-94, 116, 218, S5	Incision Measuring Caliper, Osher.....	43	Irrigation Cannula, Brauweiler.....	196, 198
Handle, Japan Model.....	85, 116	Indicator.....	255	Irrigation Cannula Handpiece.....	197-198
Handle UNO Colorline.....	S21, S23	Infusion Cannula.....	191, 209, 228, S10, S21	Irrigation Cannula Handpiece, Brauweiler..	198
Handle UNO Colorline 20G, 23G, 25G.....	218	Infusion Sleeve.....	203-205	Irrigation Cannula Handpiece, Keerl.....	198
Handpiece For Capsule Scissors.....	86	Injection Cannula.....	191, 193, 209-210, 226, 229	Irrigation Cannula, Mueller.....	194
Healon - Viscoelastic Cannula.....	193	Injection Cannula, Duisburg Model....	209, S21	Irrigation Cannula, Rentsch.....	194
Holder.....	186	Injection Cannula, Hamburg Model.....	210	Irrigation Chopper.....	202
Holder for Absorbent Strips.....	S7, S8	Injector.....	241, S4	Irrigation Chopper, R. Koch.....	202
Holder For Absorbent Strips.....	145	Intraocular Magnet.....	246, S21	Irrigation Cystotome Cannula.....	191
		Intraocular Magnet, Kuhn.....	246, S21	Irrigation Descemet Hook.....	202, S11
				Irrigation Descemet Hook,	

KEYWORD	PAGE	KEYWORD	PAGE	KEYWORD	PAGE
Althaus/Cartsburg.....	202, S10-S11	Lacrimal Probe, Wilder.....	183, S15, S17	Lens Loop, Snellen.....	161
Irrigation Descemet Scraper.....	202, S10-S11	Lacrimal Sac Retractor.....	37-38, S16-S17	Lens Loop, Spirig.....	161
Irrigation Descemet Scraper,		Lacrimal Sac Retractor, American Model.....	38	Lens Loop, Weber.....	161
Neuhann.....	202, S10-S11	Lacrimal Sac Retractor, Axenfeld.....	37, S16-S17	Lens, Machemer.....	244
Irrigation Handpiece.....	198-201, 203, 227,S3, S5	Lacrimal Sac Retractor, Goldstein.....	38	Lens Manipulation Hook.....	156, S2-S4
Irrigation Handpiece, Brauweiler.....	201	Lacrimal Sac Retractor, Mueller.....	38	Lens Manipulation Hook, Sinskey..	156, S2-S4
Irrigation Handpiece,		Lacrimal Sac Retractor, Rollet.....	37	Lens Manipulator.....	155-158
D.D. Koch - T. Kohnen.....	199	Lasek Epi Peeler.....	64	Lens Manipulator and Hook.....	156
Irrigation Handpiece, Koch-Kohnen.....	S3, S5	Lasek Epi Peeler, Lohmann.....	64	Lens Manipulator and Hook, Rentsch.....	156
Irrigation Handpiece, Schmack.....	198	LASEK Epi Peeler, Lohmann.....	S7	Lens Manipulator, Bremen Model.....	156
Irrigation Trabeculotomy Cannula.....	193	LASEK Funnel.....	243, S7	Lens Manipulator, Jacobi.....	156
Irrigation Trabeculotomy Cannula, Burk.....	193	LASEK Funnel, Lohmann.....	243, S7	Lens Manipulator, Mainz Model.....	155
Irrigation Vectis Cannula.....	193, 226, S2	LASEK Spatula.....	168	Lens, Peyman.....	244
J		LASEK Trehpine.....	243, S7	Lens Spatula.....	157
Jewelers Type Forceps.....	109	LASIK Cannula.....	196, S7	Lens Spatula, Hirschmann.....	157
K		LASIK Cannula, Bauer.....	S7	Lens, Tolentino.....	244
Karate Phaco Chopper.....	152	LASIK Cannula, W. Bauer.....	196	Lenticule Forceps.....	119, S8
Karate Phaco Chopper, Agarwal.....	152	LASIK Corneal Flap Forceps.....	119	Lenticule Forceps, Rupal Shah.....	119
Keratome.....	63	LASIK Corneal Flap Forceps, Buratto.....	119	Lenticule Micro Forceps.....	S8
Keratome, Hamburg Model.....	63	Lasik Corneal Marker.....	45, S7	Lenticule Micro Forceps, Rupal Shah.....	S8
Keratometer.....	45, S9	Lasik Corneal Marker.....	S7	Lid Plate.....	36, S15-S16
Keratometer, Maloney.....	45, S9	Lasik Corneal Marker, Chayet.....	45	Lid Plate, Jaeger.....	36, S15-S16
Keratoplasty Forceps.....	104-105, S9	LASIK Corneal Marker, Chayet.....	S7	Lid Retractor.....	35-36, S13-S16, S18
Keratoplasty Forceps, Remky.....	105	Lasik Corneal Marker, Lu.....	45	Lid Retractor, Desmarres.....	35-36,S13-S16, S18
Keratoplasty Forceps,		LASIK Corneal Marker, Lu.....	S7	Lid Retractor, Giessen Model.....	36
Thomas Neuhann.....	104, S9	Lasik Corneal Marker, Mendez.....	45	Limbal Relaxing Incision Diamond Knife.....	61
Keratoplasty Knife.....	64, S9	LASIK Corneal Marker, Mendez.....	S7	Localizer and Scleral Depressor.....	180
Keratoplasty Knife, Pauifique.....	64, S9	LASIK Corneal Protector.....	167, S7-S8	Localizer and Scleral Depressor, Meyer-Schwickerath.....	180
Keratoplasty Scissors.....	81-82, S9-S10	LASIK Corneal Protector, Zirm.....	167, S7-S8	M	
Keratoplasty Scissors, Katzin.....	S9-S10	LASIK Dissector Cannula.....	197, S7-S8	Manipulator.....	154-158
Keratoplasty Scissors, Mackensen.....	82	LASIK Dissector Cannula, Zirm.....	197, S7-S8	Manipulator, Frankfurt Model.....	156
Keratoplasty Spatula.....	163, S9	Lasik Eye Speculum.....	26, 30	Manipulator, Mainz Model.....	155
Keratoplasty Trehpine.....	69-70	LASIK Eye Speculum.....	S7	Marker.....	S7, S12, S18
Keratoplasty Trehpine, Best.....	69	Lasik Eye Speculum, Knorz.....	26	Marker, Fink.....	S18
Keratoplasty Trehpine, Franceschetti.....	69	LASIK Eye Speculum, Knorz.....	S7	Marker, Fronimopoulos.....	S12
Keratoplasty Trehpine, Hamburg Model.....	69	Lasik Eye Speculum, Maus.....	30	Marking Instrument.....	S21, S23-S24
L		Lasik Eye Speculum, Zirm.....	30	Marking Pen.....	231, S6-S7, S16, 231
Lacrimal Cannula.189- 190, 224, S3, S13, S17		LASIK Flap Forceps.....	S7, S8	Measuring Spatula.....	44
Lacrimal Cannula, Anel.....	189, S17	LASIK Flap Forceps, Buratto.....	S7, S8	Measuring Spatula, Engels.....	44
Lacrimal Cannula, Bonn Model.....	190	LASIK Irrigation Cannula.....	197, S7	Melag Metal Grid Insert.....	263
Lacrimal Intubation Set.....	223-224, S17	LASIK Irrigation Cannula, Zirm.....	197, S7	Membrane Peeler.....	175, S21, S23-S24
Lacrimal Intubation Tubing.....	223	LASIK Spatula.....	167, S8	Membrane Peeler, Morris-Kroll.....	175
Lacrimal Probe...183, 185-186, 190, S15, S17		LASIK Spatula, Neumann.....	S7	Membrane Spatula, Lucke.....	175, 179
Lacrimal Probe, Bowman.....	185, S17	LASIK Spatula, Thomas Neuhann.....	, 167	Meridional Silicone Implant (106 Style).....	235
Lacrimal Probe Cannula.....	190, S17	Lens.....	88, 244, 245, S2, S4, S21, S23, S24	Metal Stand.....	72
Lacrimal Probe Cannula, Bangerter.....	190	Lens Extraction Hook.....	155	Micro Colibri Forceps.....	104, S2-S3, S9,S16, S21, S23-S24
Lacrimal Probe Cannula, Juenemann. 190, S17		Lens Extraction Hook, Hennig / Friedrich....	155	Microfinger.....	154
Lacrimal Probe, Giessen Model.....	183, S17	Lens Handle Ring.....	245	Microfinger, Liebermann.....	154
Lacrimal Probe, Juenemann.....	186, 190, S17	Lens, Kilp.....	244, S21, S23-S24	Micro Forceps.....	116, S8
		Lens, Landers.....	244	Micro Forceps, Tobias Neuhann.....	116
		Lens Loop.....	157, 161, S2	Micro Round Knife.....	63, S12
		Lens Loop, Arlt.....	161, S2, S4		
		Lens Loop Manipulator.....	157		
		Lens Loop Manipulator, Lester.....	157		
		Lens Loop Manipulator, Zirm.....	157		

KEYWORD	PAGE	KEYWORD	PAGE	KEYWORD	PAGE
Micro Stitch Scissors.....	82	NanOEdge Trephines.....	215	Oval Silicone Sleeve (3084 Style).....	235
MICROSTO Container high.....	S21, S23-S24	NanOEdge Vacuum Corneal Punch.....	59	Oval Silicone Sponge Implant (501 Style)....	236
MICROSTOP Clip for Retainer Elements.....	250	NanOEdge Vacuum Trephine.....	216, S9	Oval Silicone Sponge Implant (506 Style)....	236
MICROSTOP Container high.....	249	NanOEdge Vacuum Trephine Punches.....	216	Oval Silicone Sponge Implant (507 Style)....	237
MICROSTOP Container High.....	S2-S3	NanOEdge Vacuum Trephine radial.....	216		
MICROSTOP Container low.....	249	Nasal Forceps.....	127		
MICROSTOP ID Label.....	250	Nasal Forceps, Troeltsch.....	127		
MICROSTOP Label.....	250	Nasal Speculum.....	242, S17	Paracentesis Knives.....	217
MICROSTOP Logistic Frame.....	251	Needle Holder.....	88, 139-144, S2-S4, S9-S10, S12-S19, S21, S23-S24	Pendulum Instrument.....	47
MICROSTOP Retainer Element.....	250, S2-S3, S21, S23-S24	Needle Holder, Arruga.....	144	PEP 2.2 mm I/A Handpiece.....	203
MICROSTOP Seal.....	251	Needle Holder, Barraquer.....	139-141, S2-S4, S9-S10, S12-S13, S19	PEP 2.2 mm Infusion Sleeve.....	203
MICROSTOP Tag.....	251	Needle Holder, Barraquer-Troutman.....	140	PEP 2.4 mm I/A Handpiece.....	204
MICROSTOP Tray Lid mini.....	250	Needle Holder, Boynton.....	144	PEP 2.4 mm Infusion Sleeve.....	204
MICROSTOP Tray mini.....	249, S2-S3, S21, S23-S24	Needle Holder, Castroviejo.....	142, S2, S12-S19, S21, S23-S24	Permanent Textile Filter.....	255
Micro Trabeculectomy Punch.....	76	Needle Holder, Heidelberg Model.....	144	PFCL Cannula.....	210, 228, S23
Micro Trabeculectomy Punch, Grehn.....	76	Needle Holder, Jacobi.....	139, S9	PFCL Cannula, Chang.....	210
MICS Capsule Forceps.....	115-116	Needle Holder, Kalt.....	143	Phaco Spatula.....	159
MICS Capsule Forceps, H.-R. Koch.....	115-116	Nucleus Dividing Hook.....	151-154, S3, S5	Phaco Spatula, Amon.....	159
MICS Capsulorhexis Forceps.....	115-116, S5	Nucleus Dividing Hook, Haefliger.....	152	Phaco Spatula, D. D. Koch.....	159
MICS Capsulorhexis Forceps, Auffarth.....	S5	Nucleus Dividing Hook, Miyoshi.....	151	Phaco Spoon.....	159
MICS Capsulorhexis Forceps, Ikeda.....	116	Nucleus Dividing Hook, Nagahara	152-153, S3, S5	Phaco Spoon, Engels.....	159
MICS Capsulorhexis Forceps, Inamura.....	115	Nucleus Dividing Hook, Thomas Neuhann.....	151, S5	Pigtail Probe.....	184, S17
MICS Capsulorhexis Forceps, Kampik.....	115	Nucleus Dividing Hook, Völcker.....	151	Prechopper.....	120, S3
MICS Double Cross Action Capsulorhexis Forceps.....	115, S3	Nucleus Dividing Hook, Zirm.....	151	Prechopper, Akahoshi.....	120, S3
MICS Double Cross Action		Nucleus Manipulator and Anterior Capsule Polisher.....	155	Prefolder.....	124
Capsulorhexis Forceps, Auffarth.....	115, S3	Nucleus Manipulator and Anterior Capsule Polisher, Drysdale.....	155	Prefolder, H.-R. Koch.....	S4
Mosquito Forceps.....	135, S2, S9, S12-S14, S16-S19, S21, S23-S24	Nucleus Rotator.....	153, 155, 169	Pressure Tube.....	211
Mosquito Forceps, Halstead.....	135, S2, S9, S12-S14, S17-S19, S21, S23-S24	Nucleus Rotator and Chopper.....	153	Probe for Stromal Channel.....	184
Mosquito Forceps, Hartmann.....	135, S16, S21, S23-S24	Nucleus Rotator and Dividing Hook.....	153	Protocol Cards.....	251
Muscle Forceps.....	126-127, S14, S18-S19	Nucleus Rotator and Dividing Hook, Thomas Neuhann.....	153	Ptosis Forceps.....	134, S16
Muscle Forceps, Bangerter.....	126, S14, S18	Nucleus Rotator, Thomas Neuhann....	153, 155	Ptosis Forceps, Berke.....	134, S16
Muscle Forceps, Sattler.....	127, S18-S19	Nucleus Rotator, Thomas Neuhann and Chopper, Tobias Neuhann.....	153	Punch.....	S9, S12
N		Nucleus Spreading Forceps.....	120		
NanOEdge Crescent Knife.....	217	Nucleus Spreading Forceps, Akahoshi.....	120		
NanOEdge Implant Knife.....	217, S2	Nucleus Sustainer.....	154		
NanOEdge MVR Stilett.....	217	Nucleus Sustainer, Akahoshi.....	154		
NanOEdge MVR Stiletto.....	S21				
NanOEdge Paracentesis Knife.....	S2-S3, S12-S13, S16, S19				
NanOEdge Phaco Slit Knife.....	217, S3				
NanOEdge Single-Use Clear Cornea Knife.....	217				
NanOEdge Single-Use Corneal Cutting Block.....	215				
NanOEdge Single-Use Knives.....	217				
NanOEdge Single-Use Paracentesis Knives	217				
NanOEdge Spoon Knife.....	217				
NanOEdge Trephine Punches.....	215				
O					
Oculopressor.....	48				
Oculopressor, Homburg Model.....	48				
Operating Scissors.....	94				
Orbital Spatula...162-163, S19, S21, S23-S24					
Orbital Spatula, Behrens-Baumann.....	162				
Orbital Spatula, Hamburg Model.....	162, S19				
Orbital Spatula, Sautter.....	162				
Oval Silicone Sleeve (3083 Style).....	235				

KEYWORD	PAGE	KEYWORD	PAGE	KEYWORD	PAGE
Round Silicone Sponge Implant (505 Style).....	236	Silicone Chamber	206, 227	Sautter	225, S5
Ruler.....	44, S19	Silicone Circling Band (40 Style).....	233	Single-Use Hydrodissection Cannula, Tuebingen Model	225
Ruler, Giessen Model	44	Silicone Circling Band (240 Style).....	233	Single-Use I/A Cannula.....	226
S		Silicone Closure	258	Single-Use I/A Cannula, Simcoe	226
Safety Seal.....	258	Silicone Insert	72	Single-Use Infusion Cannula.....	228
Sato Knife.....	65	Silicone Lens Bi-Sector	88, S4	Single-Use Infusion Tube 20G	228
Sato Knife, Dardenne	65	Silicone Lens Bi-Sector, H.-R. Koch	88, S4	UNO Colorline	222
Sato Knife, Foerster	65	Silicone Nap Strip blue	257	Single-Use Injection Cannula	226, 229
Sato Knife, Hansen	65	Silicone Perforated Tray	256	Single-Use Injection Cannula, Kirchhof	229
Scalpel Handle	231, S14-S15	Silicone Plate	257	Single-Use Iris Retractor	223, S2, S4
Scalpel Surgical Blades	231, S14-S15	Silicone Replacement Ring	244-245	Single-Use Irrigation Cannula	226
Scissors	218, 143, S2-S6, S9-S10, S12-S14, S16-S19, S21, S23-S24	Silicone Sleeve (70 Style)	234	Single-Use Irrigation Vectis Cannula	226
Scissors, Barraquer	S12	Silicone Sleeve (72 Style)	234	Single-Use Knives	217
Scissors, Westcott	89, S2, S9, S12-S13, S16, S18-S19, S21, S23-S24	Silicone Sleeve (270 Style)	234	Single-Use Lacrimal Cannula	224, S17
Scleral Depressor	179-180, S15, S19	Silicone Strip (41 Style)	233	Single-Use Mini Filter	229
Scleral Depressor, Kreissig	179	Silicone Strip (42 Style)	233	Single-Use Oil Injection System	230
Scleral Depressor, Laqua	179, S19	Silicone Strip blue	257	Single-Use Peribulbar Cannula	224
Scleral Depressor, Lucke	179	Silicone Tire (276 Style)	232	Single-Use PFCL Cannula	228, S23
Scleral Depressor, Schocket	180	Silicone Tire (277 Style)	232	Single-Use Pressure Tube	229
Scleral Fixation Ring	150	Silicone Tire (279 Style)	232	Single-Use Retrobulbar Cannula	224, S3
Scleral Fixation Ring, Flieringa	150	Silicone Tray	256	Single-Use Retrobulbar Cannula, Atkinson	224, S3
Scleral Flap Marker	44, S12	Silicone Wedge (137 Style)	235	Single-Use A. C. Cannula	225
Scleral Flap Marker, Tobias Neuhann	44, S12	Single-Use A. C. Cannula, Rycroft	225	Single-Use Silicone Chamber	227
Scleral Knife	63	Single-Use A. C. Maintainer	225	Single-Use Syringe	229
Scleral Knife, Bonn Model	63	Single-Use A. C. Maintainer, Lewicky	225	Single-Use Tano Forceps	S24
Scleral Plug	128, S21	Single-Use Backflush Cannula	S24	Single-Use Tano Forceps UNO Colorline	219
Scleral Plug Forceps	128	Single-Use Backflush Cannula	227-228, S21, S23-S24	Single-Use Three-Way Stop Cock	231
Scleral Ruler	44, S18	Single-Use Backflush Handpiece	227	Single-Use Trocar Plugs UNO Colorline	221
Scleral Ruler, Vogt	44, S18	Single-Use Bimanual I/A Handpiece Set	227	Single-Use Trocar System One-Step UNO Colorline	221
Seal	258	Single Use Cannula	228, S23	Single-Use Trocar System One-Step Single-Use Vitreous Forceps	221
Seal / Label	261	Single-Use Capsule Polishing Cannula	226	Crocodile Design	219, S21, S23
Security Seal	262	Single-Use Capsule Polishing Cannula, Kratz	226	Single-Use Vitreous Forceps	219, S23-S24
Semi-Circular Silicone Sponge Implant (510 Style)	237	Single-Use Cystotome Cannula	224, S2	Endgripping	219, S23-S24
Semi-Oval Silicone Sponge Implant (511 Style)	238	Single-Use DMEK Cartridge	S11	Single-Use Vitreous Forceps	219
Semi-Oval Silicone Sponge Implant (515 Style)	238	Single-Use DMEK-Cartridge	223	Wave Design UNO Colorline	219
Semi-Oval Silicone Sponge Implant (517 Style)	238	Single-Use DMEK Cartridge, Szurman	S11	Single-Use Vitreous Forceps	219
Serrefine	136, S2, S18-S19, S21, S23-S24	Single-Use DMEK-Cartridge, Szurman	223	Wave Design UNO Colorline, Hattenbach	219
Serrefine, Dieffenbach	136, S2, S18-S19, S21, S23-S24	Single-Use Endo-Probe C UNO Colorline	220	Single-Use Vitreous Scissors	218, S24
Silicone Ball	48	Single-Use Endo-Probe S UNO Colorline	220	UNO Colorline	218, S24
Silicone Block	71	Single-Use Eye Speculum	223	Spare Silicone Tube	210
Silicone Block for DMEK	71	Single-Use Fiber Optic S UNO Colorline	221	Spare Tip	246
Silicone Block for DMEK, Örk Linz Model	71	Single-Use Fiber Optic WS UNO Colorline	222	Spare Tubing Set	26
Silicone Block With Center Point	71	Single-Use Fiber Optic W UNO Colorline	221	Spatula	44, 153, 158-159, 162-163, 166-167, 192, S2-S3, S5, S7, S10, S12, S19, S21, S23-S24
Silicone Bulb	212	Single-Use Forceps, Tano	S24	Spatula and Dividing Hook	153
Silicone Bulb	S15	Single-Use High Speed Vitreous Cutter	MACH 2 UNO Colorline	Thomas Neuhann	153
Silicone Burl Mat	257, 262	Single-Use Hydrodissection Cannula	222	Spatula and Nucleus Dividing Hook	153
		UNO Colorline	222	Spatula and Retractor	167
		Single-Use Hydrodissection Cannula	225, S5	Spatula and Retractor, Giessen Model	167
		Single-Use Hydrodissection Cannula, Frankfurt Model	225	Spatula, Bangerter	162, S19
		Single-Use Hydrodissection Cannula, Essen Model		Spatula, Essen Model	166

KEYWORD	PAGE	KEYWORD	PAGE	KEYWORD	PAGE
Spatula, Rentsch	159, S10	Superior Rectus Forceps, Remky	125, S2, S12, S19	Triangle Scleral Depressor, Hattenbach	180
Spatula, Spirig	158	Suture Hook	172, S17	Tri-Facet Diamond Step Knife	61, S12
Spatula, Steinmetz	158	Suture Hook, Muenster Model	172, S17	Trocар Closure Forceps	129
Spatula, Thomas Neuhann and Nucleus Dividing Hook, Nagahara	153	Suture Marker	42, S9	Trocар Closure Forceps, Bartz-Schmidt	129
Speculum	S2-S3, S7-S8, S14, S23-S24	Suture Marker, Duncker	42, S9	Trocар Forceps	129, S21, S23
Spreading Forceps	120, 127, S19	Suture Marker, Thomas Neuhann	42	Tying Forceps	105-108, S2-S4, S9-S10, S12-S13, S16-S19, S21, S23-S24
Spreading Forceps, Watzke	127, S19	Suture Removing Forceps	104, S15	Tying Forceps, Bangerter	107
Sterilization Case	252-255, 257-259	Suture Removing Forceps, Barraquer	104	Tying Forceps, Bonn Model	106
Sterilization Rack	256	Suture Tool	172	Tying Forceps, Dardenne	107
Sterisafe Duro A4	259	Suturing Forceps	103-105, S3, S9, S13-S14, S16-S19, S21, S23-S24	Tying Forceps, Paufique	105
Sterisafe Duro A4 Plus	259	Suturing Forceps, Castroviejo	103, S3, S9, S13, S16, S18-S19, S21, S23-S24	Tying Forceps, Tuebingen	S9
Sterisafe Duro A4 Plus Complete Set	259, S4-S5, S7, S12, S14-S15	Suturing Forceps, Munich Model	103	Tying Forceps, Tuebingen Model	106, S2-S3, S10, S12-S13, S16-S19, S21, S23-S24
Sterisafe Duro A6	259	Suturing Forceps, Zuerich Model	104		
Sterisafe Duro Metal Tray A4-15	260	Synechiai Spatula	163	U	
Sterisafe Duro Metal Tray A4-30	260	Synechiai Spatula, Binkhorst	163	Universal Prechopper	121
Sterisafe Duro Metal Tray A6-15	260	Synechiai Spatula, Castroviejo	163	Universal Prechopper, Akahoshi	121
Sterisafe Duro Wire Basket A4-30	260	Syringe Clip	212	Universal Prefolder, H.-R. Koch	124
Sterisafe Duro Wire Basket A6-30	260	Syringe Holder	211	Universal Suturing Forcesps	105, S14, S17, S21, S23-S24
Sterisafe Filter	261			Universal Suturing Forcesps, Paufique	105, S14, S17, S21, S23-S24
Sterisafe Filter Tag Expiration Date	261	T		UNO Colorline Forceps	219, S21, S23-S24
Sterisafe Metal Grid Insert	262-263	Temporal Approach Eye Speculum	17	UNO Colorline Forceps, Tano	219, S21
Sterisafe Silicone Burl Mat	262	Temporal Approach Eye Speculum, H.-R. Koch	17	UNO Colorline Single-Use Instruments	218-222
Sterisafe Silicone Grid Insert	262	Tenotomy Scissors	89, 90, S3, S13-S14, S16-S19, S21, S23-S24	UNO Colorline Trocar System	223, S21, S23
Sterisafe Silicone Retainer Strip 75	261	Tenotomy Scissors, Bonn Model	90	UNO Colorline Trocar System, Hattenbach / Nikolic	S21, S23
Sterisafe Silicone Retainer Strip 145	261	Tenotomy Scissors, Stevens	89, S3, S13-S14, S16-S19, S21, S23-S24	UNO Colorline Vitreous Forceps	219, S24
Sterisafe Silicone Retainer Strip 230	261	Tissue Forceps	110	UNO Colorline Vitreous Forceps, Crocodile Design	219, S21, S23
Stiletto	64	Tomark Adjustable Corneal Marker	47	UNO Colorline Vitreous Forceps, Hattenbach	219, S24
Stitch Scissors	82, 92	TOMARK Adjustable Corneal Marker	56	UNO Colorline Vitreous Scissors	218, S21, S23, S24
Stitch Scissors, Spencer	92	TOMARK Pendulum Instrument	47, S6	U-Shaped Probe	186
Stopper	211	Tonometer	48, S12	U-Shaped Probe, Thomas Neuhann	186
Strabismus Hook	173-174, S13-S16, S18-S19, S21, S23-S24	Tonometer, Schioetz	48, S12		
Strabismus Hook, Bonn Model	174, S13, S19, S21, S23-S24	Towel Forceps	136, S12, S19, S21, S23-S24	V	
Strabismus Hook, Graefe	173, S13-S14, S16, S18-S19	Towel Forceps, Backhaus	136	Visual Axis Marker	45
Strabismus Hook, Hamburg Model	173	Towel Forceps, Schaedel	136, S12, S19, S21, S23-S24	Visual Axis Marker, Grene	45
Strabismus Hook, Jameson	174, S13, S15, S21, S23-S24	Trabecular Aspiration Cannula	S12	Vitreoretinal Forceps	129-132, S21, S23-S24
Strabismus Hook, Malbran	174, S18-S19	Trabecular Aspiration Cannula, Ph. Jacobi	S12	Vitreoretinal Forceps 20 gauge	130, 131
Strabismus Marker	44	Trabecular Probe	186, S12	Vitreoretinal Forceps 23 gauge	132
Strabismus Marker, Fink	44	Trabecular Probe, Harms	186, S12	Vitreoretinal Forceps 25 gauge	132
Strabismus Muscle Tucker	42, 172, S18	Trabecular Probe, Mackensen	186	Vitreoretinal Forceps, A. Mueller	130
Strabismus Muscle Tucker, Cueppers	42, 172, S18	Trabeculectomy Marker	44	Vitreoretinal Forceps, Naito	130, 133
Strabismus Scissors	90-91, S15	Trabeculectomy Marker, Fronimopoulus	44	Vitreoretinal Forceps, Reimann	129
Strabismus Scissors, Bonn Model	90	Trabeculectomy Punch	76, S12	Vitreoretinal Forceps, Tano	130, 132-133, 219, S21, S23-S24
Sub-Tenon Injection Cannula	193	Trabeculectomy Punch, Luntz-Dodick	76, S12	Vitreoretinal Scissors	95-96, S21, S23-S24
Sub-Tenon Injection Cannula, Nagata	193	Trabeculum Scissors	83	Vitreoretinal Scissors 20 gauge	95
Superior Rectus Forceps	125-126, S2, S12, S19	Trepbine	69-70, 215-216, 243, S7		
Superior Rectus Forceps, Castroviejo	126	Trepbine, Best	69		
Superior Rectus Forceps, Elschnig	126	Triangle Scleral Depressor	180		
Superior Rectus Forceps, Leydhecker	126				

KEYWORD PAGE

Vitreoretinal Scissors 23 gauge.....	96
Vitreoretinal Scissors 25 gauge.....	96
Vitreoretinal Scissors, Heimann.....	95 - 96, S21, S23 - S24
Vitreoretinal Scissors, Kroll.....	95
Vitreoretinal Scissors, Naito.....	95
Vitreoretinal Scissors, Roider.....	95
Vitreous Cutter.....	222

W

Wire Loop.....	88
----------------	----

NAME	PAGE	NAME	PAGE	NAME	PAGE
A					
Agarwal.....	152	Cauer.....	134	Graefe.....	93, 109, 173, S13-S14, S16, S18-S19
Akahoshi.....	120-21, 154, S3	Chang.....	210	Grehn.....	76, S12
Althaus / Cartsburg.....	202, S11	Chansue.....	169-170, S8	Greite.....	196
American Model.....	38	Chayet.....	45, S7	Grene.....	45
Amon.....	159	Cibis.....	245	Guell.....	170
A. Mueller.....	130	Citelli.....	75	Guthrie.....	38, S16
Anel.....	189, S17	Cook.....	23-24, S13-S14		
Anis.....	160	Cranz.....	191, S3, S13		
Arlt.....	161, S2, S4	Cueppers.....	42, 172, S18		
Arruga.....	144				
Auffarth.....	115, S3, S5	D			
Axenfeld.....	37, 52, S15-S17	Dardenne.....	65, 107, 157, S2	Haefliger.....	152
Ayer.....	134	D. D. Koch.....	159	Halstead.....	135, S2, S9, S12-S14, S17-S19, S21, S23-S24
		D.D. Koch - T. Kohnen.....	198, 199	Hamburg Model.....	63, 69, 101, 162, 173, 210, S19
B		Desmarres.....	35-36, 134, S13-S16, S18	Hansen.....	65
Backhaus.....	136	Dick.....	26, 122, 169, 202, S5	Harms.....	186, 245, S12
Bangerter.....	23, 107, 126, 166, 190, S14, S18-S19	Dieffenbach.....	136, S2, S18-S19, S21, S23-S24	Hartmann.....	135, S16, S21, S23, S24
Barraquer.....	15-19, 87, 101, 104, 109, 139-141, 145, S2-S4, S9-S10, S12-S15, S19	Dithmar.....	66	Hattenbach.....	86, 108, 117, 143, 149, 175, 180, 219, 221, 264, S4, S6, S21, S23
Barraquer Frankfurt Model.....	16	Doden.....	121, S2, S4	Hattenbach / Nikolic.....	221, S23
Barraquer-Oosterhuis.....	18	Draeger.....	102, S16, S17-S18, S21, S23-S24	Hebra.....	52
Barraquer-Troutman.....	140	Drysdale.....	155	Heidelberg Model.....	109, 111, 144, 211, 242, S17
Bartl.....	22	Duisburg Model.....	209, S21	Heimann.....	95-97, 128, S21, S23-S24
Bartz-Schmidt.....	175, S21	Duncker.....	42, S9	Hennig / Friedrich.....	155
Bauer.....	S7			Henry.....	101
Beer.....	101	E		Hess-Engelking.....	112
Behrens-Baumann.....	162	Elsa.....	243	Hirschmann.....	157
Berke.....	, 134	Elschnig.....	126, 164	Hoffer.....	41
Best.....	69	Elze.....	112	Hoffmann.....	169, S5
Beyer.....	75	Engels.....	44, 159	Homburg Model.....	48
Binkhorst.....	163, 194	Essen Model.....	127, 166	H.-R. Koch.....	17, 85, 88, 114-116, 124, 196, 200, S4
Blum.....	170, S7			Hunold.....	205
Bonn Model.....	63, 90, 92, 106, 110-111, 158, 172, 174, 190, 191, 194, S2-S4, S9, S12-S13, S15, S19, S21, S23-S24	F		I	
Bopp.....	191, S10	Fechner.....	121	Ikeda.....	116, 116
Bowman.....	185, S17	Fine.....	82, 149-150, S7-S8	Inamura.....	115
Boynton.....	144	Fine-Thornton.....	149-150, S7-S8		
Braunstein.....	43	Fink.....	44, S18	J	
Brauweiler.....	85, 196, 198, 200, 201	Flieringa.....	150, S9	Jacobi.....	156, 193, S9, S12
Bremen Model.....	156	Foerster.....	65, 102	Jaeger.....	36, S15, S16
Breyer.....	169, 170, 171, S5, S8	Franceschetti.....	69	Jameson.....	174, S13, S15, S21, S23-S24
Bunge.....	52, S14	Frankfurt Model.....	16, 156, 158, 165, 192, 225, S2, S12-S13	Japan Model.....	85, 116, S5
Buratto.....	119, S7, S8	Freiburg Model.....	102	Juenemann.....	186, 190, 265, S17
Burk.....	193	Friedmann.....	242, S17		
Busse.....	192	Fronimopoulos.....	44, S12	K	
				Kalt.....	143
C					
		Gaskin.....	108, S4	Kammann.....	29
Castroviejo.....	23, 43, 81, 83, 103, 126, 142, 145, 163, 164, S2-S3, S9, S12-S19, S21, S23-S24	Gass.....	77	Kampik.....	115, 207
		Gerten.....	46, 117, S6, S11	Katzin.....	81, S9, S10
		Giessen Model.....	36, 44, 167, 183, S17	Kaufmann.....	102
		Goldstein.....	38	Keer.....	198
		Gradle.....	101		

NAME	PAGE	NAME	PAGE	NAME	PAGE
Kellnar.....	184	Meyhoefer.....	51, S15-S16	Schaedel.....	136, S12, S19, S21, S23-S24
Kelly.....	76, S12	Meyhoeffer.....	S15	Schioetz.....	48, S12
Kelman-McPherson.....	108, S8	Miyoshi.....	151	Schmack.....	198, 199
Kerrison.....	75, S17	Mohr.....	77	Schmidt-Mumm.....	69, S12
Kershner.....	25, S9, S14, S19, S21, S23-S24	Morcher.....	241, S4	Schocket.....	180, S15
Kilp.....	244, S21, S23-S24	Morris-Kroll.....	175	Schott.....	31, 32
Kirchhof.....	229	Mueller.....	38, 130, 194, S17	Scott.....	46, S6
Klein.....	208	Muenster Model.....	172, S17	Sekundo.....	123, 160
Klöckner.....	193	Munich Model.....	103	Shaaf.....	101, S13, S15
Knorz.....	26, 118, S7, S8	Murdock.....	24-25	Shepard.....	122
Koch-Kohnen.....	S3, S5	N		Simcoe.....	194, 226, S2
Koga / Noda.....	28	Nagahara.....	152-154, S3, S5	Sinskey.....	156, 171-172, S2-S4
Kohnen.....	154, 172, 198	Nagata.....	193	Snellen.....	161
Kratz.....	17-19, 195, 226, S2	Naito.....	95, 97, 130, 133	Snyder-Osher.....	88, 125, S4, S6
Kratz-Barraquer.....	17, 19	Noda.....	27, 28, 29	Solomon.....	26
Krefeld Model.....	102	Noyes.....	86	Spaleck.....	123
Kreissig.....	179	Ohno.....	17	Spencer.....	92
Kroll.....	95, 97, 175	Örk Linz Model.....	71	Spirig.....	158, 161, 192
Krumeich.....	42, 186, 197, S10	Osher.....	43, 86, 88, 125, S4, S6	Steinmetz.....	87, 158, S4
Kuhn.....	246, S21	O		Stevens.....	89, S3, S13-S14, S16-S19, S21, S23-S24
Kuhnt.....	64, S9, S14, S16	Ohno.....	17	Strobel.....	125
Kyu.....	195	Örk Linz Model.....	71	Swiss Model.....	111, 112
L		Osher.....	43, 86, 88, 125, S4, S6	Szurman.....	64, 71, 105, 117, 168, 203, 223, S11
Lambert.....	, S16, 134	P		T	
Landers.....	244	Papst.....	195, 205, 264	Takeda.....	19
Laqua.....	179, 207, S19	Paufique.....	64, 105, S9, S14, S21, S23-S24	Tandogan.....	157
Lester.....	157	Paulmann.....	27	Tano.....	130, 132-133, 219, S21, S23-S24
Lewicky.....	225	Peyman.....	244	Teping.....	195
Leydhecker.....	126	Pfaeffl.....	118, 171, S7-S8	Thomas Neuhann.....	19, 41-42, 70, 104, 112, 118, 151, 153-155, 167-168, 186, 202, S3, S5, S8-S11
Liebermann	25, 28, 154, S2-S3, S7-S8, S12	Pfaeffl-Breyer.....	171	Thornton.....	149-150, S7-S8
Liebermann-Tennant.....	25, 28, S2-S4, S7-S8, S10, S12-S13	Pfennigwerth.....	127	T. Kohnen.....	122, 154, 198
Lindstrom.....	42	Ph. Jacobi.....	193, S12	T. Neuhan.....	123
Lohmann.....	243, S7, 64	Plange.....	55, 56	Tobias Neuhann.....	44, 116-117, 152
Lordan.....	134, S16	Pretori.....	21	Tolentino.....	244
Lu.....	45, S7	R		Troeltsch.....	127
Lucke.....	175, 179	Rau.....	168	Troutman.....	81, 140, 145, S2
Lund.....	158	Reimann.....	129	Tuebingen Model ...	84, 106, 192, 225, S2-S3, S10, S12-S13, S16-S19, S21, S23-S24
Luntz-Dodick.....	76, S12	Remky.....	105, 125, S2, S11-S12, S19		
M		Rentsch.....	156, 159-160, 194, S6, S10		
Machemer.....	244	R. Koch.....	201, 202		
Mackensen.....	82, 186	Rohrschneider.....	144		
Mainz Model.....	155	Roider.....	95, 210		
Malbran.....	174, S18-S19	Rollet.....	37		
Maloney.....	45, S9	Rupal Shah.....	119, S8		
Masket.....	113	Rycroft.....	225		
Maus.....	30	S			
Mead.....	242, S17	Sattler.....	127, S18-S19		
Mellinger.....	20	Sauer.....	19		
Mellinger-Beard.....	20	Sautter.....	81, 162, 192, 225, S3-S5, S9, S12		
Mendez.....	45-46, S7	W			
Meyer-Schwickerath.....	180	Walser.....	245, S14		

NAME INDEX

NAME PAGE

Watzke	127, S19
W. Bauer.....	196
Weber.....	161, 210, S21
Wecker.....	86, 87, 165
Weiss.....	21, 22
Wells.....	167, S14
West.....	243
Westcott.....	89, S2, S9, S12-S13, S16, S18-S19, S21, S23-S24
Wilder.....	183, S15, S17
Wright.....	243

Y

Yae	85, S5
Yörük.....	117

Z

Zirm	30, 151, 157, 167, 197, S7-S8
Zuerich Model.....	104



Geuder®
Precision made in Germany

GEUDER AG
Hertzstrasse 4
69126 Heidelberg
Germany

Phone: +49 6221 3066
Fax: +49 6221 303122
info@geuder.de
www.geuder.com