

 **LABTICIAN™**
LAB™EDGE

Labtician Ophthalmics **Lab™EDGE** Catalogue



LAB™EDGE KNIVES & BLADES PRODUCTS

Our **Lab™EDGE** portfolio of knives and blades are offered in a wide variety of sizes to meet all your needs during cataract surgery. **Lab™EDGE** is designed to offer best-in-class sharpness, precise cutting dimensions, and consistent cutting force through the use of electrochemical sharpening and grindless technology. This ensures the cutting edges of the blades are smooth, absent of residual burrs, and creates a true incision size that can be replicated consistently from blade to blade.

FEATURES OF THE Lab™EDGE SINGLE-USE KNIVES	1
STAB	2
Stab And Paracentesis	2
Controlled Depth	2
CRESCENT	3
Straight	3
Angled	3
PHACO SLIT	4
Single Bevel Angled	4
Double Bevel Angled	4
MVR/SIDEPORT	5
Straight	5
Angled	5
IMPLANT	5
Angled	5

Continued next ▶

OTHER 6

Scleral 6

Lamellar 6

Needle 6

MINI BLADES 7

Chisel Tip 7

Spear Tip 7

Round Tip 7

Lance Tip 8

Curved Tip 8

Round Tip Straight 8

Round Tip Angled 8

HANDLES 9

Cross-cut collet 9

Round 9

Hex 10

www.labtician.com/docs


















FEATURES OF THE LAB[™]EDGE SINGLE-USE KNIVES

The Lab[™]EDGE knives have best-in-class sharpness* that reduces the force required for penetration, creating a precise wound architecture and minimizing surface damage.

Lab[™]EDGE knives are manufactured using electrochemical sharpening and grindless technology to eliminate residual burrs on the cutting edges, producing a true incision size that can be replicated consistently from blade to blade.



For greater incisional accuracy, tip size is indicated in the description and on the handle.

* Testing was performed on a Mark-10 Force Testing gauge and motorized stand that measures force over time.

STAB STAB AND PARACENTESIS



Full Handle Knife


Item No.	Description
LAB961501*	15 °
LAB962201*	22.5 °
LAB963001*	30 °
LAB964501*	45 °

Screw-in Blade

Item No.	Description
LAB961502*	15 °

6 pcs. per package

STAB CONTROLLED DEPTH



Full Handle Knife

Item No.	Description
LAB7514*	15° 3.0 mm depth Light blue
LAB7516*	15° 5.0 mm depth Green
LAB7523*	22.5° 4.0 mm depth Grey
LAB7531*	30° 3.5 mm depth Black
LAB7536*	30° 5.0 mm depth Black

Screw-in Blade

Item No.	Description
LAB7511	15° 1.5 mm depth Orange
LAB7513	15° 3.0 mm depth Light blue
LAB7515	15° 5.0 mm depth Green
LAB7522*	22.5° 4.0 mm depth Grey
LAB7530*	30° 3.5 mm depth Black

6 pcs. per package

*The lead times to get this product are variable. Please contact your sales representative for more information.

CRESCENT STRAIGHT



Full Handle Knife

Item No. Description

LAB960021* 2.0 mm

6 pcs. per package

CRESCENT ANGLED



Full Handle Knife

Item No. Description


LAB940002* 2.5 to 2.0 mm Tapered

LAB960022 2.0 mm

6 pcs. per package

*The lead times to get this product are variable. Please contact your sales representative for more information.

PHACO SLIT SINGLE BEVEL ANGLED




Full Handle Knife

Item No.	Description	Item No.	Description
LAB961861*	1.8 mm	LAB962751	2.75 mm
LAB962061*	2.0 mm	LAB962861*	2.80 mm
LAB962261	2.2 mm	LAB962851*	2.85 mm
LAB962461*	2.4 mm	LAB962961*	3.0 mm
LAB962561*	2.5 mm	LAB963261	3.2 mm
LAB962661*	2.65 mm	LAB963561*	3.5 mm

6 pcs. per package

PHACO SLIT DOUBLE BEVEL ANGLED




Full Handle Knife

Item No.	Description	Item No.	Description
LAB961868*	1.8 mm	LAB962758*	2.75 mm
LAB962068*	2.0 mm	LAB962868*	2.80 mm
LAB962268*	2.2 mm	LAB962858*	2.85 mm
LAB962468*	2.4 mm	LAB962968*	3.0 mm
LAB962568*	2.5 mm	LAB963268*	3.2 mm
LAB962668*	2.65 mm		

6 pcs. per package

*The lead times to get this product are variable. Please contact your sales representative for more information.

MVR/SIDEPORT STRAIGHT



Item No.	Description	Item No.	Description
LAB965564*	25 ga, 0.6mm incision	LAB965560*	20 ga, 1.2mm incision
LAB965562*	23 ga, 0.8 mm incision	LAB965561	19 ga, 1.4 mm incision
LAB965563*	21 ga, 1.1mm incision	LAB965550*	1.0mm incision

6 pcs. per package


MVR/SIDEPORT ANGLED



Item No.	Description	Item No.	Description
LAB965569*	25 ga, 0.6mm incision	LAB965565*	20 ga, 1.2mm incision
LAB965567*	23 ga, 0.8 mm incision	LAB965566*	19 ga, 1.4 mm incision
LAB965568*	21 ga, 1.1mm incision	LAB965570*	1.0mm incision

6 pcs. per package

IMPLANT ANGLED




Item No.	Description
LAB965261*	5.2 mm

6 pcs. per package

*The lead times to get this product are variable. Please contact your sales representative for more information.

OTHER SCLERAL



Full Handle Knife

Item No.	Description
LAB5710*	Standard US style
LAB5710S*	Small European style

Screw-in Blade

Item No.	Description
LAB5700	Standard US style
LAB5700S*	Small European style

6 pcs. per package

OTHER LAMELLAR



Screw-in Blade

Item No.	Description
LAB6600	Angled, bevel up

6 pcs. per package

OTHER NEEDLE



Screw-in Blade

Item No.	Description
LAB5910	3.0 mm
LAB5920*	5.0 mm

6 pcs. per package

*The lead times to get this product are variable. Please contact your sales representative for more information.

MINI BLADES CHISEL TIP



Screw-in Blade

Item No.	Description
LAB6100*	Chisel 1.5 mm
12 pcs. per package	

MINI BLADES CHISEL TIP



Screw-in Blade

Item No.	Description
LAB6200*	Chisel 2.7 mm
12 pcs. per package	

MINI BLADES SPEAR TIP



Screw-in Blade

Item No.	Description
LAB6300*	Spear tip
6 pcs. per package	

MINI BLADES ROUND TIP



Screw-in Blade

Item No.	Description
LAB6400	Round tip, sharp one side
12 pcs. per package	

*The lead times to get this product are variable. Please contact your sales representative for more information.

MINI BLADES LANCE TIP



Screw-in Blade

Item No. Description

LAB6500 Lance tip

12 pcs. per package

MINI BLADES CURVED TIP



Screw-in Blade

Item No. Description

LAB6700 Curved tip, sharp one side

12 pcs. per package

MINI BLADES ROUND TIP STRAIGHT




Screw-in Blade

Item No. Description

LAB6900 Round tip, sharp all around, straight

6 pcs. per package

MINI BLADES ROUND TIP ANGLED



Screw-in Blade

Item No. Description

LAB6910* Round tip, sharp all around, angled

6 pcs. per package

*The lead times to get this product are variable. Please contact your sales representative for more information.

HANDLES CROSS-CUT COLLET



1305

Screw-in Blade

Item No.	Description
LAB1305*	Cross-cut collet
1 pcs. per package	

HANDLES HEX AND ROUND



1311



1312



1310



1350



1360



1320



1330

Screw-in Blade

Item No.	Description
LAB1311*	5.0 cm round
LAB1312*	7.5 cm round
LAB1310*	10.0 cm round
LAB1350*	13.0 cm round
LAB1360*	13.0 cm round heavy duty
LAB1320*	15.5 cm round
LAB1330*	10.0 cm hex
1 pcs. per package	

*The lead times to get this product are variable. Please contact your sales representative for more information.

Labtician Ophthalmics, Inc.

2140 Winston Park Drive, Unit 6, Oakville, Ontario, Canada L6H 5V5

Telephone: **905.829.0055** Toll free: **1.800.265.8391** (Canada Only)

labtician.com

Order Desk & Customer Service

Toll free: **1.855.651.4934** (Canada & USA only)

Toll free fax: **1.855.390.9790** (Canada & USA only)

Email: orders@labtician.com

